

Product Information

Fluorescamine

≥ 98% (TLC), powder, used for detection of primary amines

F9015

Product Description

CAS Number: 38183-12-9

Synonym: 4-Phenylspiro-[furan-2(3H),1-phthalan]-

3,3'-dione

Molecular Weight: 278.26 Molecular Formula: C₁₇H₁₀O₄ Melting Point: 154-155 °C¹

 $\lambda_{max} \colon 306 \text{ nm}, \, 284 \text{ nm}, \, 276 \text{ nm}, \, 235 \text{ nm}$

Extinction Coefficient: $E^{mM} = 3.8, 4.1, 3.9, 25.9$

(ether)

Excitation Wavelength: 390 nm (borate buffer,

pH 8.5)^{2,3}

Emission Wavelength: 465 nm;2 475 nm3

Fluorescamine, a heterocyclic dione, reacts with primary amines to form a fluorescent product. Free NH₃ yields a non-fluorescent product. The fluorescence of a solution that contains protein plus fluorescamine is proportional to the quantity of free amine groups present.² This is the basis of a fluorescent protein assay.^{4,5} This product has been used to label casein, to give a derivatized casein substrate for measuring protease activity.⁶

Fluorescamine is used in many sensitive detection methods, such as:

- Characterization of poly-L-lysine (pLL) / DNA complexes post-modified with a multivalent hydrophilic polymer⁷
- Spectrofluorimetric analysis of procaine⁸
- Detection and quantitation of residual aminopenicillins by HPLC after fluorescamine derivatization⁹
- Determination of lisinopril in human plasma and urine by HPLC with fluorescence detection¹⁰
- Determination of sulfonamides in honey by HPLC with fluorescence detection¹¹

Several theses¹²⁻¹⁶ and dissertations^{17,18} have cited use of product F9015 in their protocols.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

This product may be stored at room temperature.

Solubility

1

This product is tested for solubility in acetone at 50 mg/mL. Fluorescamine is also soluble in DMSO^{19,20} and other anhydrous, aprotic solvents.

A stock solution of 7.5 mg of fluorescamine in 25 mL of acetone is stable at room temperature, if kept free of moisture. Fluorescamine hydrolyzes quite rapidly in water to give non-fluorescent products. The half-life for the reaction with peptides is 10-100 milliseconds, with hydrolysis taking 1-10 seconds.



References

- 1. The Merck Index, 12th ed., Entry# 4193 (1996).
- 2. Gore, M. (ed.), Spectrophotometry and Spectrofluorimetry: A Practical Approach. Oxford University Press (New York, NY), p. 63 (2000).
- Rickwood, D. (ed.), Centrifugation: A Practical Approach. Oxford University Press (Oxford, UK) (1984).
- 4. Böhlen, P. et al., Arch. Biochem. Biophys., **155(1)**, 213-220 (1973).
- Udenfriend, S. et al., Science, 178(4063), 871-872 (1972).
- 6. Sogawa, K., and Takahashi, K., *J. Biochem.*, **83(6)**, 1783-1787 (1978).
- Read, M.L. et al., FEBS Lett., 461(1-2), 96-100 (1999).
- 8. Segura Carretero, A. *et al.*, *J. Pharm. Biomed. Anal.*, **21(5)**, 969-974 (1999).
- 9. Hong, C., and Kondo, F., *J. Food Prot.*, **63(10)**, 1421-1425 (2000).
- 10. Sagirli, O., and Ersoy, L., J. Chromatogr. B Analyt. Technol. Biomed. Life Sci., **809(1)**, 159-165 (2004).
- 11. Maudens, K.E. *et al.*, *J. Chromatogr. A*, **1047(1)**, 85-92 (2004).
- 12. Gentes, Marie-Line, "Health Assessment of Tree Swallows (*Tachycineta bicolor*) Nesting on the Athabasca Oil Sands, Alberta". University of Saskatchewan, M.Sc. thesis, p. 112 (2006).
- Olsgard, Mandy Lee, "Toxicological evaluation of inhalation exposure to benzene and toluene in a raptorial bird, the American kestrel, *Falco sparverius*". University of Saskatchewan, M.Sc. thesis, p. 112 (2007).
- 14. Esposito, Tullio Vito Francesco, "Assessment of the Biodistribution of Aluminum-Based Vaccine Adjuvants using 111/115 In-AlO(OH)". University of British Columbia, M.Sc. thesis, p. 11 (2016).
- 15. Bhattarai, Rewati Raman, "Effect of Food Structure on Macronutrients Digestion and Fermentation Kinetics". University of Queensland, M.Sc. thesis, p. 36 (2018).
- 16. Ekmyr, Elina Toss, "Total proteolysis in bovine raw bulk milk of Northern Sweden". Swedish University of Agricultural Sciences, M.Sc. thesis, p. 29 (2020).

- 17. Brammell, Benjamin Frederick, "Characterization of pollutant response in teleosts with varying degrees of pollutant sensitivity". University of Kentucky, Ph.D. dissertation, p. 152 (2005).
- 18. Drzewiecki, Kathryn Emily, "Collagen Methacrylamide A Photocrosslinkable, Thermoreversible Collagen-Based Biomaterial: Characterization and Applications". Rutgers, The State University of New Jersey, Ph.D. dissertation, p. 52 (2016).
- 19. Froehlich, P.M., and Cunningham, T.D., *Analytica Chim. Acta*, **84(2)**, 427-430 (1976).
- 20. Toro, T.B., and Watt, T.J., *Protein Sci.*, **24(12)**, 2020-2032 (2015).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at <u>SigmaAldrich.com/terms</u>.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

F9015pis Rev 07/22 NDH,PHC,GCY

