

Stick with PureProteome™ Kappa & Lambda Ig-binding Magnetic Beads For Fab-ulous human immunoglobulin capture!

Producing Fab fragments for therapeutic development?

Immunoprecipitating protein:Ig complexes?

Isolate more immunoglobulins, faster, with PureProteome™ magnetic beads. Our newest beads bind either the Kappa or Lambda light chains of antibodies, offering:

- High binding capacity (>2 mg/mL slurry)
- Workflow improvements: save time, no centrifugation, easier handling
- Affinity for ALL Ig subtypes – bind IgG, IgA, IgM, IgE, and IgD (Proteins A and G only bind IgG)

	% Ig Captured
IgA	99.79
IgG ₁	99.94
IgG ₂	99.91
IgG ₃	99.98
IgG ₄	99.83
IgD	99.99
IgE	99.11
IgM	99.92


Table 1. Efficient capture of all Ig subtypes from human serum. A blend of 150 µL of each of PureProteome™ Kappa magnetic bead slurry and PureProteome™ Lambda magnetic bead slurry was used to capture (deplete) Ig from 25 µL human serum (pooled). The % Ig captured was determined by measuring Ig in the input and unbound serum samples by ELISA (IgD and IgE) and MILLIPLEX® MAP assays (IgA, IgG and IgM, Catalogue No. HGAMMAG-301K).

SPECIAL INTRODUCTORY OFFER

Save 25%
on PureProteome™ Kappa or Lambda
Ig-Binding Magnetic Beads


Use Promotion Code: ZBIOS543
Offer valid until November 30, 2012.

*See reverse side for special offer details



Purify the Fab portion away from the Fc portion of Immunoglobulins.

- PureProteome™ Kappa & Lambda magnetic beads specifically bind to Fab fragment light chains, enabling you to separate Fab and/or F(ab')₂ from Fc fragments.
- Benefits of separating Fab and Fc fragments include: reduced immunogenicity in *in vivo* studies and less interference in samples containing Fc-binding molecules.



High purity, Fab-ulous recovery.




Figure 1. High purity isolation of human Lambda Fab from hybridoma medium. 100 μ L of PureProteome™ Lambda Ig-binding magnetic bead slurry was used to capture the Fab fragment from 100 μ L of Hybridoma Medium spiked with human Lambda Fab. The captured protein was eluted 3 times using 100 μ L of 100 mM triethylamine pH 11.5. 10 μ L of each fraction was resolved by SDS-PAGE and Coomassie stained.

Ordering Information

Description	Qty/PK	Catalogue No.
PureProteome™ Kappa Ig-Binding Magnetic Beads	2 mL	LSKMAGKP02
PureProteome™ Lambda Ig-Binding Magnetic Beads	2 mL	LSKMAGLM02



www.merckmillipore.com/offices

Merck Millipore, the M mark and PureProteome are trademarks of Merck KGaA, Darmstadt, Germany. MILLIPLEX is a registered trademark of Merck KGaA, Darmstadt, Germany. Trademarks belonging to third parties are the properties of their respective owners.
Lit No. PR2277ENEU LS-SBU-12-06697 6/2012 Printed in the USA.
© 2012 EMD Millipore Corporation, Billerica, MA USA. All rights reserved.

Save 25%
on PureProteome™ Kappa or Lambda
Ig-Binding Magnetic Beads

To take advantage of this offer, contact your local Customer Service office and reference promotional code: ZBIOS543*.

*Promotion code ZBIOS543 must be used when placing your order to be eligible for discount. May not be combined with any other offers. Void where prohibited by law or company policy. Offer valid until November 30, 2012.

NOTE: For human immunoglobulins only.

To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 01805 045 645

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: www.merckmillipore.com/offices

For Technical Service visit:

www.merckmillipore.com/techservice

Get Connected!

Join Merck Millipore Bioscience on your favorite social media outlet for the latest updates, news, products, innovations, and contests!

 facebook.com/MerckMilliporeBioscience

 twitter.com/Merck4Bio