

KAPA Taq HotStart PCR Kit

KR0355_S - v3.20

Product Description

KAPA Taq DNA Polymerase is the single-subunit *Taq* DNA polymerase of the thermophilic bacterium *Thermus aquaticus*, purified from recombinant *Escherichia coli*. KAPA Taq DNA Polymerase has $5'\rightarrow 3'$ polymerase and $5'\rightarrow 3'$ exonuclease activity, but no $3'\rightarrow 5'$ exonuclease (proofreading) activity. The enzyme system has an error rate of approximately 1 error per 2.2 x 10^5 nucleotides incorporated.

In the HotStart formulation, the enzyme is combined with a proprietary antibody that inactivates the enzyme until the first denaturation step. This prevents nonspecific amplification during reaction setup, increases sensitivity, and improves reaction efficiency. PCR products generated with KAPA Taq HotStart are A-tailed and may be cloned into TA cloning vectors.

KAPA Taq HotStart Buffer is a uniquely-formulated buffer to facilitate specific primer annealing. This translates to higher yields of specific product when compared to traditional *Taq* buffers, and improved amplification of GC- and AT-rich templates. However, KAPA Taq HotStart DNA Polymerase may be used in combination with any standard *Taq* buffer with a pH of 8.3 or higher.

Product Applications

The KAPA Tag HotStart PCR Kit is ideally suited for:

- Routine PCR
- Multiplex PCR
- Amplification of DNA for Sanger sequencing
- Any standard PCR application for which a hotstart formulation of a high-quality thermostable DNA polymerase is required.

Product Specifications

Shipping and Storage

KAPA Taq HotStart PCR kits are shipped on dry ice or ice packs, depending on the country of destination. Upon arrival, store kit components at -15°C to -25°C in a constant-temperature freezer. When stored under these conditions and handled correctly, full activity of the kit is retained until the expiry date indicated on the kit label. KAPA Taq HotStart Buffer contains isostabilizers and may not freeze solidly, even when stored at -15°C to -25°C. This will not affect the shelf-life of the product.

Handling

Always ensure that the product has been fully thawed and mixed before use. Reagents may be stored at 2°C to 8°C for short-term use (up to 1 month). Return to -15°C to -25°C for long-term storage.

Kit Codes and Components			
KK1508 (250 U)			
KK1510 (500 U)	KAPA Taq HotStart DNA Polymerase (5 U/μL) KAPA Taq HotStart Buffer (5X) MgCl ₂ (25 mM)		
KK1513 (2500 U)	go., ₂ (=0)		
KK1509 (250 U)	KAPA Taq HotStart DNA Polymerase (5 U/μL) KAPA Taq HotStart Buffer (5X)		
KK1511 (500 U)	MgCl ₂ (25 mM) KAPA dNTP Mix (10 mM each)		
KK1512 (2500 U)	KAPA Taq HotStart DNA Polymerase (5 U/μL)		

Quick Notes

- KAPA Taq DNA Polymerase can replace any commercial Taq DNA polymerase in an existing protocol. The final MgCl₂ concentration may need to be optimized to account for differences in buffer formulation.
- KAPA Taq Buffers contain MgCl₂ at a final concentration of 1.5 mM.
- Buffer A is recommended as first approach and for applications requiring high yields.
- Buffer B is recommended for applications where high sensitivity is required (e.g. when the template is limiting).
- Both buffers may be evaluated to determine the buffer most suitable for a specific application.
- The KAPA Taq PCR system is suitable for the amplification of fragments up to 3.5 kb from genomic DNA or 5 kb from less complex targets.

Quality Control

Each batch of KAPA Taq HotStart DNA Polymerase is confirmed to contain <2% contaminating protein (Agilent Protein 230 Assay). KAPA Taq HotStart PCR kits are subjected to stringent quality control tests, are free of contaminating exo- and endonuclease activity, and meet strict requirements with respect to DNA contamination levels.

Safety Information

Precautions

- Handle all samples as if potentially infectious, using safe laboratory procedures. As the sensitivity and titer of potential pathogens in the sample material can vary, the operator must optimize pathogen inactivation and follow the appropriate measures according to local safety regulations.
- Do not eat, drink, or smoke in the laboratory area.
- Do not pipette by mouth.
- Wear protective disposable gloves, laboratory coats, and eye protection, when handling samples and kit reagents.
- Wash hands thoroughly after handling samples and reagents.

Waste handling

- Discard unused reagents and waste in accordance with country, federal, state, and local regulations.
- Safety Data Sheets (SDS) are available online at <u>www.sigmaaldrich.com</u>, or upon request from <u>www.sigma-aldrich.com/techservice</u>.

KAPA Tag HotStart PCR Protocol

KAPA Taq HotStart DNA Polymerase can be used to replace any commercial *Taq* DNA polymerase in an existing protocol. To allow the most seamless integration of KAPA Taq HotStart into existing protocols, be sure to match reaction conditions, particularly the MgCl₂, primer and enzyme concentrations, as closely as possible.

Step 1: Prepare the PCR master mix

- Ensure that all reagents are properly thawed and mixed.
- Prepare a PCR master mix containing the appropriate volume of all reaction components common to all or a subset of reactions to be performed.
- Calculate the required volumes of each component based on the following table:

Component	25 μL reaction¹	Final conc.
PCR-grade water	Up to 25 μL	N/A
5X KAPA Taq HotStart Buffer	5.0 μL	1X
25 mM MgCl ₂	1.5 μL	1.5 mM²
10 mM dNTP Mix	0.5 μL	0.2 mM each
10 μM Forward Primer	0.5–1.25 μL	0.1–0.25 μM
10 μM Reverse Primer	0.5–1.25 μL	0.1–0.25 μM
5 U/µL KAPA Taq HotStart DNA Polymerase ³	0.1 μL	0.5 U
Template DNA ⁴	As required	As required

 $^{^1}$ Reaction volumes of 10–50 μL are recommended. For volumes other than 25 μL , scale reagents proportionally.

 $\mbox{NOTE:}$ For GC-rich or other difficult templates or amplicons, include DMSO at a final concentration of 5%.

 $^{^2\}text{A}$ final MgCl $_2$ concentration of 1.5 mM is sufficient for most standard applications. For assays that do not perform well with 1.5 mM MgCl $_2$, the optimal MgCl $_2$ concentration for each primer–template combination should be determined empirically.

 $^{^3}$ For GC-rich and other difficult templates, higher enzyme concentrations (up to 2.5 U per 25 μL reaction) may be required.

 $^{^{4}}$ $\leq\!250$ ng for genomic DNA; $\leq\!25$ ng for less complex DNA (e.g. plasmid, lambda).

Step 2: Set up individual reactions

- Transfer the appropriate volume of PCR master mix, template and primer to individual PCR tubes/wells of a PCR plate.
- Cap or seal individual reactions, mix and centrifuge briefly.

Step 3: Run the PCR

• Perform PCR with the following cycling protocol:

Step	Temperature	Duration	Cycles
Initial denaturation	95°C	3 min¹	1
Denaturation	95°C	30 sec	
Annealing ²	T _m – 5°C	30 sec	35³
Extension	72°C	1 min/kb	
Final extension (optional) ⁴	72°C	1 min/kb	1
Hold	4–10°C	8	1

 $^{^{\}rm 1}$ Initial denaturation for 3 min at 95°C is recommended for most assays. For GC-rich targets (>65% GC content), 5 min at 95°C may be used.



Manufacturing, R&D Cape Town, South Africa Tel: +27 21 448 8200

Fax: +27 21 448 6503

Technical Support sigma-aldrich.com/techservice

© 2020 KAPA is a trademark of Roche. All other product names and trademarks are the property of their respective owners.

 $^{^2}$ An annealing temperature 5°C lower than the calculated melting temperature ($T_{\rm m}$) of the primer set is recommended as first approach. If low yields and/or nonspecific amplification is obtained, an annealing temperature gradient PCR is recommended to determine the optimal annealing temperature for the primer set empirically.

 $^{^{\}rm 3}$ 35 cycles are sufficient for most assays. A higher number of cycles may be necessary for assays requiring a higher level of sensitivity.

 $^{^{\}rm 4}$ Final extension should be included if PCR products are to be cloned into TA cloning vectors.