

STAINS-ALL Molecular Biology Reagent

ProductInformation

Product No. **E7762** Lot 89H3429

[7423-31-6]

(1-Ethyl-2-(3-[1-ethylnaphtho(1,2-d)-thiazolin-2-ylidene]-2-methylpropenyl)naphtho(1,2-d)-thiazolium bromide) Mol. Wt. 559.6

Product Description

Suitable for differential staining of nucleic acids. Detected as little as 3 ng (123 BP fragment) of pBR322/Hae III DNA and 90 ng tRNA.

Preparation Instructions

Stock Stain: 0.1% (w/v) Stains-All in formamide

Staining Solution:

0.005% (w/v) Stains-All 10% (v/v) formamide 25% (v/v) isopropanol 15 mM Trizma 7 -HCl, pH 8.8 65% (v/v) H $_{2}$ O

Storage

Store at 2-8EC. Store solutions containing Stains-All protected from light.

Procedure

pBR322/Hae III DNA and tRNA were run on 10-20% polyacrylamide gradient gels. The gels were stained in the dark overnight at room temperature with freshly prepared staining solution. Destaining was accomplished by removing the gel from the staining solution and exposing it to the light from a light box until sufficient destaining had occurred. The gels were then immediately photographed. **Note**: Upon continued exposure to light the nucleic acid-dye complex will fade and disappear.

References

- 1. ISCO Staff, ISCO Tables, 10th edition, p. 35 (1989)
- 2. Dahlberg, A.E., et al., J. Mol. Biol. 41, 139 (1969)