



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### UDP-Glucuronosyltransferase 2B4 Supersomes™ human, recombinant

Product Number **U 3133**  
Storage Temperature  $-70^{\circ}\text{C}$

EC 2.4.1.17

Synonyms: UGT2B4; UDP-glycosyltransferase

#### Product Description

This UDP-glucuronosyltransferase product is a microsomal fraction of insect cells (BTI-TN-5B1-4) infected with a baculovirus strain containing the cDNA for human UGT2B4.

Glucuronidation, conjugation with glucuronic acid, plays an important role in the metabolic fate of many drugs and other xenobiotics. Examining glucuronidation by the human UDP-glucuronosyltransferase (UGT) family of enzymes is essential when investigating the metabolism of new therapeutic agents or chemicals present in the environment. This biosynthetic reaction also has a role in the conjugation and excretion of endogenous substrates, such as steroids, bilirubin, and bile acids.<sup>1</sup> UGT activity results in the conjugation of glucuronic acid to substrates containing sulfhydryl, hydroxyl, aromatic amino, or carboxylic acid moieties. The glucuronides formed are more polar (water soluble) than the parent organic substrate and are generally excreted through the kidney.

The UGT enzymes comprise a superfamily of integral membrane proteins of the endoplasmic reticulum that have been subdivided into two families, UGT1 and UGT2, based on the evolutionary divergence of their genes.<sup>2</sup> The enzymes of the UGT1A family play an important role in the metabolism of dietary constituents, phenols, and therapeutic drugs, and also the glucuronidation of bilirubin and iodothyronines. The enzymes of the UGT2B family are involved in the metabolism of bile acids, phenol derivatives, catechol-estrogens and steroids.<sup>1,3</sup> Although it is widely recognized that the liver is the major site of glucuronidation, it is now clear that UGT enzymes are also found in extra-hepatic tissues.

The product is supplied as a solution containing 100 mM Tris, pH 7.5.

UGT2B4 Specific Activity: minimum 0.3 unit/mg of protein (Lowry).

Unit Definition: One unit will transfer 1 nanomole of glucuronic acid from uridine-5'-diphosphoglucuronic acid to 7-hydroxy-4-trifluoromethylcoumarin per minute at pH 7.5 at  $37^{\circ}\text{C}$ .

#### Storage/Stability

This UGT2B4 product ships on dry ice and should be stored at  $-70^{\circ}\text{C}$ . If not using the entire contents, aliquot to minimize freeze-thaw cycles. A loss of less than 15% of the catalytic activity has been observed after the product has gone through 8 freeze/thaw cycles.

#### Procedure

##### 7-Hydroxy-4-trifluoromethylcoumarin Glucuronidation Activity:

Reactions were performed in 0.2 ml volumes containing 0.25 mg/ml of UGT2B4 protein, 10 mM magnesium chloride, 0.025 mg/ml of alamethicin, 1 mM uridine diphosphoglucuronic acid, and 50  $\mu\text{M}$  7-hydroxy-4-trifluoromethylcoumarin in 50 mM Tris, pH 7.5. After incubation for 20 minutes at  $37^{\circ}\text{C}$ , the reaction was stopped with 100  $\mu\text{l}$  of 94% acetonitrile/6% acetic acid with subsequent centrifugation at 10,000  $\times g$ .

From the supernatant, 60  $\mu\text{l}$  was injected on a 5  $\mu\text{m}$  C18 HPLC column at  $45^{\circ}\text{C}$ . There were three mobile phases: phase A (10% methanol), phase B (100% methanol), and phase C (30% acetonitrile and 1 mM perchloric acid). The initial gradient was 80% phase A, 10% phase B, and 10% phase C. Over a 15 minute period, phase B was increased to 90%, while phase C remained constant at 10%. The reaction product, 7-hydroxy-4-trifluoromethylcoumarin glucuronide (Product Code T 6410) was detected by measuring absorbance at 325 nm and quantitated against a standard.

Notes: With respect to enzyme concentration, catalysis is linear up to at least 0.8 mg of UDP-glucuronosyltransferase isozyme per ml. The glucuronidation of 7-hydroxy-4-trifluoromethylcoumarin is approximately linear for 60 minutes. Other substrates may not exhibit similar linearity.

### References

1. Tephly, T. R., et al., Metabolism of endobiotics and xenobiotics by UDP-glucuronosyltransferase. *Adv. Pharmacol.*, **42**, 343-346 (1998).
2. Mackenzie, P. I., et al., The UDP glycosyltransferase gene superfamily: recommended nomenclature update based on evolutionary divergence. *Pharmacogenetics*, **7**, 255-269 (1997).
3. Burchell, B., et al., Specificity of human UDP-glucuronosyltransferases and xenobiotic glucuronidation. *Life Sci.*, **57**, 1819-1831 (1995).
4. Matern, H., et al., Radioassay of UDP-glucuronosyltransferase activities toward endogenous substrates using labeled UDP-glucuronic acid and an organic solvent extraction procedure. *Anal. Biochem.*, **219**, 182-188 (1994).
5. Bansal, S. K., and Gessner, T., A unified method for the assay of uridine diphosphoglucuronyl transferase activities toward various aglycones using uridine diphospho[U-<sup>14</sup>C] glucuronic acid. *Anal. Biochem.*, **109**, 321-329 (1980).
6. Strassburg, C. P., et al., Expression of the UDP-glucuronosyltransferase 1A Locus in Human Colon: Identification and Characterization of the Novel Extrahepatic UGT1A8. *J. Biol. Chem.* **273**, 8719-8726 (1998).

Supersomes is a trademark of BD Gentest.

JWM/LM/MDR/MAM 8/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.