MultiScreen™ Permeability Plates

The evaluation of the reproducibility of passive, transcellular drug permeability assays

Daniel Schmidt and John Lynch, PhD Millipore Corporation, Life Science Division, Danvers, MA

Abstract

The reproducibility and precision of a high throughput method based on a published method¹ using a Millipore plate (MultiScreen Cat. MAPBMN310) for predicting passive, transcellular compound permeability was assessed. The permeability of six drugs (propranolol, methotrexate, warfarin, carbamazepine, furosemide, and testosterone) was measured on different days using five different lots of (MultiScreen) plates. In addition several protocol variations (incubation times, volumes, temperature, etc.) were tested to simulate the analytical variability that might be encountered from lab to lab and from operator to operator. Our results show that under typical conditions, this permeability assay is robust and generates reproducible data. We also observed that small protocol variations can have an effect on the apparent permeability rates of some drugs but in general the rank order (by log P_{e}) of the compounds is unaffected.

Introduction

The permeability assay is a non-cell based assay designed to predict passive, transcellular permeability of drugs in early drug discovery. The assay is carried out in a 96-well MultiScreen Permeability plate (MAPBMN310) and measures the ability of compounds to diffuse from a Donor to an Acceptor compartment separated by a hexadecane liquid layer on a polycarbonate membrane support.


After the artificial membrane has been applied to the polycarbonate membrane in the filter plate (known as the Donor plate), the 96-well Donor plate is filled with buffer solutions containing the compounds to be tested. The Donor plate is then placed upon a 96-well Acceptor plate filled with sufficient buffer so that there is liquid contact between the liquid in the Acceptor plate and the polycarbonate membrane. The Acceptor plate

can be either the plate provided with this device or an equivalent such as MSSACCEPTOR (a PTFE 96-well plate). The Donor and Acceptor plates are incubated together for 5-7 hours after which time the Donor plate is removed from the Acceptor plate. The 96 wells in the Acceptor plate can then be analyzed by LC/MS or transferred to a UV compatible 96-well plate and analyzed immediately in a UV/Vis spectrophotometer. At the end of the incubation time, the integrity of the artificial membrane layer can be measured using electrical resistance. The permeability method can also be used to determine the effect of pH on compound permeability by adjusting the pH of the solutions used in the analysis. These plates are particularly recommended for use in pre-ADME or Discovery programs requiring compound rank ordering or profiling.

application note



MILLIPORE



General Protocol Considerations

In addition to the MultiScreen Permeability plate, a UV/Vis spectrophotometer capable of analyzing 96-well plates and UV compatible, 96-well sample plates are also required to run the assay. Preparation of the hexadecane layer in a fume hood is recommended.

Experimental Materials

Methotrexate (cat. A-7019), propranolol (cat. P-0884), warfarin (cat. A-2250), carbamazepine (cat. C-8981), furosemide (cat. F-4381), testosterone (cat. T-1500), hexadecane (cat. H670-3), hexane (cat. 27050-4), dimethyl sulfoxide (cat. D-8779) and phosphate buffered saline (cat. P-3813) were purchased from Sigma Chemical Co. (St. Louis, MO). MultiScreen Permability assay plates (cat. MAPBMN310) and PTFE Acceptor plate (cat. MSSACCEPTOR) are available from Millipore Corporation (Bedford, MA). Spectramax® Plus microtiter plate reader, SoftMax[®] Pro and UV compatible quartz plate (part no. R8024) were purchased from Molecular Devices (Sunnyvale, CA). Polypropylene reagent reservoirs (cat. 175-RBAS-000) were purchased from ELKay laboratory consumables (Shrewsbury, MA). Finnpipette® electronic pipettor (cat. 21377232) was purchased from Thermolab Systems (Helsinki, Finland). Biohit™ 8 channel electronic pipettor and polypropylene tips (cat. W67-710-800 and W16-160045) were purchased from Vangard International (Neptune, NJ). An ohm meter and 96 well Trans Epithelial Electrical Resistance (TEER) testing tray (model #'s EVOMX-G and MULTI-96) were purchased from World Precision Instruments (Sarasota, FL).

Electrical Resistance Testing

To ensure that hexadecane layers were intact, electrical resistance measurements were made both before and after permeability assays were conducted. Intact hexadecane layers exhibit extremely high electrical resistance (normally exceeding 25 k Ω). Data from wells with electrical resistance measurements below 5 k Ω were excluded.

Methods

The following protocol was used to determine log P_e for methotrexate, propranolol, warfarin, carbamazepine, furosemide and testosterone.

- a) Prepare a 5% solution (v/v) of hexadecane in hexane (~3 mL/plate).
- b) Pipette 15 µL of the hexadecane/hexane mixture each well carefully, avoiding pipette tip contact with the membrane. Note: use polypropylene reservoir.
- c) Allow the plates to dry for 1 hour in a fume hood to ensure complete evaporation of the hexane resulting in a uniform layer of hexadecane.
- d) Add 300 µL of buffer (5% DMSO in phosphate buffered saline pH 7.4) to each well of the PTFE Acceptor plate (MSSACCEPTOR).
- e) Place the hexadecane treated Donor plate into the Acceptor plate making sure the underside of the membrane is in contact with the buffer.
- f) Dissolve drugs of interest in 5% DMSO, PBS to the desired concentration and add 150 µL to each well in the Donor plate (for the following experiments testosterone Donor concentration equals 100 µM, all other drugs at 500 µM).

- g) Replace plate lid and incubate at room temperature for 5 hours.
- h) After incubation, measure UV/Vis absorption from 250 to 500 nm for 100 μL/well of the Donor solution and 250 μL/well of the Acceptor solution.
- Make up drug solutions at the theoretical equilibrium (i.e., the resulting concentration if the Donor and Acceptor solutions were simply combined) and measure UV/Vis absorption from 250 to 500 nm for 250 µL/well of each.
- j) Calculate log P_e using the equations provided (see below).

Drug Standard Curves

Standard curves were prepared to determine limits of quantification for each drug when using UV/Vis spectroscopy as the detection method. The limit of quantification is defined as the concentration of compound whose absorbance (optical density, OD) is five times areater than the standard deviation of the background (noise) absorbance. Two-fold serial dilutions starting at 100 µM of each drug were made in 5% DMSO/PBS. UV/Vis absorbance for each dilution was measured across 250 - 500 nm in 10 µm steps using the Spectramax Plus plate reader and quartz sample plate. Peak maximum and area under the curve data were then plotted for each drug (see Table 1 and Figures 1a - e).

Results of Log Pe Studies

To assess the analytical variability (precision) within a plate, the permeability of testosterone and propranolol were measured in all 96 wells of multiple plates (n = 1 plate for testosterone, n = 3 plates for propranolol). Each value in Table 2 is the average of the log P_e calculated from 96 readings (wells).

Equations

Log Pe can be calculated from the following equation as reported by Faller et al.¹

 $\log P_e = \log \{C \bullet -1n(1 - \frac{[drug]_{Acceptor}}{[drug]_{equilibrium}})\} \text{ where } C = (\frac{V_D \bullet V_A}{(V_D + V_A) \text{ Area } \bullet \text{ time}})$

The five drugs listed in Table 3 were assayed for permeability using four different plates over two days to measure the reproducibility of assay performance from plate to plate. Each value in the table is the average of 16 wells per drug for each plate.

To determine the assay reproducibility from day to day, the six drugs in Table 4 were assayed for permeability on eight different days, in two different plates using 16 wells per plate. [Note: Testosterone data on days 1, 2 and 3 were not included because of problems with solubility. Beginning on day 4 the Donor concentration was dropped to 100 µM.]

Five different lots of membrane were tested to determine lot to lot reproducibility. All assays were performed using the standard protocol described. These results are listed in Table 5. [Note: Lot 2 was tested before the solubility issue with testosterone was resolved.]

Summarized in Table 6 are the results of several small protocol changes that were tested to determine the impact of protocol variability. The impact of changes in how the hexadecane layer is formed was a major focus of the tested changes. The average and standard deviation for the standard protocol includes data from at least 5 lots, 8 days and 2 plates with 16 wells for each drug per plate (i.e., at least 1280 assays). For each of the protocol changes, 1 plate was evaluated with 16 wells for each drug. Diffusion of the tested compounds from the Donor to Acceptor compartments was monitored over the course of 48 hours. Drugs were added to successive rows of a plate for a set time (48, 32, 28, 24, 8, 6, 4, 2 hours) prior to plate separation and analysis by UV/Vis. Each row of the plate had 2 replicates for each of the 6 drugs and two plates were assayed resulting in a total of 4 replicates of each drug per time point. The average of these replicates is plotted in Figures 2 and 3 for the Donor and Acceptor plates, respectively. The log Pe calculated for each time point is listed in Table 7.

Table 1: UV/Vis Limits of Quantification (LOQ)

Drug	LOQ Concentration	Typical Acceptor Concentration After Incubation
Propranolol	5.0 µM	82.5 μM
Warfarin	2.89 µM	12.3 µM
Methotrexate	1.36 µM	0.07 µM
Carbamazepine	4.51 µM	59.8 µM
Furosemide	1.6 µM	0.09 µM
Testosterone	0.91 µM	24.1 µM

Figure 1a: Propranolol Standard Curve




Figure 1b: Warfarin Standard Curve




Figure 1c: Methotrexate Standard Curve




Figure 1d: Carbamazepine Standard Curve




Figure 1e: Furosemide Standard Curve





Figure 1f: Testosterone Standard Curve



Discussion

The permeability of the 6 drugs was determined using a Donor plate concentration of $500 \ \mu$ M. For testosterone, the Donor plate concentration was reduced to $100 \ \mu$ M after initial testing revealed its marginal solubility at the higher concentration. Overall, there were no significant (< 0.2 log units) variations due to well-to-well,

Table 2: Well-to-well Reproducibility

plate-to-plate, day-to-day and lot-to-lot changes for the medium and high permeability drugs. Log P_e values obtained for furosemide and methotrexate were more variable than for the other four drugs because of analytical limitations in detecting the small amounts present in the Acceptor plate (using the standard, 5 to 7 hour incubation period). The results obtained here agree well with the ranking reported by Faller *et al.*¹ with the exception of warfarin. This variation is probably due to the difference in pH (7.4 versus 6.8) and the fact that warfarin contains an ionizable group with a pKa in this range.

In addition to testing assay reproducibility, several small protocol changes were tested to determine the effect that typical protocol variability might have on the calculated log P_e . A table of the tested protocol changes can be found in the results section (Table 6). By and large, the impact of any of these protocol modifications is quite small and does not result in any

Drug	$Log P_e \pm 1 S.D.$	$Log P_e \pm 1 S.D.$	$Log P_e \pm 1 S.D.$					
Testosterone	-3.7 ± 0.09							
Propranolol	-4.3 ± 0.07	-4.1 ± 0.10	-4.0 ± 0.09					

Table 3: Plate-to-plate Reproducibility

Day	Plate	Data	Drug Carbamazepine	Furosemide	Methotrexate	Warfarin	Propranolol
1	1	Ave.	-4.2 ± 0.00	-6.8 ± 0.48	-6.9 ± 0.52	-4.9 ± 0.03	-4.1 ± 0.06
	2	Ave.	-4.1 ± 0.03	-6.9 ± 0.62	-6.7 ± 0.56	-4.8 ± 0.00	-4.0 ± 0.05
2	1	Ave.	-4.3 ± 0.00	-6.6 ± 0.24	-6.7 ± 0.30	-5.1 ± 0.03	-4.1 ± 0.05
	2	Ave.	-4.3 ± 0.03	-6.5 ± 0.43	-7.0 ± 0.36	-5.1 ± 0.04	-4.0 ± 0.04

Table 4: Day-to-day Reproducibility

Day	Data	Carbamazepine	Furosemide	Methotrexate	Testosterone	Warfarin	Propranolol
1	Ave.	-4.1 ± 0.05	-6.3 ± 0.78	-6.3 ± 0.61	Sol. Issue	-4.8 ± 0.09	-4.0 ± 0.07
2	Ave.	-4.2 ± 0.05	-6.8 ± 0.55	-6.8 ± .054	Sol. Issue	-4.8 ± 0.05	-4.1 ± 0.06
3	Ave.	-4.3 ± 0.02	-6.5 ± 0.34	-6.9 ± 0.35	Sol. Issue	-5.1 ± 0.04	-4.0 ± 0.05
4	Ave.	-4.1 ± 0.00	-7.2 ± 0.63	-7.6 ± 0.56	-3.7 ± 0.06	-4.9 ± 0.00	-3.9 ± 0.05
5	Ave.	-4.1 ± 0.04	-6.8 ± 0.42	-7.1 ± 0.51	-3.7 ± 0.05	-4.9 ± 0.00	-4.0 ± 0.05
6	Ave.	-4.2 ± 0.05	-7.0 ± 0.37	-7.4 ± 0.42	-3.8 ± 0.07	-5.0 ± 0.03	-4.1 ± 0.05
7	Ave.	-4.2 ± 0.04	-7.0 ± 0.42	-7.3 ± 0.85	-3.8 ± 0.05	-5.0 ± 0.06	-4.0 ± 0.02
Tote	al Ave.	-4.2 ± 0.07	-6.8 ± 0.59	-7.0 ± 0.69	-3.7 ± 0.12	-4.9 ± 0.11	-4.0 ± 0.07

Table 5: Lot-to-lot Reproducibility

Lot	Data	Drug Carbamazepine	Furosemide	Methotrexate	Testosterone	Warfarin	Propranolol
1	Ave.	-4.2 ± 0.08	-6.7 ± 0.49	-6.8 ± 0.45	-3.7 ± 0.15	-5.0 ± 0.13	-4.1 ± 0.06
2	Ave.	-4.1 ± 0.05	-6.3 ± 0.78	-6.3 ± 0.61	Sol. Issue	-4.8 ± 0.09	-4.0 ± 0.09
3	Ave.	-4.1 ± 0.05	-7.1 ± 0.58	-7.4 ± 0.60	-3.7 ± 0.09	-4.9 ± 0.05	-4.0 ± 0.05
4	Ave.	-4.2 ± 0.08	-7.0 ± 0.41	-7.4 ± 0.69	-3.7 ± 0.10	-5.0 ± 0.06	-4.0 ± 0.07
5	Ave.	-4.2 ± 0.05	-6.9 ± 0.40	-7.2 ± 0.56	-3.7 ± 0.09	-4.9 ± 0.05	-4.0 ± 0.06
То	tal Ave.	-4.2 ± 0.07	-6.8 ± 0.59	-7.0 ± 0.69	-3.7 ± 0.12	-4.9 ± 0.11	-4.0 ± 0.09

change in the rank ordering of compound permeability.

As described, the typical orientation of the assay is top down (starting with drug on the top), although the assay should perform comparably if the format is inverted (i.e., by adding the compounds to the bottom and analyzing the solution on top of the artificial membrane). Included in Table 6 are the results of performing the assay in the inverted orientation. The values for high-permeable compounds are reduced while low permeable compounds are unaffected and this is because the unstirred water layer is thicker in the inverted orientation. It is not recommended to use this orienta-

Table 6: Tested Protocol Variations

tion because it reduces the dynamic range of the assay. In general the rank order was unaffected.

A 48-hour time course experiment was performed to monitor the concentration changes of the Donor and Acceptor compartments for all six drugs. Graphs of the OD readings for each drug can be seen in the results section (see figures 2 & 3). As can be seen in the data, permeability rates are not constant and are governed by many factors. Short incubation times (e.g., 2 hours) result in observed permeability rates that are elevated relative to the standard protocol (5 hours). This observation may be due to concentration gradients driving the initial rate. At the end of the 48 hour incubation some drugs are approaching equilibrium, indicated by the flattening out of the curve.

As with any assay, there is an initial period of learning and trial and error before reproducible results are obtained. Experience and familiarity with an assay helps in spotting potential problems and making adjustments accordingly. Some potential pitfalls are not obvious, such as a drug's ability to be assayed by UV/Vis spectroscopy. Many drugs are not good UV/Vis chromophores so the proper analytical techniques must be chosen accordingly. Refer to Table 1 for limits of detection determined from the stan-

Assay Changes	Propranolol	Methotrexate	Warfarin	Carbamazepine	Furosemide	Testosterone
Standard Protocol	-4.1 ± 0.02	-7.3 ± 0.53	-5.0 ± 0.02	-4.3 ± 0.02	-7.0 ± 0.32	-3.8 ± 0.05
Inverted assay	-4.5 ± 0.03	-6.8 ± 0.34	-5.4 ± 0.04	-4.7 ± 0.02	-6.1 ± 0.19	-4.4 ± 0.03
Hex. side of well	-4.1 ± 0.03	-7.2 ± 0.59	-5.2 ± 0.03	-4.3 ± 0.03	-6.4 ± 0.16	-3.8 ± 0.02
4 °C incubation	-4.6 ± 0.04	-7.1 ± 0.61	-5.6 ± 0.14	-4.9 ± 0.04	-6.1 ± 0.22	-4.2 ± 0.04
37 °C incubation	-3.9 ± 0.02	-7.3 ± 0.36	-4.5 ± 0.12	-4.0 ± 0.05	-6.2 ± 0.06	-3.5 ± 0.10
2hr incubation	-4.0 ± 0.02	-6.7 ± 0.31	-5.2 ± 0.04	-4.3 ± 0.03	-6.0 ± 0.05	-3.7 ± 0.05
16hr incubation	-4.2 ± 0.03	-7.6 ± 0.51	-5.1 ± 0.01	-4.3 ± 0.02	-6.6 ± 0.26	-3.8 ± 0.10
10 µL Hex. mixture (5%)	-4.0 ± 0.03	-7.0 ± 0.30	-5.1 ± 0.01	-4.2 ± 0.02	-6.1 ± 0.23	-3.8 ± 0.08
20 µL Hex. mixture (5%)	-4.0 ± 0.01	-7.7 ± 0.40	-5.1 ± 0.02	-4.2 ± 0.01	-6.2 ± 0.11	-3.7 ± 0.02
3% Hex. (15 µL applied)	-3.8 ± 0.04	-4.9 ± 0.61	-4.5 ± 0.25	-4.0 ± 0.11	-4.4 ± 0.50	-3.6 ± 0.05
7% Hex. (15 µL applied)	-3.9 ± 0.02	-7.2 ± 0.46	-5.0 ± 0.01	-4.2 ± 0.01	-6.1 ± 0.08	-3.6 ± 0.06

Table 7: Average log Pe for Different Incubation Times

Time (hrs)	Drug Carbamazepine	Furosemide	Methotrexate	Propranolol	Testosterone	Warfarin
2	-4.2	-5.5	-5.3	-3.9	-3.7	-4.9
4	-4.2	-6.0	-6.3	-4.0	-3.7	-5.0
6	-4.2	-6.1	-6.2	-4.0	-3.7	-5.0
8	-4.3	-6.2	-6.3	-4.1	-3.8	-5.0
24	-4.3	-6.7	-7.0	-4.3	-3.8	-5.0
28	-4.3	-6.9	-7.2	-4.3	-4.0	-5.0
32	-4.3	-6.8	-6.9	-4.4	-4.0	-5.1
48	-4.4	-7.1	-5.5	-4.3	-4.2	-4.9

dard curves of the six drugs assayed. Carbamazepine, for example, exhibits a comparatively high limit of detection by UV/Vis, however due to its relatively high permeability; sufficient quantities are present in the Acceptor compartment to permit the use of this detection method. Methotrexate, on the other hand, is a good UV/Vis chromophore and exhibits a much lower limit of detection, however due to its slow permeability rate, very little compound is present in the Acceptor volume to be measured. Consequently, greater variability in the calculated log P_{e} of methotrexate is observed.

Another significant source of method error may be related to compound solubility. Initial attempts to measure testosterone permeability were conducted using a Donor solution concentration of 500 μ M – a concentration at which testosterone is only marginally soluble. This resulted in the drug partitioning onto the plate and membrane surfaces – thereby rendering accurate log P_e calculations impossible.

Lastly, environmental factors such as the temperature in the lab can affect the reproducibility of permeability experiments. The permeability assay

Figure 2: Ratio of OD Donor/Initial Donor Concentration





Figure 3: Ratio of OD Acceptor/Equilibrium OD



protocol calls for room temperature incubations. If the ambient temperature in the lab tends to fluctuate it could affect the measured permeability rates of some drugs. The melting temperature of hexadecane is 18 °C and even small temperature changes have the potential to affect the viscosity/permeability of the artificial membrane layer.

Reference

 (a) Faller, B.; Wohnsland, F. "Physicochemical Parameters as Tools in Drug Discovery and Lead Optimization." in: Testa, B., van de Waterbeemd, H., Folkers, G., Guy, R. (Eds.), Pharmacokinetic Optimization in Drug Research, Verlag Helvetica Chimica Acta: Zürich and Wiley – VCH: Weinheim, (2001) pp. 257 – 274.

(b) Wohnsland, F.; Faller, B. "High-throughput Permability pH profile and High-throughput alkane/water log P with artificial membranes." *J. Med. Chem.* (2001), 44, pp. 923 – 930.

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office. In the U.S. and Canada, call toll-free **1-800-MILLIPORE** (1-800-645-5476) In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439) Internet: www.millipore.com

Tech Service: www.millipore.com/techservice

pure e-commerce

Now you can buy Millipore products online @



www.millipore.com/purecommerce

Millipore Worldwide

AUSTRALIA Tel. 1 800 222 111 or (02) 9888 8999 Fax (02) 9878 0788 AUSTRIA Tel. (01) 877 89 26 Fax (01) 877 16 54 BALTIC COUNTRIES Tel. +358 9 804 5110 Fax +358 9 256 5660 BELGIUM AND LUXEMBOURG Tel. +32 2 726 88 40 Fax +32 2 726 98 84 BRA7II Tel. (011) 5548-7011 Fax (011) 5548-7923 CANADA Tel. 1-800-645-5476 Fax 1-800-645-5439 CHINA, PEOPLE'S REPUBLIC OF Beiiina Tel. (86-10) 8519 1250 (86-10) 6518 1058 Fax (86-10) 8519 1255 Guangzhou: Tel. (86-20) 8752 0187 (86-20) 8752 0173 Fax (86-20) 8752 0172 Hong Kong: Tel. (852) 2803 9111 Fax (852) 2513 0313 Shanghai: Tel. (86-21) 5306 9100 Fax (86-21) 5306 0838

CZECH REPUBLIC Tel. 02-2051 3841 Fax 02-2051 4298 DENMARK IRELAND

ITALY

Roma:

JAPAN

KOREA

MALAYSIA

MEXICO

Life Sciences and

Laboratory Water Divisions:

Tel. +44 1923 816375

Fax +44 1923 818297

Tel. (021) 883 666

Fax (021) 883 048

Vimodrone (Milano):

Tel. (02) 25.07.81

Fax (02) 26.50.324

Tel. (06) 52.03.600

Fax (06) 52.95.735

Tel. (03) 5442-9711

9736 Life Sciences

9737 BioPharm.

9734 Lab Water

Tel. (822) 551-0230

Fax (822) 551-0228

Tel. 03-7957-1322

Fax 03-7957-1711

Tel. (55) 5576 9688

Fax (55) 5576 8706

THE NETHERLANDS

Tel. 0900 7645645

Fax 0900 7645644

Fax (03) 5442-

BioPharmaceutical Division:

Tel. 70 10 00 23 Fax 70 10 13 14 EASTERN EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF

Life Sciences Division: Tel. +33 3 88 38 9536 Fax +33 3 88 38 9539 BioPharmaceutical Division: Tel. +43 1 877-8926 Fax +43 1 877-1654 Laboratory Water Division: Tel. +33 1 30 12 7000 Fax +33 1 30 12 7180 FINLAND

Tel. (09) 804 5110 Fax (09) 256 5660 FRANCE

Tel. 01 30 12 7000 Fax 01 30 12 7180 GERMANY

```
Tel. (06196) 494-0
Fax (06196) 43901
```

HUNGARY Tel. 01 381 0433

```
01 381 0434
01 209 3232
Fax 01 209 0295
```

INDIA Tel. (91) 80-839 46 57 Fax (91) 80-839 63 45

NORWAY Tel. 22 67 82 53 Fax 22 66 04 60

POLAND

Tel. 22-669 12 25 22-663 70 31 Fax 22-663 70 33

PUERTO RICO Tel. (787) 273-8495 Fax (787) 747-6553

SINGAPORE Tel. 6842-1822 Fax 6842-4988

SPAIN AND PORTUGAL Tel. +34 917 283 960

Fax +34 917 292 909 SWEDEN

Tel. 08-628 6960 Fax 08-628 6457

SWITZERLAND Tel. (01) 908-30-60

Fax (01) 908-30-80

TAIWAN

Tel. 886-2-2792-9333 Fax 886-2-2792-6555

U.K. Tel. 01923 816375 Fax 01923 818297

U.S.A. Tel. 1-800-645-5476

Fax 1-800-645-5439

COUNTRIES

Millipore Intertech (U.S.A.) Tel. +1 (781) 533-8622 Fax +1 (781) 533-8630

Millipore is a registered trademark of Millipore Corporation. MultiScreen is a trademark of Millipore Corporation. Spectramax and SoftMax are registered trademarks of Molecular Devices Corporation. Finnpipette is a registered trademark of Thermolab Systems Oy. Biohit is a trademark of Biohit Oy. Lit. No. AN1725EN00 Rev. A 12/02 Printed in the U.S.A. 02-305 © 2002 Millipore Corporation, Bedford, MA 01730 U.S.A. All rights reserved.

MILLIPORE