

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

#### Xanthine

BioUltra

Catalog Number X4002

CAS RN 69-89-6

Synonyms: 2,6-Dihydroxypurine Molecular Formula: C₅H₄N₄O₂ Molecular Weight: 152.1

pK<sub>a</sub>:1 7.7, 11.9

 $\lambda_{max}$ :1 277 nm (pH = 10)

Extinction Coefficient:  $E_{mM} = 9.3 \text{ (pH = 10)}$ 

### **Product Description**

Xanthine is a catabolic product of purine nucleotide metabolism, produced from several different precursors in the purine metabolic pathway, such as:

- Deamination of guanine by guanine deaminase
- Conversion of hypoxanthine by xanthine oxidoreductase

Xanthine oxidase converts xanthine to uric acid.

Natural sources of xanthine or its derivatives include animal organs, yeast, potatoes, coffee beans, and tea.<sup>2</sup> Xanthine and xanthine oxidase are used to generate superoxide radicals used to measure the activity of superoxide dismutase.<sup>3</sup>

Various publications have cited use of this product in different systems, including:

- Yeast culture<sup>4</sup>
- Human cell line culture<sup>5</sup>
- Recombinant viruses cloned as bacterial artificial chromosomes in E. colf<sup>6</sup>

Trace elemental analyses have been performed on the BioUltra xanthine. The Certificate of Analysis (CofA) provides lot-specific results. BioUltra xanthine is for applications which require tight control of elemental content.

#### **Precautions and Disclaimer**

For R&D use only. Not for drug, household, or other uses.

#### **Preparation Instructions**

Xanthine is soluble in sodium hydroxide solutions and in acidic solutions. It is soluble in 1 M NaOH (50 mg/mL), with sonication for less than 5 minutes, yielding a clear solution. It is slightly soluble in water (1 g/14.5 L, 16 °C) and in ethanol.

## Storage/Stability

The decomposition of solutions of xanthine in either  $0.5 \text{ M H}_2\text{SO}_4$  or in 10 M NaOH is <10% after one hour at 100 °C.<sup>7</sup> Stock solutions of at least 10 mM xanthine in NaOH can be stored at 2-8 °C for one week.<sup>8</sup>

#### References

- Specifications and Criteria for Biochemical Compounds, Third Edition. National Academy of Sciences (Washington, DC), p. 182 (1972).
- 2. The Merck Index, 13th ed., Entry# 10116.
- 3. McCord, J.M., and Fridovich, I., *J. Biol. Chem.*, **244(22)**, 6049-6055 (1969).
- Castoria, R. et al., Phytopathology, 93(5), 564-572 (2003).
- 5. Pillay, N. et al., Cancer Cell, **35(3)**, 519-533 e8 (2019).
- 6. Borst, E.-M. et al., Methods Mol. Biol., **256**, 221-239 (2004).
- 7. Data for Biochemical Research, 3rd ed. (Dawson, R.M.C. et al., eds.). Oxford University Press (Oxford, UK), pp. 94-95 (1986).
- Heinz, F., and Reckel, S., in *Methods Of Enzymatic Analysis*, Vol. III, 3rd ed. (Bergmeyer, H.U., ed.). Verlag-Chemie (Weinheim, Germany), pp. 210-216 (1983).

ARO, RXR, GCY 04/20-1