3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

SIGMA QUALITY CONTROL TEST PROCEDURE

Enzymatic Assay of ADENOSINE DEAMINASE (EC 3.5.4.4)

PRINCIPLE:

Adenosine + H₂O Adenosine Deaminase > Inosine + NH₃

CONDITIONS: T = 25°C, pH = 7.5, A_{265nm} , Light path = 1 cm

METHOD: Continuous Spectrophotometric Rate Determination

REAGENTS:

A. 100 mM Potassium Phosphate Buffer, pH 7.5 at 25°C (Prepare 150 ml in deionized water using Potassium Phosphate, Monobasic, Anhydrous, Sigma Prod. No. P-5379. Adjust to pH 7.5 at 25°C with 1 M KOH.)

- 1.35 mM Adenosine Solution, pH 7.0 at 25°C¹ B. (Prepare 10 ml in deionized water using Adenosine, Sigma Prod. No. A-9251. Adjust to pH 7.0 at 25°C with 10 mM NaOH.)
- C. 0.1% (w/v) Bovine Serum Albumin Solution (BSA) (Prepare 50 ml in Reagent A using Albumin, Bovine, Sigma Prod. No. A-4503.)
- D. Adenosine Deaminase Enzyme Solution (Immediately before use, prepare a solution containing 0.20 - 0.40 unit/ml of Adenosine Deaminase in cold Reagent C.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable quartz cuvettes:

	<u>Test</u>	Blank
Deionized Water	1.30	1.30
Reagent A (Buffer)	1.50	1.50
Reagent B (Adenosine)	0.10	0.10

SPADEN01

Enzymatic Assay of ADENOSINE DEAMINASE (EC 3.5.4.4)

PROCEDURE: (continued)

Mix by inversion and equilibrate to 25°C. Monitor the A_{265nm} until constant, using a suitably thermostatted spectrophotometer. Then add:

	<u>Test</u>	<u>Blank</u>
Reagent C (BSA)		0.10
Reagent D (Enzyme Solution)	0.10	

Immediately mix by inversion and record the decrease in the A_{265nm} for approximately 5 minutes. Obtain the ΔA_{265nm} /minute using the maximum linear rate for both the Test and Blank.

CALCULATIONS:

Units/ml enzyme =
$$\frac{(\Delta A_{265nm}/min \text{ Test - } \Delta A_{265nm}/min \text{ Blank})(3)(df)}{(8.1) (0.1)}$$

3 = Total volume (in milliliters) of assay

df = Dilution Factor

8.1 = Millimolar extinction coefficient of adenosine at 265 nm

0.1 = Volume (in milliliters) of enzyme used

UNIT DEFINITION:

One unit will deaminate 1.0 µmole of adenosine to inosine per minute at pH 7.5 at 25°C.

FINAL ASSAY CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 53.3 mM potassium phosphate, 0.045 mM adenosine, 0.003% (w/v) bovine serum albumin and 0.02 - 0.04 unit of adenosine deaminase.

Revised: 08/05/95

SPADEN01

Enzymatic Assay of ADENOSINE DEAMINASE (EC 3.5.4.4)

REFERENCE:

Murphy, J., Baker, D.C., Behling, C. and Turner, R.A. (1982) Anal. Biochemistry 122, 328-337

Bergmeyer, H.U. (1983) in *Methods of Enzymatic Analysis*, 3rd ed., Vol. 2, pp 135-136, Verlag Chemie, Deerfield Beach, FL

NOTES:

- 1. Substrate concentration should not be significantly increased above that listed in the procedure. Higher concentrations result in a deviation from Beers Law leading to errors in the rate determinations according to Murphy et al. (1982).
- 2. This assay is based on the cited references.
- 3. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.

SPADEN01 Page 3 of 3 Revised: 08/05/95