

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

Pyronin Y

Product Number **P 9172** Store at Room Temperature

Product Description

Molecular Formula: C₁₇H₁₉CIN₂O

Molecular Weight: 302.8 CAS Number: 92-32-0 Melting Point: 250-260 °C λ_{max} : 548 nm (50% ETOH)

Synonyms: Pyronine Y, Pyronin J, and Pyronin G,

and cationic xanthene dye

Pyronin Y solutions are red in transmitted light, displaying a yellow flouresence in reflected light. Pyronin Y is used widely in combination with methyl green to selectively and differentially stain nucleic acids. In the dichromatic dye solution, Pyronin Y stains RNA red, while methyl green stains the DNA green. ¹ Transfer RNA (tRNA) will complex with Pyronin Y at a 1:10 molar ratio when analyzed by polyacrylamide gel electrophoresis.² tRNA

(4 to 8 μ g) complexed with Pyronin Y is visible as a red band during electrophoresis. The stained tRNA is also detectable by UV fluoresence at 302 nm. Pyronin Y can also be used in correlating RNA content by flow cytometry. In addition, Pyronin Y solutions may be utilized in staining liver cells that have been fixed in Carnoy's fluid. 5

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water 1 mg/ml), yielding a clear dark red solution.

Procedure

Suitability Test for RNA Staining:

- Prepare a tRNA (Product No. R 4251) solution at 1 μg/μl.
- Prepare Pyronin Y solution: Dissolve 25 mg in 2.5 ml of 15% acetic acid in water. Once fully dissolved dilute the solution 5-fold with water.
- 3. Mix 220 μl of tRNA solution (step 1) with 2.5 μl of Pyronin Y solution (step 2). Equilibrate at 22 °C for 1.5 hours. To a 50 μl aliquot, add 12.5 μl of gel loading solution (G 2526). Load 10 μl of the tRNA/gel loading solution into the wells of a 20% TBE polyacrylaminde gel. Run until the blue dye front is about 1 cm from the bottom of the gel. Carefully remove the gel from the gel holder and photograph.

References

- 1. The Sigma-Aldrich Handbook of Stains, Dyes & Indicators, Green, F. J., ed., Aldrich Chemical Co. (Milwaukee, WI: 1990), pp. 601-602.
- Hassur, S. M., and Whitlock, H. W. Jr., A gel electrophoresis method for determining the relative binding constants of biologically active, intercalating fluorescent stains. Anal. Biochem., 95(2), 329-339 (1979).
- Crissman, H. A., et al., Normal and perturbed Chinese hamster ovary cells: correlation of DNA, RNA, and protein content by flow cytometry. J. Cell. Biol., 101(1), 141-147 (1985).
- 4. Clinical Flow Cytometry: Principles and Applications, Bauer, K. D., et al., Williams and Wilkins, (Philadelphia, PA: 1993), p. 377.
- 5. Conn's Biological Stains, 9th ed., Lillie, R. D., Williams and Wilkins (Baltimore, MD: 1977), p. 590.

TMG/RXR 10/03