Instructions for Use

Reticulum Stain Kit

Procedure No. HT102

Intended Use

The Reticulum Stain Kit is intended for use as a general laboratory histochemical stain to demonstrate reticular fibers. Reticulum Stain Kit reagents are for "In Vitro Diagnostic Use." For professional use only. The data obtained from this manual, qualitative procedure demonstrates reticular fibers in tissue samples of human specimens. Histological visualization of reticular fibers is a unique technique currently commonly used in medicine which may be used as an aid to diagnosis of human disease, particularly some liver diseases.

The main function of reticular fibers is to provide support. They are normally found throughout the body, particularly in liver, lymph node, spleen and kidney.¹ Ammoniacal silver stains are the most commonly used methods for demonstration of reticular fibers. In the procedure of Gordon and Sweets, tissue sections are oxidized by potassium permanganate with oxalic acid removing the excess potassium permanganate. Ferric ammonium sulfate acts as the sensitizer. After the silver impregnation, formalin is used to reduce the silver to its visible metallic form. Gold chloride tones the sections and any unreduced silver is removed by sodium thiosulfate. A counterstain may be used, if desired.²

Reagents

Sodium Hydroxide Solution (Cat. No. HT1021-100ML)
Sodium Hydroxide, 3% aqueous solution

Potassium Permanganate Solution (Cat. No. HT1022-100ML)
Potassium Permanganate, 1% aqueous solution

Oxalic Acid Solution (Cat. No. HT1023-100ML)
Oxalic Acid Dihydrate, 1% aqueous solution

Ferric Ammonium Sulfate Solution (Cat. No. HT1024-500ML)
Ferric Ammonium Sulfate, 2.5% aqueous solution

Silver Nitrate Solution (Cat. No. HT1025-50ML)
Silver Nitrate, 10% aqueous solution. Warning. Causes skin irritation. Causes serious eye irritation. Very toxic to aquatic life. Avoid release to the environment.

Gold Chloride Solution (Cat. No. HT1026-100ML)
Gold Chloride Trihydrate, 0.2% aqueous solution

Sodium Thiosulfate Solution (Cat. No. HT1027-500ML)
Sodium Thiosulfate, 5% aqueous solution

Special Materials Required but Not Provided

- Positive control slides, such as Sigma Reticulum TISSUE-TROL™ slides (Cat. No. TTR010) should be included in each run
- Counterstain (optional) Eosin Y solution or Nuclear Fast Red Solution (Cat. No. N3020), nuclear fast red 0.1% in 5% aluminum sulfate
- Ammonium Hydroxide Solution, concentrated
- Ethanol, absolute
- Formalin Solution, 10%
- Reagent Alcohol
- Xylene or xylene substitute
- Forceps, plastic or paraffin coated
- Coplin jars, chemically clean

Storage and Stability

Store unopened Reticulum Stain Kit in refrigerator (2–8°C). After opening, Sodium Hydroxide, Potassium Permanganate, Oxalic Acid, Ferric Ammonium Sulfate and Sodium Thiosulfate Solutions may be stored in the refrigerator or at room temperature (2–26°C). Store Silver Nitrate Solution and Gold Chloride Solution in refrigerator (2–8°C). Reagents are stable until the expiration dates shown on the labels.

Silver Nitrate Working Solution should be used once and then discarded.

Nuclear Fast Red Solution (Cat. No. N3020) nuclear fast red 0.1% in 5% aluminum sulfate, should have limited exposure to air. Material oxidizes when exposed to air. Repeat exposure can cause the solution to become more alkaline over time, causing aluminum sulfate to fall out of solution.

In dry form ammoniacal salts may present an explosive hazard. Do not allow ammoniacal silver solution to dry out. Store used ammoniacal silver solution in plastic bottles – do not store in glass. Wallington recommended the inactivation of the ammoniacal silver solution by the addition of dilute hydrochloric acid or sodium chloride solution.³

Deterioration

Potassium Permanganate Solution should be purple. Solution may be used several times, but should be discarded if it turns brown.

Preparation

To prepare Ammoniacal Silver Nitrate Solution:

- Pipet 5 mL Silver Nitrate Solution in an Erlenmeyer flask.
- In a hood, while shaking or swirling the flask continuously, add concentrated ammonium hydroxide, drop by drop, until the precipitate it formed is completely dissolved. Do not add

excess ammonium hydroxide.

- Add 5 mL Sodium Hydroxide Solution to the flask. Solution will turn black and precipitate will form. Continuously swirl the flask and add concentrated ammonium hydroxide, drop by drop, until the precipitate just dissolves. At this stage the solution should not be completely clear. NOTE: If no cloudiness remains, add Silver Nitrate Solution drop by drop, until one drop causes the solution to become permanently cloudy. Only a faint cloudiness is desirable.
- Dilute the resulting solution to 50 mL with distilled or deionized water. Filter into a chemically clean Coplin jar. Use once and then discard.

10% Formalin Solution is prepared by diluting 5 mL formaldehyde, 37–40%, with 45 mL distilled or deionized water. Prepare fresh daily.

The Potassium Permanganate Solution and the Oxalic Acid Solution should be divided into two 50-mL aliquots. Each aliquot can be used 5 times and then discarded. Do not mix used reagent with unused solution in original bottle.

Other reagents are supplied ready for use.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Reticulum TISSUE-TROL™ control slides are paraffin embedded human tissue containing reticulum and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Fix specimen in 10% neutral buffered formalin, process and embed in paraffin. Cut paraffin sections at 4 to 5 microns. Incorporate appropriate controls.

Procedure

- Deparaffinize sections and hydrate to distilled water.
- Oxidize sections in Potassium Permanganate Solution for **5 minutes**.
- Rinse slides in tap water for **2 minutes**.
- Bleach in Oxalic Acid Solution for **2 minutes** or until sections are colorless
- Wash slides in tap water for **2 minutes**.
- Sensitize sections in Ferric Ammonium Sulfate Solution for **15 minutes**.
- Wash slides in several changes of distilled water.
- Impregnate sections with Ammoniacal Silver Nitrate Solution for **2 minutes**.
- Rinse slides well with distilled water.
- Reduce sections for **2 minutes** in 10% Formalin Solution.
- Wash slides in tap water for **3 minutes**.
- Tone sections in Gold Chloride Solution for **10 minutes**.
- Rinse slides in distilled water.
- Place slides in Sodium Thiosulfate Solution for **1 minute**.
- Wash slides in tap water for **2 minutes**.
- Counterstain, if desired, with Nuclear Fast Red Solution for **3–5 minutes** or an Eosin Y Solution for **1–2 minutes**. Generally, all sections except those from liver are counterstained. Wash well in water.
- Dehydrate in two changes each of 95% ethanol and absolute ethanol.
- Clear in xylene and mount with synthetic resin.

Performance Characteristics

Target Structure	Staining Result
Reticulum	Black
Background	Pink to rose (if counterstained with Nuclear Fast Red Solution)

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT1021	Sodium Hydroxide Solution	Reticulum	3 of 3	3 of 3	3 of 3	3 of 3
HT1022	Potassium Permanganate Solution	Reticulum	3 of 3	3 of 3	3 of 3	3 of 3
HT1023	Oxalic Acid Solution	Reticulum	3 of 3	3 of 3	3 of 3	3 of 3
HT1024	Ferric Ammonium Sulfate Solution	Reticulum	3 of 3	3 of 3	3 of 3	3 of 3
HT1025	Silver Nitrate Solution	Reticulum	3 of 3	3 of 3	3 of 3	3 of 3
HT1026	Gold Chloride Solution	Reticulum	3 of 3	3 of 3	3 of 3	3 of 3
HT1027	Sodium Thiosulfate Solution	Reticulum	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT102A:

H290: May be corrosive to metals.

H314: Causes severe skin burns and eye damage.

H410: Very toxic to aquatic life with long lasting effects.

P234: Keep only in original packaging.

P273: Avoid release to the environment.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P310: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed, CV Mosby Co., St. Louis, MO, 1980, pp 181–182
2. Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 144-155
3. Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 172

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 4.0	2020
Rev. 5.0	2022
Rev. 6.0	2022
	Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.
Rev. 7.0	2025
	Updated naming conventions. Revised Inteded Use statement. Added CH-REP and importer information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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MDSS GmbH
MDSS CH GmbH
Chemie GmbH