

## Product Information

# Anti-Mouse IgG (Whole Molecule)-Peroxidase

Produced in Goat, IgG Fraction of Antiserum

**A5278**

## Product Description

Anti-Mouse IgG is developed in goat using purified mouse IgG as the immunogen. Whole antiserum is fractionated and then further purified by ion exchange chromatography to provide the IgG fraction of antiserum. This fraction is essentially free of other goat serum proteins. Goat anti-mouse IgG is conjugated to Sigma Horseradish Peroxidase, Type VI (Cat. No. P8375) by a modification of the periodate method of Wilson and Nakane.<sup>1</sup>

Specificity of the goat anti-mouse IgG antibody is determined by Ouchterlony Double Diffusion (ODD) prior to conjugation. The antibody preparation reacts with mouse IgA, IgG1, IgG2a, IgG2b, IgG3, and IgM, myeloma proteins.

Identity and purity of the antibody is established by immunoelectrophoresis (IEP), prior to conjugation. Electrophoresis of the product followed by diffusion against anti-goat IgG and the anti-goat whole serum results in single arcs of precipitation in the gamma region.

## Reagent

The product is provided as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% BSA with preservative.

## Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Consult the SDS for information regarding hazards and safe handling practices.

## Storage/Stability

Store at -20 °C. The antibody may be stored at 2-8 °C for up to one month. For extended storage, the solution should be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

## Product Profile

We are now reporting lot specific information as a titer by direct ELISA (minimum 1:10,000) rather than as a working dilution. Titer is defined as the dilution of conjugate sufficient to give a change in absorbance of 1.0 at 450 nm after 30 minutes of substrate conversion at 25 °C. Microwell plates are coated with purified mouse IgG at a concentration of 5 µg/mL in 0.05 M carbonate/bicarbonate buffer, pH 9.6 (Carbonate/Bicarbonate Buffer Capsules are available as Cat. No. C3041).

Substrate: *o*-Phenylenediamine Dihydrochloride (OPD, Cat. No. P8287), 0.4 mg/mL in 0.05 M phosphate-citrate buffer, pH 5.0, containing 0.03% sodium perborate (Phosphate-Citrate Buffer Capsules with Sodium Perborate are available as Cat. No. P4922).

Working dilution should be determined by titration assay. Due to product improvement and changes in the assay procedure, we now list a lot specific titer by direct ELISA for this product. Due to differences in assay systems, this titer may not reflect the user's actual working dilution.

## References

1. Wilson, M. and Nakane, P., Immunofluorescence and Related Staining Techniques, Elsevier/North Holland BioMedical Press, Amsterdam, p. 215 (1978).
2. Voller, A., et al., Bull. World Health Organ., **53**: 55 (1976).

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.  
A5278dat Rev 05/21

