



# Regulatory Changes Affecting Pyrogen Testing with the Monocyte Activation Test



## Gaetan Bour

*R&D Scientist Expert  
Millipore SAS, Molsheim, France  
an affiliate of Merck, KGaA, Darmstadt, Germany*

Pyrogens represent a broad family of molecules that, when present in parenteral drugs, may induce adverse reactions in humans, ranging from fever to septic shock. This is why pyrogenicity testing became mandatory for injectable products, initially only with the Rabbit Pyrogen Test (RPT). However, the European Pharmacopoeia (Ph. Eur.) has now effectively banned the RPT and revised 57 monographs to remove the reference to it. A new general chapter on Pyrogenicity was introduced, and the texts on MAT were revised. The animal-based Limulus Amebocyte Lysate (LAL) test remains available, but it detects only Endotoxins (Etx), pyrogens from Gram negative bacteria, not Non-Endotoxin Pyrogens (NEP). The Ph. Eur. now promotes the Monocyte Activation Test (MAT), which uses human monocytic blood cells to mimic our immune reaction, as an alternative for the full spectrum of pyrogens. Launched in 2018, the PyroMAT® system was the first commercially available MAT solution based on a human monocytic cell line, offering, in a convenient format, ready-to-use monocytic cells, a ready-to-use IL-6 ELISA kit and validated analysis protocols for data processing.

## Endotoxin or full-spectrum pyrogen testing?

Chapter 5.1.13:2025 on pyrogenicity states that the choice of pyrogen test should depend on the possible sources of contamination. If the presence of NEP can be ruled out, then tests for bacterial endotoxins can be used; if not, the recommendation is to perform the MAT. To determine if NEP may be present, the MAT should be carried out alongside a Bacterial Endotoxin Test (BET) on the same product batches, as early as during development of the manufacturing process,

and again if manufacturing changes lead to the identification of new contamination risks. A risk assessment must be performed, taking into consideration all potential sources of pyrogen contamination. The risk can be influenced by many factors, including the production process, the bioburden of the API, the type and origin of the excipients or raw materials, the capability of the process to remove pyrogenic substances and the final use of the drug.

## The advantages of continuous cell lines

The new 2024 version of Ph. Eur. chapter 2.6.30 on the MAT now regards continuous monocytic cell lines, once qualified, to be suitable for both Etx and NEP detection. Continuous cell lines represent a virtually infinite source for large scale production and a long-term supply of cells with identical properties, offering standardization that the previously used whole blood or peripheral blood mononuclear cells could not achieve. The qualification tests require that cell performance is monitored, notably for intact receptor signaling in addition to growth criteria, signal and background noise. There is also clear guidance on qualitative criteria to demonstrate the stability of the cell line over the number of its passages.

For the Mono-Mac-6 (MM6) cells that our PyroMAT® system uses to detect endotoxin and non-endotoxin pyrogens, such tests are performed routinely for every batch, in particular to verify the absence of viral and microbial contamination, cell identity, cell performance (dose response curve, background level) as well as the activation of pattern recognition receptors by various TLR agonists targeting TLR 1/2/4/5/6/7/8.

## Changes to product-specific validation

Although revalidation of the MAT method is not required, its suitability for a specific substance in a specific environment must be determined. The new Ph. Eur. chapter 2.6.30:2024 now defines a series of criteria to verify that monocytic continuous cell lines are qualified to detect pyrogens. These tests must show that the criteria for the endotoxin curve are satisfied, that there is no interference between the test solution and the MAT cells for IL-6 expression nor with the IL-6 quantification by ELISA and, finally, that the test is able to detect non-pyrogen contaminants. The criteria revised in 2024 are summarized in the table.

Table 1.		
Topic	Criteria 2.6.30:2017	Criteria 2.6.30:2024
Standard curve: Limit of Detection (LOD)/ Test sensitivity	Value corresponding to the cut-off value	Renamed "Test Sensitivity", corresponding to the lowest point of the standard curve exceeding the cut-off value
Maximum Valid Dilution (MVD)	CLC x C/LOD	CLC x C/Test sensitivity
Standard curve: Blank OD value	OD < 0.1 in ELISA	OD as low as possible
Standard curve: Effect of dose	The regression of responses shall be statistically significant	Criteria removed
Standard curve: Coefficient of determination ( $r^2$ )	No criteria	Standard curve is valid only if $r^2 > 0.975$
Standard curve: Number of concentrations for regression	At least 4 concentrations	6 concentrations for a 5-parameter logistic model
Spiked NEP recovery	No criteria	Valid test: 50% < NEP recovery < 200% Invalid test: NEP recovery < 50% Synergism: NEP recovery > 200%
Dilution factors for the sample to test	2-fold dilution from the lowest dilution determined during PSV (not exceeding the MVD)	Geometric dilution from the lowest dilution determined during PSV (not exceeding the MVD)
Spiked Etx recovery	Valid test: 50% < Etx recovery < 200%	Valid test: at least one of the tested dilutions 50% < Etx recovery < 200%
CLC	1 valid dilution < CLC is the minimum	Failed test if mean pyrogen concentration > CLC for any concentration

To account for the changes in the Ph. Eur. on the MAT, the PyroMAT<sup>®</sup> analysis software was updated and fully revalidated according to good engineering practices so it complies with the new criteria and ensures conformity of the results with Ph. Eur. 2.6.30:2024.


Learn more about the PyroMAT<sup>®</sup> system for compliant detection of endotoxin and non-endotoxin pyrogens  
[SigmaAldrich.com/Sirius](https://SigmaAldrich.com/Sirius)

**MERCK**

## THE END OF RABBIT PYROGEN TEST

Effective as of July 2025 in the Ph. Eur.

Anticipate the Health Authority requirements & make the move to the Monocyte Activation Test (MAT)



[SigmaAldrich.com/pyromat](https://SigmaAldrich.com/pyromat)



The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Millipore and PyroMAT are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

05/2025

**Millipore<sup>®</sup>**

Preparation, Separation, Filtration & Monitoring Products