

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

ACETIC ACID, GLACIAL Sigma Prod. Nos. A0808 24,285-3 is an exact replacement for A0808

CAS NUMBER: 64-19-7

Boiling point: 118°C

SYNONYMS: acetic acid, methanecarboxylic acid, ethanoic acid, vinegar acid, ethylic acid¹

PHYSICAL DESCRIPTION: 1,2

Appearance: Clear colorless liquid Molecular formula: CH₃COOH Molecular weight: 60.05 Density: 1.049 g/mL Melting point: 16.7°C

pK_a = 4.75 at 25°C³; K_a has a maximum value of 1.754 x 10⁻⁵ at 25°C, but decreases in value at

temperatures above and below than 25°C. (1.700 x 10⁻⁵ at 5°C; 1.703 x 10⁻⁵ at 40°C³)

Glacial acetic acid has a very pungent odor, and is a corrosive liquid. It is also flammable, with a flash point of 40°C.¹

STORAGE / STABILITY AS SUPPLIED:

Acetic acid is extremely stable stored at room temperature. It should be stored tightly sealed to prevent vapors form escaping.

SOLUBILITY / STABILITY OF SOLUTIONS:

Acetic acid is infinitely miscible in water, also miscible with alcohol, glycerol, ether, carbon tetrachloride. It is practically insoluble in carbon disulfide.

Aqueous solutions are stable at room temperature, if well sealed.

GENERAL REMARKS:

By specification A0808 ("ACS Reagent") is ≥99.7%. In addition, A0808 is subjected to additional testing as required by the American Chemical Society.⁴

The name "glacial" is often used to emphasize the use of 100% acetic acid; frozen acetic acid looks like water ice.

ACETIC ACID, GLACIAL Sigma Prod. Nos. A0808

GENERAL REMARKS: (continued)

Acetic acid is widely used in manufacturing either as starting material or solvent; in the food industry, dilute solutions are used for preparation and preservation. Because acetic acid is a "weak acid", it is effective as buffer at pH 4.75 ± 1 . This buffer can be prepared by titrating acetic acid with hydroxide or by mixing equal volumes of equimolar acetic acid and sodium acetate solutions (standard mixing tables have been published).^{5,6}

METHODS OF PREPARATION:

Although acetic acid has historically been produced form bacterial oxidation of ethanol, on the industrial scale, it is usually obtained by destructive distillation of wood or oxidation of acetaldehyde.¹

REFERENCES:

- 1. Sigma Material Safety Data Sheet.
- 2. *Merck Index*, 12th ed., #52 (1996).
- 3. Handbook of Chemistry and Physics, 74th ed., D.R. Lide, ed. (CRC Press, 1993), p. 8-45 and 8-47.
- 4. Reagent Chemicals, 8th ed. (American Chemical society, 1993), pp. 98-101.
- 5. *Methods in Enzymology*, 182, 32 (1990).
- 6. Data for Biochemical Research, 3rd ed., Dawson et al., ed. (Oxford Press, 1987), p. 429.

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.