

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

## **ProductInformation**

# Pyruvate Kinase preparation from rabbit muscle

Product Number **P 1506** Storage Temperature 2-8 °C

**Product Description** 

Molecular Weight: 237 kDa<sup>1</sup> CAS Number: 9001-59-6

Enzyme Commission (EC) Number: 2.7.1.40

Molecular Weight: 237 kDa<sup>1</sup> Extinction Coefficient: E<sup>0.1%</sup> = 0.54<sup>1</sup>

Pyruvate kinase from rabbit muscle is a tetramer consisting of four equal subunits of molecular weight 57 kDa.<sup>1</sup>

Pyruvate kinase catalyzes the following reaction:

ATP + Pyruvate → ADP + Phophoenolpuruvate

Reported K<sub>M</sub> values are ATP (0.86 mM), pyruvate (10 mM), ADP (0.3 mM), and PEP (0.07 mM).<sup>2</sup>
Pyruvate kinase can also utilize other dinucleotide phosphates as substrates including GDP, IDP, dADP, UDP, CDP, dCDP.<sup>1</sup> The pH optimum for pyruvate kinase is 7.5<sup>3</sup> and the mechanism of the reaction catalyzed by pyruvate kinase has been described in the literature.<sup>4</sup>

Both Mg<sup>2+</sup> and K<sup>+</sup> are required metal cofactors for optimal activity and the enzyme is inhibited by Ca<sup>2+</sup>.1

Pyruvate kinase is often used as a coupling enzyme in conjunction with lactic dehydrogenase in quantifying ADP and the activity of enzymes that catalyze the formation of ADP.

#### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

### Storage/Stability

This enzyme is offered as a crystalline suspension in 3.2 M ammonium sulfate, pH 6.0. Dilute solutions should not be stored.

#### References

- Kayne, F. J., in The Enzymes, 3rd ed., Vol. 8, pt. A, Boyer, P. D., ed., Academic press (New York, NY: 1973), pp. 353-382.
- 2. Methods of Enzymatic Analysis, 2nd English ed., Vol. 1, Bergmeyer, H. U., ed., Academic Press (New York, NY: 1974), pp. 509-511.
- 3. Bucher, T., and Pfleiderer, G., Pyruvate kinase from muscle. Meth. Enzymol., 1, 435-440 (1955).
- Ainsworth, S., et al., The regulatory properties of rabbit muscle pyruvate kinase. The influence of substrate concentrations. Biochem. J., 209(2), 401-411 (1983).

TMG/RXR 1/03