





Can I use Benzonase® Nuclease to make my protein experimental results perfectly clear?

Yes! By completely digesting nucleic acids,
Benzonase® Nuclease reduces the viscosity
of your samples – yet is perfectly safe for proteins.
So you can enjoy sharper SDS-PAGE gels,
better purifications, faster chromatography flow rates...
and even a clearer mindset.

That's what's in it for you. EMD Chemicals

### QUICK REFERENCE

Product Description	3
Enzyme Characteristics	4
Applications	6
Frequently Asked Questions	10
Related Products	11
User Citations	14

There is **only one effective biochemical method** to remove DNA and RNA. It is called **Benzonase® Nuclease**.

**First commercialized** by us over 15 years ago, Benzonase Nuclease is still manufactured only by **Merck KGaA Darmstadt, Germany**. Benzonase Nuclease has **proven its value** in the laboratory.

Why risk it with others...
....accept no substitute!

## Benzonase® Nuclease

Works so well that there is no need for any other nuclease.

nacica

You don't need to add a lot of extra 'protein' to your prep to get the results you want.

Free of detectable proteolytic activity For safe protein production processes and everyday applications.

Effective over a wide range of

Exceptionally high specific activity

operating conditions

Allows greater flexibility of reaction conditions to accommodate your protein needs.

Benzonase® Nuclease is a genetically engineered endonuclease from *Serratia marcescens*. It degrades all forms of DNA and RNA (single-stranded, double-stranded, linear and circular) while having no proteolytic activity. The enzyme consists of two subunits of 30kDa each. It is effective over a wide range of conditions and possesses an exceptionally high specific activity. It completely digests nucleic acids to 5'-monophosphate terminated oligonucleotides 3 to 5 bases in length (below the hybridization limit), which is ideal for removal of nucleic acids from recombinant proteins enabling compliance with FDA guidelines for nucleic acid contamination.

Product	Size	Cat. No.
Benzonase® Nuclease, Purity >99%	10 KU	70664-3
Benzonase® Nuclease HC, Purity >99%	25 KU	71206-3
Benzonase® Nuclease, Purity >90%		70746-4 70746-3
Benzonase® Nuclease HC, Purity >90%	25 KU	71205-3



For Bulk and Custom options, please contact Technical Services:

Phone 800.628.8470 bioscienceshelp@emdchemicals.com

Ideal for a variety of applications including:

- 1 | Elimination of nucleic acid and viscosity from recombinant proteins
- 2 | Enhanced protein purification via complete cell lysis and viscosity reduction
- 3 | Increased gel resolution
- 4 | Prevention of cell clumping



For current pricing on Benzonase® Nuclease or to access additional information, visit...






## **Applications**

# 1 | Elimination of nucleic acids and viscosity from recombinant proteins

"Traditional" methods for the removal of interfering nucleic acids, such as extraction, sonication or precipitation, are not specific. These methods can be harsh and result in a loss in yield of protein due to denaturation and precipitation effects.

Many recombinant proteins have requirements concerning residual nucleic acids. This is a typical field of application for the use of Benzonase® Nuclease. In particular, Benzonase Nuclease completely digests nucleic acids (all kinds of DNA and RNA) to 5'-monophosphate-terminated oligonucleotides 3 to 5 bases in length. This is significantly below the hybridization limit, thus enabling recombinant proteins to meet the FDA guidelines for nucleic acid contamination.

Also, the ability of Benzonase Nuclease to rapidly hydrolyze nucleic acids makes the enzyme an ideal choice for reducing cell lysate viscosity consistently. This results in reduced processing times, increased protein yields, and increased efficiency of chromatographic purification steps.



Destroy the

nucleic acids,

not your protein!


# 2 | Enhanced protein purification via complete cell lysis and viscosity reduction

When purifying proteins for functional or structural proteomics studies, the first step is to disrupt the cells or tissue samples and extract the relevant protein fraction. The purification process should be carried out without any time delay in order to minimize losses of the target protein by proteolytic attack.

However, cell extracts often show high viscosity caused by the release of nucleic acids during the disruption of cells. The addition of Benzonase® Nuclease rapidly hydrolyzes nucleic acids resulting in viscosity reduction thus reducing processing times and increasing protein yields.

For further purification enhancement and protein yield from *E. coli* cultures, the enzyme is compatible with BugBuster® and PopCulture® Protein Extraction Reagents enabling complete cell lysis, nucleic acid digestion and viscosity reduction all in the same step. This provides a simple, rapid, low-cost alternative to mechanical methods such as French press or sonication for releasing expressed target protein in preparation for purification.

Benzonase Nuclease can also be added to eukaryotic extracts obtained with CytoBuster<sup>TM</sup> and Insect PopCulture Protein Extraction Reagents to achieve the same protein purification enhancements seen with  $E.\ coli$  extracts.



#### Viscosity reduction by Benzonase Nuclease

E. coli BL21(DE3) cells containing a pET construct were suspended in BugBuster Reagent (5 ml/g wet weight). Aliquots of the suspension were treated with the indicated amounts of Benzonase Nuclease for 10 min at room temperature, centrifuged at 350 x q for 3 min and photographed.

Enhanced IMAC purification results using the combination of Benzonase® Nuclease and BugBuster® Protein Extraction Reagent					
Media	Extraction	Flow rate ml/min	Recovery total mg	Media capacity mg/ml	Purity
His•Bind® Resin¹	BugBuster only	0.7	32.1 mg	8.0	84%
His•Bind Resin¹	BugBuster + Benzonase	7.5	37.1 mg	9.3	88%
Ni-NTA His•Bind Resin²	BugBuster only	0.5	31.6 mg	7.9	85%
Ni-NTA His•Bind Resin²	BugBuster + Benzonase	5.1	44.4 mg	7.9	90%

- <sup>1</sup> IDA Ni<sup>2+</sup> chelation. Can be regenerated many times with no loss in performance. Note that chelating or other reducing agents can produce poor purification results due to higher metal leaching rates vs. NTA resins.
- 2 NTA Ni<sup>2+</sup> contain additional chelation site which minimizes leaching and is compatible with 20 mM β-mercaptoethanol for reduction of disulfide bonds.

After expression of a recombinant protein in E. coli, the culture was split in half. One half of the cells were lysed using BugBuster® Protein Extraction Reagent and the other half was lysed with BugBuster plus Benzonase® Nuclease following the recommend protocol. Following the extraction, the extracts were equally split into two and the protein was purified using an either a HiseBind® Resin or a Ni-NTA HiseBind Resin column while measuring the flow rate of the columns. Subsequently, the yields and purities were determined.



For current pricing on Benzonase® Nuclease or to access additional information, visit...








## Applications (continued)

## 3 | Increased gel resolution

Two-dimensional gel electrophoresis is a powerful technique in proteomics for the high resolution separation of complex mixtures of proteins. Nucleic acids are negatively charged molecules. They tend to form complexes by electrostatic interactions with positively charged domains on the surface of proteins. The formation and shape of these adducts usually cannot be predicted. These complexes behave completely different when migrating in an electric field than would be expected for a pure protein. Besides other effects, this may lead to band shifts in the typical protein patterns and cause poor resolution in two-dimensional gel electrophoresis. Pre-treatment of the sample with Benzonase® Nuclease will significantly improve the resolution of electrophoretic separation as demonstrated in the figure below.



Silver-stained two dimensional gel electrophoretic spearation of bacterial cells (Proteus vulgaris) a. without use of Benzonase Nuclease

b. by use of Benzonase Nuclease

Note: The presence of Benzonase Nuclease could not be detected on the gel.

See more clearly, so you don't get cloudy results!



### **4** | Prevention of cell clumping

A novel application of Benzonase® Nuclease is its incorporation into cell culture media to prevent cell clumping, especially when thawing frozen cell samples. Not only is Benzonase Nuclease free of protease activity, but it poses no threat to healthy cells making it ideal for such a role.

Recent studies demonstrate utility of Benzonase Nuclease in developing vaccines against some of the world's most serious infectious diseases. Peripheral Blood Mononuclear Cells (PBMCs) isolated from whole blood have important applications in vaccine research, for example, during quantification of vaccine-induced T-cell responses. T-cell response assay methods initially require freshly isolated cells for optimal signal detection, posing a serious practical limitation for sample handling during large clinical trials. However, frozen PBMCs (particularly PBMCs prepared from stored blood) tend to clump together upon thawing, preventing further analysis.

In 2001, Smith et al. (1) demonstrated that inclusion of Benzonase Nuclease in PBMC thawing buffer prevented cell clumping, allowing implementation of PBMC cryopreservation. This method has since been widely adopted in vaccine evaluation studies (2,3).



For current pricing on Benzonase® Nuclease or to access additional information, visit...

www.emdbiosciences.com/ Benzonase

When your cells just want to be alone...





<sup>1.</sup> Smith, J.G., et al. 2001. Clin. Diagn. Lab. Immunol. 8, 871.

<sup>2.</sup> Huaman, M.C., et al. 2008. J. Immunol. 180, 1451.

<sup>3.</sup> Bull, M., et al. 2007. J. Immunol. Meth. 322, 57

## **Enzyme Characteristics**

We are able to meet the widest possible range of processing requirements because we are the only manufacturer of Benzonase® Nuclease. This enables us to offer the most options for your Benzonase Nuclease needs. Benzonase Nuclease is available in two different grades of purity and concentrations:

Specific

activity

1 x 10<sup>6</sup> units/

1 x 10<sup>6</sup> units/

1 x 106 units/

1 x 10<sup>6</sup> units/

mg protein

mg protein

mg protein

mg protein

Comments

SDS-PAGE

Produced by chromatographic

purification; produces a single

Produces the dominant band by

band by SDS-PAGE

**Endotoxins** Concentration levels Benzonase Nuclease HC, Ultra Pure 250 U/μI <0.25 Purity >99% EU/1000 units Ultra Pure 25 U/μI <0.25 Benzonase Nuclease. EU/1000 units Purity >99% Benzonase Nuclease HC. Pure 250 U/µI Not tested Purity >90% Benzonase Nuclease, Pure 25 U/μI Not tested Purity >90%

### **Reaction Conditions**

Benzonase Nuclease retains its activity in a wide range of reaction conditions, as specified in the following table and figures 1 to 3. The enzyme is also active in the presence of ionic and non-ionic detergents, reducing agents, PMSF (1 mM), EDTA (1 mM) and urea (relative activity depends on specific conditions). Activity is inhibited by >300 mM monovalent cations, >100 mM phosphate, >100 mM ammonium sulfate, or >100 mM guanidine HCl.

Condition	Optimal*	Effective**
Mg <sup>2+</sup> concentration	1 – 2 mM	1 – 10 mM
pH	8.0 - 9.2	6.0 – 10.0
Temperature	37 °C	0 – 42 °C
Dithiothreitol (DTT)	0 – 100 mM	> 100 mM
B-Mercaptoethanol	0 – 100 mM	> 100 mM
Monovalent cation concentration (Na+, K+, etc.)	0 – 20 mM	0 – 150 mM
PO <sub>4</sub> <sup>3-</sup> concentration	0 – 10 mM	0 – 100 mM

<sup>\* &</sup>quot;Optimal" is defined as the condition in which Benzonase Nuclease retains ≥90 % of its activity.


BUIK

For Bulk and Custom options, please contact Technical Services:

Phone 800.628.8470 bioscienceshelp@emdchemicals.com

<sup>\*\* &</sup>quot;Effective" is defined as the condition in which Benzonase Nuclease > 15 % of its activity.

### Effect of detergents (Triton® X-100, sodium deoxycholate) on Benzonase® Nuclease activity\*




- Concentrations of Triton® X-100 < 0.4% have no effect on the activity of Benzonase Nuclease.
- At concentrations of sodium deoxycholate < 0.4%, Benzonase Nuclease retains at least 70% of its activity.




For current pricing on Benzonase® Nuclease or to access additional information, visit...

www.emdbiosciences.com/ Benzonase


### Effect of different SDS concentrations on Benzonase Nuclease activity<sup>†</sup>






- Benzonase Nuclease retains 100% of its activity in SDS concentrations up to 0.05%.
- At SDS concentrations between 0.1% and 1%, Benzonase Nuclease remains active for a short period of time before being denatured. This can be partially compensated by increasing the concentration of Benzonase Nuclease.

### Effect of different concentrations of urea on Benzonase Nuclease activity\*†





- Benzonase Nuclease is activated by urea of concentrations up to approximately 6M.
- At 6M urea, enzyme activity first increases, then decreases over time.
- At 7M urea, Benzonase Nuclease denatures after 15 min, and activity is lost. However, significant degradation of nucleic acids occurs before the enzyme is inactivated.
- Higher concentrations of Benzonase Nuclease can partially compensate the effects of 7M urea.

Reaction conditions defined by your needs, not your enzymes.





\* Standard assay conditions = 100% activity. Substrate is DNA. +Benzonase Nuclease concentrations are given in units/ml.



Kari Weber, Ph.D. EMD Chemicals – Novagen

# Frequently Asked Questions What you need to know about Benzonase® Nuclease

Which quality of Benzonase Nuclease will be adequate for certain applications? What does HC mean?

There are several parameters that influence the activity of Benzonase Nuclease, such as the incubation time, incubation temperature, and the quantity used. Hence, the optimal purity/ amount of the nuclease and the optimal conditions will vary from process to process and needs to be determined experimentally. For viscosity reduction, Benzonase Nuclease, Purity >90% will often be sufficient, especially if *E. coli* is used as the expression system. In addition to the standard concentration of 25 U/ $\mu$ l, Benzonase Nuclease is also provided in a High Concentration (HC) formulation of 250 U/ $\mu$ l.

What type of nucleic acids does Benzonase Nuclease work on? Can I use it when isolating RNA? What is the end result of the digestion?

Benzonase Nuclease is a promiscuous endonuclease that attacks and degrades all forms of DNA and RNA (single stranded, double stranded, linear, circular, and super-coiled). Upon complete digestion, all nucleic acids present in the solution are reduced to 5'-monophos-phate-terminated oligonucleotides 3 to 5 bases in length.

### For what purposes would I want to use Benzonase Nuclease?

Some applications include: elimination of nucleic acids, reduction of viscosity to facilitate processing (e.g., recombinant protein purification, viscosity-sensitive applications using mammalian cell lysates, etc.); reduction of cell clumping in stored peripheral blood mononuclear cell (PBMC) samples prior to processing (see Perrin, 2008 *inNovations* 28, 21); preparation of inclusion bodies to allow successful protein renaturation; and removal of negatively charged nucleic acids from samples prior to two-dimensional SDS-PAGE.

How can Benzonase Nuclease be inactivated? How can it be removed? Reversible inhibition can be achieved using EDTA to chelate essential metal ions. Irreversible inactivation can only be accomplished with extreme conditions (100 mM NaOH at 70°C for 30 min). The enzyme can easily be separated from the target product using chromatography. Please note that Benzonase Nuclease consists of two subunits of 30 kDa each. Because of the robust nature of this endonuclease, we recommend that Benzonase Nuclease

not be used if a nuclease-free end product is required. For more information, visit www.novagen.com/benzonase.

My Benzonase Nuclease was left out on the bench all weekend. Is it still good? We have done extensive stability testing on Benzonase Nuclease, and find that it is extremely stable. Even with extended incubations at 37°C, Benzonase Nuclease maintained >90% activity for several months. For more information on stability testing, go to www.novagen.com/benzonase.

I want to use a different buffer. What conditions are absolutely required for full activity of Benzonase® Nuclease? What will reduce its activity?

Benzonase Nuclease requires 1-2 mM Mg²+ for activity. It is inhibited (approximately 50% activity) by monovalent cation concentrations >300 mM, phosphate concentrations >100 mM, and by ammonium sulfate concentrations >100 mM. In addition, concentrations of >1mM EDTA will also inhibit Benzonase Nuclease activity.

Is Benzonase Nuclease compatible with protease inhibitor cocktails?

Benzonase Nuclease is compatible with protease inhibitors. However, caution should be exercised, since many protease inhibitor cocktails include EDTA. As mentioned previously, concentrations of >1 mM EDTA will inhibit Benzonase Nuclease activity.

My protein is insoluble and I need to perform purification under denaturing conditions. Will Benzonase Nuclease still work in urea?

Benzonase Nuclease activity actually increases in presence of urea at concentrations up to 6 M (see page 9). At 6 M urea, enzyme activity first increases, then decreases over time. At 7 M urea, Benzonase Nuclease denatures after 15 min, and activity is lost. However, significant degradation of nucleic acids occurs before inactivation. Higher initial concentrations of Benzonase Nuclease can partially compensate the effects of 7 M urea.

### How much Benzonase Nuclease should I use?

One unit of Benzonase Nuclease is defined as the amount of enzyme that results in an  $\Delta A_{260}$  of 1.0 in 30 min, which corresponds to complete digestion of 37  $\mu g$  DNA. When used in conjunction with BugBuster® or CytoBuster™ Protein Extraction Reagents, we recommend the addition of 25 U Benzonase Nuclease per 1 ml cell lysate (assuming BugBuster or CytoBuster reagents are used at standard concentrations). Conditions may need to be altered when working at lower temperature, as is the case with most mammalian cell lysis protocols. For example, when using the Benzonase Nuclease at 4°C, we recommend increasing the duration of the digestion.



For current pricing on Benzonase® Nuclease or to access additional information, visit...





## **Related Products**

Being the gold standard in protein research, we also offer many other valuable tools to accommodate your protein purification needs. To propel your research forward, don't forget to ask your local technical sales specialist about how else we can help!

### BugBuster® Master Mix

BugBuster Master Mix combines BugBuster Protein Extraction Reagent with Benzonase® Nuclease and rLysozyme™ Solution into one convenient reagent. BugBuster® Protein Extraction Reagent is formulated to gently disrupt the cell wall of *E. coli* and liberate soluble proteins. It provides a simple, rapid, low-cost alternative to mechanical methods such as French press or sonication for releasing expressed target protein in preparation for purification or other applications. The proprietary formulation utilizes a detergent mix that is capable of cell wall perforation without denaturing soluble protein. BugBuster Master Mix allows for maximum recovery of active, soluble protein from both Gram-negative and Gram-positive bacteria. BugBuster Master Mix combines the activities of BugBuster Reagent, rLysozyme solution and Benzonase Nuclease, thus significantly increasing protein extraction efficiency and facilitating downstream processing of protein extracts. The two available package sizes provide sufficient reagents for protein extraction from 20 g and 100 g cell paste.


### Lysonase<sup>™</sup> Bioprocessing Reagent

Lysonase™ Bioprocessing Reagent is an optimized, ready-to-use blend of rLysozyme Solution and Benzonase Nuclease. rLysozyme Solution contains a highly purified and stabilized recombinant lysozyme with specific activity 250 times greater than that of chicken egg white lysozyme. The combined activities of rLysozyme and Benzonase Nuclease significantly increase protein extraction efficiency and facilitate downstream processing of protein extracts. For efficient protein extraction with BugBuster Protein Extraction Reagent, use 10 µl Lysonase per 1 g cell paste. In addition, Lysonase can be used to enhance the effectiveness of non-detergent based cell lysis procedures.

### ProteoExtract® Complete Protein Extraction Kit

The ProteoExtract® Complete Proteome Extraction Kit (C-PEK) is designed for fast and easy extraction of total proteins from mammalian cells and tissues, without the need for sonication or precipitation. C-PEK provides a straightforward two-step isolation of complete proteomes in a single microcentrifuge tube. It uses optimized extraction reagents to provide improved protein solubilization, resulting in an increased number of resolved spots on 2D gels. The procedure uses proprietary Benzonase Nuclease, leading to clear, nonviscous protein solutions and improved resolution on 2D gels. Using the C-PEK procedure, protein concentration is not necessary. Extracted proteins are ready

for immediate use in standard downstream proteomics applications. Each kit contains all the reagents needed for 20 reactions  $(1-2 \times 10^8 \text{ cells per reaction})$ .



Different samples show high resolution in 2D gel electrophoresis.

Product	Size	Cat. No.
BugBuster® Master Mix	100 ml	71456-3
	EOO ml	71.4EC 2

Product	Size	Cat. No.
Lysonase™ Bioprocessing	0.2 ml	71230-3
Reagent	1 ml	71230-4
	5 x 1 ml	71230-5

Product	Size	Cat. No.
ProteoExtract® Complete	1 kit	539779
Mammalian Proteome		
Extraction Kit		

### ProteoExtract® Subcellular Protein Extraction Kit

ProteoExtract Subcellular Proteome Extraction Kit (S-PEK) is designed for fast and reproducible extraction of subcellular proteomes from adherent and suspension-grown mammalian cells. The S-PEK takes advantage of the different solubilities of certain subcellular compartments in the four selected reagents. The extraction procedure utilizes four extraction buffers prepared with ultra-pure chemicals, a protease inhibitor cocktail to prevent protein degradation during the extraction procedure, and Benzonase® Nuclease to degrade contaminating nucleic acids. In the case of adherent cells, the procedure is performed directly in the tissue culture dish without the need for cell removal. Cells or the parts of the cells remain attached to the plate during sequential extraction of subcellular compartments, until the appropriate extraction reagent is used. Thus, the early destruction of the cellular structure by enzymatic or mechanical detachment of cells from the tissue culture plate and any mixing of different subcellular compartments is prevented. For suspension-grown cells, extraction starts with gentle sedimentation and washing of the cells. For tissues, fragmentation is required before proceeding with the extraction protocol. The stepwise extraction delivers four distinct protein fractions from one sample: cytosolic fraction; nucleic protein fraction; membrane/organelle protein fraction; and cytoskeletal fraction.

Proteins are obtained in the native state (except the cytoskeletal fraction), making the S-PEK suitable for many downstream applications such as 1D and 2D gel electrophoresis, immunoblotting, enzyme activity assays, and protein microarrays. Each kit contains all reagents needed for 20 reactions (Cat. No. 539790) or 10 reactions (Cat. No. 539791). Sample size:  $1 \times 10^6$  or 25-50 mg tissue.

Product

ProteoExtract®

Subcellular Proteome Extraction Kit ProteoExtract® Subcellular Proteome

Extraction Kit, Mini

Cat. No.


539790

539791

1 kit

1 kit


Fraction 1: Cytosol



Fraction 2: Membrane/Organelle

Fraction 4: Cytoskeleton

		700	13162	11/22	riĝi
	1.00	3402		~ 457.) ·	-
-		2.55		ter.	-
-		Company	100	-	25
		. 312	411		15
					Ŧ
			114	137.5	10.7
		- 22			



Protein profiling by 2D gel electrophoresis

The comparison of fractions 1 through 4 by two-dimensional SDS-PAGE indicates a large number of protein spots that are specific to each respective subcellular fraction.



For current pricing on Benzonase® Nuclease or to access additional information, visit...







As discussed earlier, Benzonase® Nuclease is ideal for a variety of applications including reduction of viscosity caused by nucleic acids, purification of proteins, sample preparation for electrophoresis and chromatography and prevention of cell clumping in the cell culture media. Below are selected citations that demonstrate the use of Benzonase Nuclease for each of these applications.

## Elimination of Nucleic Acids and Viscosity from Recombinant Proteins

Shi, H., et. al. (2009) Mitogen-activated protein kinase signaling is necessary for the maintanence of skeletal muscle mass. *Am. J. Physiol. Cell Physiol.* **296**, C1040.

Hervas-Stubbs , S., et al. (2007) Insect baculoviruses strongly potentiate adaptive immune responses by inducing type I IFN. *J. Immuno.* **178**, 2361.

Kercher, L, et al. (2007) Prion protein (PrP) expression differences in microglia and astroglia influence scrapie-induced neurodegeneration in retina and brain of transgenic mice. *J. Virology* **81**, 10340.

Cho, S., et al. (2005) Glucocorticoid receptor ligand binding domain is sufficient for the modulation of glucocorticoid induction properties by homologous receptors, coactivator transcription intermediary factor 2, and Ubc9. *Mol. Endocrinology* **19**, 290.

Han, W., et al. (2005) C-terminal ECFP fusion impairs synaptotagmin 1 function: crowding out synaptotagmin 1. *J. Biol. Chem.* **280**, 5089–5100.

Meyn, M., et al. (2005) Src family kinase activity is required for murine embryonic stem cell growth and differentiation. *Mol. Pharmacol.* **68**, 1320.

Murphy, M., et al. (2003) An improved method for the *in vitro* evolution of aptamers and applications in protein detection and purification. *Nucleic Acids Res.* 31, e110.

Persad, S., et al. (2001) Regulation of protein kinase B/ Akt-serine 473 phosphorylation by integrin-linked kinase. J. Bio. Chem. 276, 27462.

## Enhanced Protein Purification via Complete Cell Lysis and Viscosity Reduction

Collin, R., et. al. (2009) The DEXD/H-box RNA Helicase DDX19 is regulated by an  $\alpha$ -helical switch. *J. Biol. Chem.* **284**, 10296.

Gustin, K., et. al. (2009) Bovine coronavirus nsp1 (p28) is an RNA binding protein that binds terminal genomic cis-replication elements. *J. Viol.* **10.1128/JVI**, 00160.

Wanyiri, J., et. al. (2009) Role CpSUBI, a subtillisin-like protease, in *Cryptosporidium parvum* infection *in vitro*. *Eukaryot*. *Cell.* **8**, 470.

Shen, X., et. al. (2009) Identification of novel SHPS-I-associated proteins and their roles in regulation of insulin-like growth factors-dependent responses in vascular smooth muscle cells. *Mol. Cell. Proteomics.* 10.1074/MCP, M800543.

Heikkilä, O., et. al. (2009) Integrin  $\alpha V \beta 6$  is a high-affinity receptor for coxsackievirus A9. *J. Gen. Virol.* **90**, 197.

Masip, L., et al. (2008) Laboratory evolution of *Escherichia coli* thioredoxin for enhanced catalysis of protein oxidation in the periplasm reveals a phylogenetically conserved substrate specificity determinant. *J. Bio. Chem.* **283**, 840.

Suzuki, C., et al. (2008) PDCD4 inhibits translation initiation by binding to eIF4A using both its MA3 domains. *Proc. Natl. Acad. Sci.* **105**, 3274.

Lloyd-Burton, S., et al. (2007) Regulation of Ins(1,4,5)P3 3-kinases by calcium and localization in cells. J. Biol. Chem. 282, 9526.

Park, M., et al. (2007) Glyceraldehyde-3-phosphate ferredoxin oxidoreductase from *Methanococcus maripaludis*. *J.Bacteriol.* **189**, 7281.

Kobayashi, S., et al. (2007) Serine 64 phosphorylation enhances the antiapoptotic function of Mcl-1. *J. Biol. Chem.* **282**, 18407.

Miyazaki, K., et al. (2006) Thermal stabilization of *Bacillus subtilis* family-11 xylanase by directed evolution. *J. Biol. Chem.* **281**, 10236.

Rangarajan, E., et al. (2005) Crystallographic trapping of the glutamyl-CoA thioester intermediate of family I CoA transferases. *J. Biol. Chem.* **280**, 42919.

Frankel, P., et al. (2005) RalA interacts with ZONAB in a cell density-dependent manner and regulates its transcriptional activity. *Eur. Mol. Bio. Org. J.* **24**, 54–62.

Johnsen, L., et al. (2005) 1.6 Å crystal structure of EntA-Im: A bacterial immunity protein conferring immunity to the antimicrobial activity of the pediocin-like bacteriocin enterocin A. *J. Biol. Chem.* **280**, 19045.

Kumaran, D., et al. (2005) Structure and mechanism of ADP-ribose-1"-monophosphatase (Appr-1"-pase), a ubiquitous cellular processing enzyme. *Prot. Sci.* **14**, 719.

Li, L., et al. (2005) Clarification of cinnamoyl co-enzyme a reductase catalysis in monolignol biosynthesis of aspen. *Plant Cell Physiol.* **46**, 1073.

Yamakoshi, Y., et al. (2005) Porcine dentin sialoprotein is a proteoglycan with glycosaminoglycan chains containing chondroitin 6-sulfate. *J. Biol. Chem.* **280**, 1552.

Zanghi, C., et al. (2005) A simple method for displaying recalcitrant proteins on the surface of bacteriophage lambda. *Nucleic Acids Res.* **33**, e160.

Kasim, V., et al. (2004) Control of siRNA expression using the Cre-loxP recombination system. *Nucleic Acids Res.* **32**, e66.

Lunin, V., et al. (2004) The structure of the MAPK scaffold, MP1, bound to its partner, p14 A complex with a critical role in endosomal MAP kinase signaling.

J. Biol. Chem. 279, 23422.

Rizzo, M., et al. (2004) An improved cyan fluorescent protein variant useful for FRET. *Nat. Biotech.* 22, 445.

Chan, W., et al. (2002) Human phosphatidylcholine transfer protein: purification, crystallization and preliminary X-ray diffraction data. *Biochim. Biophys Acta* **1596**, 1.

Huhtinen, K., et al. (2002) The peroxisome proliferator-induced cytosolic type I acyl-CoA thioesterase (CTE-I) is a serine-histidine-aspartic acid  $\alpha/\beta$  hydrolase. *J. Biol. Chem.* **277**, 3424.

Lobel, L., et al. (2002) Bacterial expression of a natively folded extracellular domain fusion protein of the hFSH receptor in the cytoplasm of Escherichia coli. *Protein Expr. Purif.* **25**, 124.

Peitz, M., et al. (2002) Ability of the hydrophobic FGF and basic TAT peptides to promote cellular uptake of recombinant Cre recombinase: a tool for efficient genetic engineering of mammalian genomes. *Proc. Natl. Academy of Science* **99**, 4489.

Suarez, T., et al. (2002) Biochemical defects in retina-specific human ATP binding cassette transporter nucleotide binding domain 1 mutants associated with macular degeneration. *J. Biol. Chem.* 277, 21759.

Glick, E., et al. (2001) Mutations in human DNA polymerase  $\eta$  motif II after bypass of DNA lesions. *Eur. Mol. Bio. Org. J.* **20**, 7303

Lobel, L., et al. (2001) High-level bacterial expression of a natively folded, soluble extracellular domain fusion protein of the human luteinizing hormaone/ chorionic gonadotropin receptor in the cytoplasm of *Escherichia coli. Endocrine* 14, 205

### Increased Gel Resolution

Milovnik, P., et. al. (2009) Selection and characterization of DARPins specific for the neurotensin receptor I. *Protein Eng. Des. Sel.* 10.1093, 1.

Oliveira, J., et al. (2007) NRAS mutation causes a human autoimmune lymphoproliferative syndrome. *Proc. Natl. Acad. Sci.* **104**, 8953.

Kwade, Z., et al. (2005) Identification of four adenosine kinase isoforms in tobacco by-2 cells and their putative role in the cell cycle regulated cytokinin metabolism. *J. Biol. Chem.* **280**, 17512.

Poirier, M., et al. (2005) A structure-based analysis of huntingtin mutant polyglutamine aggregation and toxicity: evidence for a compact beta-sheet structure. *Human Mol. Gen.* 14. 765.

Mercer, K.L.N. and Weiss, D. (2002) The *Escherichia coli* cell divisionn protein FtsW is required to recruit its cognate transpeptidase, Ftsl (PBP3), to the division site. *J. Bacteriol.* **184.** 904.

### | Prevention of Cell Clumping |

Doria-Rose, N., et. al., (2009) Frequency and phenotype of human immunodeficiency virus envelope-specific B cells from patients with broadly cross-neutralizing antibodies. *J. Virol.* **83**, 188.

Huaman, M., et al. (2008) "Ex vivo cytokine and memory T cell responses to the 42-kDa fragment of *Plasmodium falciparum* merozoite surface protein-1 in vaccinated volunteers." *J. Immuno.* **180**, 1451.

Mingozzi, F., et al. (2007) Modulation of tolerance to the transgene product in a nonhuman primate model of AAV-mediated gene transfer to liver. *Blood* 110, 2334.

Keating, S., et al. (2005) Durable human memory T cells quantifiable by cultured enzyme-linked immunospot assays are induced by heterologous prime boost immunization and correlate with protection against malaria. *J. Immunol.* 175,

Vuola, J., et al. (2005) Differential immunogenicity of various heterologous prime-boost vaccine regimens using DNA and viral vectors in healthy volunteers. *J. Immunol.* **174**, 449.

For current pricing on Benzonase® Nuclease or to access additional information, visit...

www.emdbiosciences.com/ Benzonase





Prices and availability are subject to change. ©Copyright 2008 EMD Chemicals, an affiliate of Merck KGaA, Darmstadt, Germany. All rights reserved. Each product is sold with a limited warranty which is provided with each purchase. Each product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. EMD Chemicals products may not be resold, modified for resale, or used to manufacture commercial products without written approval of EMD Chemicals. His-Pesin®, BugBuster®, PopCulture® and Novagen® are registered trademarks of EMD Chemicals Inc. in the United States and in certain other jurisdictions. CytoBuster® and Lysonase® are trademarks of EMD Chemicals, Inc. Benzonase® and ProteoExtract® are registered trademarks of Merck KgaA, Darnstadt, Germany. Triton® is a registered trademark of Dow Chemical Company.

The use of Benzonase Nuclease is covered under U.S. Patent 5,173,418 and European Patent 0,229,866. Nycomed Pharma A/S, Denmark holds worldwide patent rights to Benzonase Nuclease and Merck KGaA holds an exclusive license on all patents.

For more information or to order Novagen products, contact EMD

Orders

Phone 800 854 3417 Fax 800 776 0999

*Technical Support*Phone 800 628 8470

E-mail bioscienceshelp@emdchemicals.com

Visit our website

www.emdbiosciences.com

You can also order Novagen products through our distribution partner, VWR International

Phone 800 932 5000 Web www.vwr.com

Printed in the USA



10003362 USD