3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**



# **Western Blot Procedure**

#### **Product Description**

Some Prestige Antibodies<sup>®</sup> are verified as primary reagents for Western blot analysis using the described procedure.

## **Preparation Instructions**

Wash Buffer – 10 mM Tris-HCl with 150 mM NaCl and 0.05% (v/v) TWEEN<sup>®</sup> 20, pH 7.5

Blocking Buffer - 10 mM Tris-HCl with 150 mM NaCl, 5% non-fat dried milk, and 0.5% (v/v) TWEEN 20, pH 7.5

### Electrophoresis and Electroblotting

Protein samples to be examined, such as human plasma and lysates from selected human tissues and cell lines, are separated by SDS-PAGE on an appropriate gel, e.g., Tris-HCl gel, 10–20%. An average of 20 µg of total protein is loaded per lane onto the gel. The protein is transferred onto a PVDF or nitrocellulose membrane by electroblotting. After transfer, the membrane is completely dried to increase the protein retention prior to immunoblotting.

#### **Procedure**

All washes are performed at room temperature (RT) on a shaker.


Antibody dilution: 1:250-1:500

<u>Note</u>: The specified working dilution of the primary antibody is to be considered as a guideline only. Optimal dilution must be determined by the user.

- 1. PVDF membranes are soaked in methanol.
- Non-specific sites are blocked on both PVDF or nitrocellulose membranes with Blocking Buffer for 45 minutes at RT or overnight at 4 °C.
- 3. The membrane is quickly rinsed in Wash Buffer followed by incubation for 1 hour at RT in the primary antibody diluted in 5 ml of Blocking Buffer.

- 4. The membrane is rinsed 4 times, 5 minutes each wash, with large volumes of Wash Buffer.
- The secondary antibody horseradish peroxidase conjugate is diluted 3,000-fold in 5 ml of Blocking Buffer and the membrane is incubated for 1 hour at RT.
- 6. The membrane is washed as in step 4.
- Excess Wash Buffer is drained from the membrane, which is then placed in a chemiluminescent substrate solution and incubated for 30 seconds at RT in darkness.
- 8. Excess detection reagent is drained off and a CCD-camera is used for detection and capture of a digital image

Figure 1.



Prestige Antibodies is a registered trademark of Sigma-Aldrich® Biotechnology LP and Sigma-Aldrich Co. TWEEN is a registered trademark of Croda International PLC.

KB,LPG,MAM 06/09-1

Sigma brand products are sold through Sigma-Aldrich, Inc.