

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Polynucleotide phosphorylase from *Escherichia coli*, histidine tagged recombinant, expressed in *Escherichia coli*

Catalog Number **N9664** Storage Temperature –70 °C

EC 2.7.7.8

Synonyms: Polyribonucleotide nucleotidyltransferase, PNPase

Product Description

Polynucleotide phosphorylase (PNPase) is a bifunctional enzyme with a 3'-terminal oligonucleotide polymerase activity and a phosphorolytic 3' to 5' exoribonuclease activity. It is also involved in mRNA processing and degradation in bacteria, plants, and humans. The subunit molecular mass is 84–86 kDa. The native protein is present as a trimeric homopolymer with a molecular mass of 230 kDa. 3.4

Polynucleotide phosphorylase (PNPase) is present in the chloroplasts and mitochondria of some eukaryotic cells. The enzyme is a functional part of the "degradosome", a multienzyme complex (molecular mass ~500 kDa). ^{5,6} PNPase was shown to protect *E. coli* against oxidative stress by specifically binding to RNA molecules that were oxidatively damaged. ⁷ The *E. coli* PNPase enzyme was shown to polymerize all ribonucleotides in the presence of divalent cations such as Mg²⁺ or Mn²⁺ for GDP. ^{7,8} The affinity of *E. coli* PNPase to poly(G) sequences is very low and thus, this polynucleotide can be used as an effective barrier to exonuclease activity in yeast and chloroplasts. ⁵

It is important to mention the enzyme activity on different nucleotides may require different assay conditions, e.g., maximal polymerization activity was observed with ADP or IDP at 55 $^{\circ}$ C and pH 9.1; whereas, CDP and UDP were maximal at 37 $^{\circ}$ C. GDP polymerization is very sensitive to pH changes and gave the lowest activity of any of the nucleotides tested at 37 $^{\circ}$ C.

The product is supplied as a solution in 20 mM HEPES buffer, pH 7.9, with 0.1 mM EDTA, 2 mM DTT, 12.5 mM MgCl₂, 200 mM KCl, and 21.4% (w/v) glycerol.

Purity: ≥90% (SDS-PAGE)

Specific activity: ≥250 units/mg-protein

Unit definition: One unit will polymerize 1.0 μ mole of ADP releasing 1.0 μ mole of inorganic phosphate in 15 minutes at pH 9.1 at 37 °C.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the product at -70 °C. The product is stable for at least 2 years as supplied. Avoid repeated freezing and thawing. After initial thawing, the enzyme may be kept for up to 2 months at 4 °C or refrozen in aliquots at -70 °C.

References

- Yehudai-Resheff, S., et al., Polynucleotide phosphorylase functions as both an exonuclease and a Poly (A) polymerase in spinach chloroplasts. *Mol. Cell Biol.*, 21, 5408-5416 (2001).
- Liou, G.G., et al., RNA degradosomes exist in vivo in Escherichia coli as multicomponent complexes associated with the cytoplasmic membrane via the N-terminal region of ribonuclease E. Proc. Natl. Acad. Sci. USA, 98, 63-68 (2001).
- Soreq, H., and Littauer, U.Z., Purification and characterization of polynucleotide phosphorylase from *Escherichia coli*. Probe for the analysis of 3' sequences of RNA. *J. Biol. Chem.*, 252, 6885-6888 (1977).
- Kimhi, Y., and Littauer, U.Z., Purification and properties of polynucleotide phosphorylase from *Escherichia coli*, *J. Biol. Chem.*, 243, 231-240 (1968).

- Lisitsky, I., and Schuster, G., Preferential degradation of polyadenylated and polyuridinylated RNAs by the bacterial exoribonuclease polynucleotide phosphorylase. *Eur. J. Biochem.*, 261, 468-474 (1999).
- 6. Carpousis, A.J., The *Escherichia coli* RNA degradosome: structure, function and relationship in other ribonucleolytic multienzyme complexes. *Biochem. Soc. Trans.*, **30**, 150-155 (2002).
- 7. Wu, J., et al., Polynucleotide phosphorylase protects *Escherichia coli* against oxidative stress. *Biochemistry*, **48**, 2012-2020 (2009).
- 8. Thang, M.N., et al., Synthesis of Polyguanylic acid by polynucleotide phosphorylase of *Escherichia coli. Biochim. Biophys. Acta*, **108**, 125-131 (1965).

RBG,EB,MAM 12/09-1