

AXIS[™]: Axon Investigation System Spatially Controlled Outgrowth

Spatially Controlled Outgrowth Platform

www.millipore.com/AXIS

Introduction

Nervous system development results in a network of synaptic connections between participating neurons. Understanding the process of formation of this network is crucial to improving therapeutic treatments for patients suffering from nervous system developmental disorders and from neurodegenerative diseases, such as Parkinson's disease, Huntington's disease, and Lewy body disease¹. Neural connections are made when neuron filaments grow outwards in response to axon guidance cues, which include extracellular chemoattractant gradients and intracellular signaling molecules². Historically, neurite outgrowth assays have been able to detect only gross correlations between signaling events, molecular gradients and axon growth. From these studies, it has been shown that proteins such as F-actin and microtubules accumulate in the neuron's growth cone in response to axon growth cues³.

Recent research has focused on studying the causes of directionality of axon growth. The polarity of axon growth coincides with gradients of extracellular signals, but it is not yet clear how extracellular gradients translate into asymmetric distribution and function of intracellular proteins driving neurite outgrowth⁴. To fully understand the biochemical mechanisms by which axons grow in response to signaling, one should assay each axon, in isolation from somas and other neurons.


Spatial studies of directional neuronal outgrowth have been hampered by a lack of a solid, flexible, neurite culturing platform. Traditionally, outgrowth experiments are conducted either by culturing neural stem cells (NSCs), primary neurons, or tissue slices with standard cultureware. This leads to random neurite outgrowth or outgrowth in particular directions guided by specific growth factors as dictated by the experimental design. Even in growth-factor-guided situations, neurite outgrowth is often crowded and somewhat haphazard. The ability to spatially isolate and study individual neurites is difficult in these preparations due to erratic growth, clumping and bifurcation. A simple, inexpensive, repeatable way to grow and isolate neurites would greatly enhance qualitative and quantitative studies.

The AXIS platform is Millipore's most advanced tool for the study of neurite outgrowth. This slide-mounted microfluidic chamber system enables the deposition and culture of neural cells and the spatially controlled addition of growth factors, toxins, and other reagents. Neurite outgrowth is restricted to narrow, parallel channels, and the resultant outgrowth or collapse behavior is easily observed under a microscope. The result is a powerful platform for the study of somas, neurites or synaptic formation.

Cell Body Growth Cues


How the AXIS Isolation Device Works

AXIS Axon Isolation Device is a two chamber system, each composed of two wells and an interconnected channel, separated by a set of microgrooves. The hydrostatic pressure formed by volume differential between chambers induces fluidic isolation of the solution on the low volume side of the device. The microfluidic design of an AXIS device allows for development and maintenance of a fluidic gradient of chemoattractants, toxins or other molecules of interest, facilitating controlled exposure and differentiation of axons.



Features, Benefits, Advantages:

- Organize, visualize, and characterize neuronal cell culture.
- o Detect protein expression with better spatial resolution.
- o Isolate cell bodies from axons through fluidics.
- Reduce time and expense through optimized protocols and QC validated products.
- o Attain superior performance over in-house protocols.
- Optically clear transparent, inert, non-toxic, and nonflammable polymer mold.
- $\circ\,$ Available in 150 $\mu m,$ 450 $\mu m,$ 900 $\mu m,$ or 6-well.



N1E-115 cells clearly demonstrate neurite outgrowth through the AXIS channels (150 µm) using the Milli-Mark™ FluoroPan neuronal marker (MAB2300X) shown in green, versus DAPI (blue).

ORDERING INFORMATION

Description	Catalogue No.
AXIS Axon Isolation Device (150 µm)	AX15010
AXIS Axon Isolation Device (450 μm)	AX45005, AX45010
AXIS Axon Isolation Device (900 µm)	AX43010 AX90010
AXIS Axon Isolation Device (6-well)	AX50010


Related Products

Description	Catalogue No.
Neurite Outgrowth Kit	NS220
Neurite Outgrowth Plus Kit,	NS230
including antibody staining controls	
Milli-Mark™ Pan-Neuronal Marker,	MAB2300
monoclonal antibody blend	
Milli-Mark Chromapan Neuronal Marker	NS420
Milli-Mark Chromapan Neuronal Marker - Mouse Open	NS330
Milli-Mark Chromapan Neuronal Marker - Rabbit Open	NS340
Milli-Mark FluoroPan-Neuronal Marker,	MAB2300X
Alexa Fluor® 488-conjugated	
Anti-Actin	MAB1501
Anti-Actin, Alexa Fluor 488-conjugated	MAB1501X
Anti-α-Tubulin	05-829
Anti-β-III-Tubulin	MAB1637
Anti-β-III-Tubulin, Alexa Fluor 488-conjugated	CBL412X



www.millipore.com/offices

ADVANCING LIFE SCIENCE TOGETHER™ Research. Development. Production.



Triple staining of E18 rat hippocampal cell bodies, axons and growth cones using DAPI (blue), Map2 (AB15452, green) and axon specific staining with anti- β -III tubullin (MAB1637, red) using the 450 μm AXIS device.

References:

- Yaron, A. and Zheng, B. (2007). Navigating their way to the clinic: emerging roles for axon guidance molecules in neurological disorders and injury. Dev. Neurobiol. 67:1216-1231.
- 2. Chilton JK. Molecular mechanisms of axon guidance. Dev Biol. 2006 Apr 1; 292(1):13-24.
- Lin CH and Forscher P. Cytoskeletal remodeling during growth cone-target interactions. J Cell Biol. 1993 Jun;121(6):1369-1383.
- Quinn CC and Wadsworth WG. Axon guidance: asymmetric signaling orients polarized outgrowth. Trends Cell Biol. 2008 Oct; 18(12):597-603.

TO PLACE AN ORDER OR RECEIVE TECHNICAL ASSISTANCE

In the U.S. and Canada, call toll-free 1 800-Millipore (1-800-645-5476)

In **Europe**, please call Customer Service: France: 0825.045.645 Spain: 901.516.645 Option 1 Germany: 01805.045.645 Italy: 848.845.645 English UK: 0870.900.46.45

For other countries across Europe and the world, please visit **www.millipore.com/offices**.

For Technical Service, please visit www.millipore.com/techservice.

Millipore is a registered trademark of Millipore Corporation. The M mark, Advancing Life Science Together, and Milli-Mark are trademarks of Millipore Corporation. Alexa Fluor is a registered trademark of Life Technologies, Inc. Lit. No. PB2876EN00 Rev. - Printed in U.S.A. 09/09 LS-SBU-09-02370 © 2009 Millipore Corporation, Billerica, MA 01821 U.S.A. All rights reserved.