

## Product Information

### bcl-2 PCR Primers Set

Product No. **B9179**  
 Store at  $-20^{\circ}\text{C}$

#### Product Description

In recent years, several genes have been linked to apoptosis. The bcl-2 family of genes regulates PCD either positively or negatively. Bcl-2 and members of its family have been found to block apoptotic cell death. Bcl-2 protein heterodimerizes with Bax (Bcl-2 Associated X protein), which is a potent mediator of programmed cell death. The Bcl-2/Bax ratio appears to determine whether some cells live or die.<sup>1-4</sup>

The bcl-2 PCR Primers Set contains both sense and antisense primers for the amplification of the bcl-2 $\alpha$  gene. It is designed for PCR<sup>†</sup> detection of human, rat and mouse cDNA levels (representing mRNA expression) of the bcl-2 $\alpha$  apoptotic gene. No amplification of the genomic DNA has been observed.

The size of the amplified product resulting from the use of the bcl-2 PCR Primers Set is 127 bp.

#### Component

- bcl-2 PCR Primers Set, Product No. B9179 1 vial

#### Equipment and Reagents Required but Not Provided

(Sigma product numbers have been given where appropriate)

- Thermal cycler
- Taq DNA polymerase, Product No. D4545 or equivalent
- Deoxynucleotide mix, 10 mM, Product No. D7295 or equivalent
- Agarose
- Ethidium bromide, 500  $\mu\text{g}/\text{ml}$ , Product No. E1385
- PCR 100 bp low ladder, Product No. P1473
- Gel loading solutions, Product No. G2526 or G7654
- PCR grade water, Product No. W1754
- Mineral oil, Product No. M8662
- PCR microtubes, Product No. Z37,487-3 or Z37,496-2

#### Storage

Store the vial at  $-20^{\circ}\text{C}$ .

#### Preparation Instructions

The bcl-2 PCR Primers Set contains 1 nmoles of each primer (sense and antisense). Centrifuge the tube briefly in order to collect the tube contents. For the following procedure, resuspend the primers set in 100  $\mu\text{l}$  deionized water to a final concentration of 10 pmole/ $\mu\text{l}$ . Mix until the solution is homogenous. Once suspended, store the solution at  $-20^{\circ}\text{C}$ . To avoid repeated freeze-thaw cycles, aliquot the primer solution for long-term storage.

#### Procedure

**Note:** Use aseptic techniques and use aerosol barrier tips while performing PCR experiments.

1. Thaw the bcl-2 PCR Primers Set on ice, being sure that the solution is homogenous.
2. Add the following reagents to a PCR microcentrifuge tube in the following order:

|  | Amount for 50 $\mu\text{l}$ single PCR reaction | Final concentration in the PCR reaction |
|--|---|---|
| Water  | To 50 $\mu\text{l}$                             | ----                                    |
| 10X PCR Buffer                                 | 5 $\mu\text{l}$                                 | 1X                                      |
| 2 mM dNTP solution                             | 5 $\mu\text{l}$                                 | 0.2 mM of each dNTP                     |
| 25 mM $\text{MgCl}_2^*$ solution               | 3 $\mu\text{l}$                                 | 1.5 mM                                  |
| bcl-2 PCR Primers Set, 10 pmole/ $\mu\text{l}$ | 2 $\mu\text{l}$                                 | 0.4 $\mu\text{M}$                       |
| cDNA**   | 2 $\mu\text{l}$                                 | $\sim$ 30 ng                            |
| Taq DNA Polymerase, 5 units/ $\mu\text{l}$     | 1 $\mu\text{l}$                                 | 0.1 units/ $\mu\text{l}$                |
| Total volume                                   | 50 $\mu\text{l}$                                | -                                       |

\* When using the bcl-2 PCR Primers Set for the first time, you may set two additional reaction tubes with a higher and a lower  $\text{MgCl}_2$  concentrations (see Note at the end of this section).

\*\* Optimize this parameter with your own cDNA.

3. Mix gently by finger tapping and centrifuge briefly to collect the mixture in the bottom of the tube. Overlay the reaction mixture with 2 drops (~30 µl) of mineral oil to cover the surface of the reaction mixture if not using a thermal cycler with a heated lid. Place the tube in the thermal cycler when the thermal cycler reaches 95 °C, and run the following PCR program.

95 °C for 2 min  
 94 °C for 45 sec  
 53 °C for 45 sec  
 72 °C for 1.5 min  
 72 °C for 7 min

} x 30 cycles

The amplified DNA can be evaluated by agarose gel electrophoresis.

Note: Using different thermal cyclers:

For a better detection of the amplified product you may increase the number of amplification cycles. In case you do not see differences in the amount of the amplified DNA fragments, decrease the number of cycles to verify your results.

In rare cases, some of the parameters should be optimized for the specific thermal cycler or cDNA samples. The most frequently adjusted factors are MgCl<sub>2</sub> concentration and annealing temperature. You may prepare three different reactions using MgCl<sub>2</sub> at a concentration of 0.5-3 mM (e.g., 0.5-0.8 mM, 1.5 mM and 3 mM). Optimize the MgCl<sub>2</sub> and/or the annealing temperature on your instrument using the positive control cDNA provided before using your own cDNA.

Troubleshooting Guide

| Problem  | Cause  | Solution   |
|--|--|--|
| No PCR products                                | A PCR component may be missing or degraded.                    | Try to isolate the problematic reagent by replacing it with a fresh one. A checklist is also recommended when assembling reactions.  |
|  | cDNA or MgCl <sub>2</sub> concentration is not optimal.        | Optimize the cDNA and MgCl <sub>2</sub> concentrations.  |
| High background, smearing or nonspecific bands |  | Increase the annealing temperature or decrease the MgCl <sub>2</sub> concentration. Another solution for avoiding high background is to decrease the amount of cDNA template used for amplification. |
| Amplified products are not the correct size    | Contamination with other DNA<br>cDNA quality is not sufficient | Use sterile techniques while performing PCR experiments.<br>Use a different cDNA preparation.  |
|  | Non-optimal PCR conditions                                     | Optimize PCR conditions especially cDNA and MgCl <sub>2</sub> concentrations and annealing temperature.  |
| Poor resolution of products in agarose gel     |  | Use 2% agarose gel and increase run time.  |

**References**

- Oltavi, Z.N., *et al.*, Cell, **74**, 609 (1993)
- Korsmeyer, S.J., Cancer Research (Suppl), **59**, 1693s (1999)
- Agarwal, N. and Metha, K., Biochem. Biophys. Res. Commun., **230**, 251 (1997)
- Aggarwal, S. and Gupta, S., J. Immunol., **160**, 1627 (1998)

†The PCR process is covered by patents owned by Hoffman-LaRoche, Inc. Purchase of this product does not convey a license under these patents.

ya 6/00

Sigma brand products are sold through Sigma -Aldrich, Inc.  
 Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma -Aldrich publications.  
 Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.  
 Please see reverse side of the invoice or packing slip.