CHEMICON now part of Millipore

### **B6-WHITE™ MURINE ES CELL LINE**

			6		
CATALOG NUMBER:	SCR011	QUANTITY:	5 x 10° cells, passage 10		
LOT NUMBER:		CONCENTRATION:	$2.5 \times 10^6$ cells / vial		
BACKGROUND:	The generation of gene-modified mice, created by homologous recombination in embryonic stem (ES) cells, has become a fundamental tool for analyzing gene function. The influence of genetic background on phenotype has been shown to be an important consideration in selection of a mouse model. Furthermore, the time required to achieve congenic status, which can be 18-24 months has further stimulated demand for targeting in pure inbred strains including C57BL/6 (B6) mice.				
DESCRIPTION:	Millipore's B6-White Murine ES cell line is the first commercially available C57BL/6 tyr <sup>c-2J</sup> albino line that allows for rapid coat-color determination of successful chimerism in the C57BL/6 mouse strain. When B6-White ES cells are injected into C57BL/6 blastocysts, chimeric mice are easily identified by their coat color (a mix of black and white patches), while non-chimeric littermates are black. These cells allow for the efficient generation of gene-targeted mice in a pure C57BL/6 genetic background, thus providing more experimental flexibility. Additionally the use of these cells can lower production costs by eliminating the need to maintain an albino blast donor colony, in order to assess chimerism by coat color when targeting in the C57BL/6 strain.				
	C57BL/6J-Tyr <sup>C-2J</sup> /J - available from Jackson Labs, stock number 000058).				
STRAIN:	C57BL/6-tyr <sup></sup> albino mice				
KARYOTYPE:	40,XY [20] Apparently NORMAL Male Mou	se Karyotype			
	Banding Technique: GTL Band Resolution: Good Metaphases Counted: 20 Analyzed: 6 Karyotyped: 6 Interpretation: Cytogenetic anal from mouse cell line SCR011 p male karyotype. No abnormal of	ysis was performed on twenty 11 and all twenty cells demon cells with trisomy 8 and/or trisc	y G-banded metaphase cells strated an apparently normal omy 11 were detected.		
PATHOGEN TESTING:	Cells tested negative for mycop	lasma and other pathogens b	y infectious disease PCR.		
PRESENTATION:	Cells are supplied frozen in basal medium supplemented with 10% FCS and 10% DMSO				
STORAGE/HANDLING:	Place vials in the vapor phase of liquid nitrogen storage immediately upon receipt until it is convenient to proceed to subculture.				

### CHEMICON now part of Millipore

#### B6-WHITE ES CELL MEDIUM:

Mix all ingredients in the top of a 250 mL filter (0.22 µm PES unit, Cat. No. SCGPUO2RE), and filter sterilize. Store at 4°C. Discard unuse d media 7-10 days after preparation.

### *IMPORTANT NOTE:* For successful B6-White ES cell culture, the following medium formulation is highly recommended.

Component	250 mL Final Volume	Millipore Cat. No.	
IMDM, with 25mM HEPES, 3,024mg/L NaHCO <sub>3</sub> & L-Glutamine, w/out Alpha-thioglycerol & beta- mercaptoethanol	185 mL	SLM-063-B	
FCS- <b>ES cell qualified</b>	50 mL	ES-009-B, ES-010-B, ES-011-B	
L. glutamine (100x), 200mM	2.5 mL	TMS-002-C	
Non-Essential Amino Acids (100x)	2.5 mL	TMS-001-C	
Pen/Strep	2.5 mL	TMS-AB2-C	
Na Pyruvate (100x), 100mM	2.5 mL	TMS-005-C	
2-mecaptoethanol (100x), 0.2mM final conc.	5.0 mL	ES-007-E	
ESGRO <sup>®</sup> mLIF Medium Supplement (1000 U/mL final conc.)	25 μL ESG1107 or 250 μL ESG1106	ESG1107 (1 x 10 <sup>7</sup> units) or ESG1106 (1 x 10 <sup>6</sup> units)	

### SUBCULTURING:

B6-White ES cells are fast growing and should be maintained as small dense, but not confluent cultures. The colonies will be tightly packed with phase bright borders. Optimal results are achieved when the cells are maintained as small colonies at high density, fed daily and passaged 1:3 to 1:5 every other day. It is recommended that B6-White ES cells are cultured on a monolayer of mitotically inactivated primary mouse embryonic fibroblast (PMEF) feeder cells in the presence of 1000U/mL ESGRO mLIF Medium Supplement (ESG1106, ESG1107).

### **Culture Notes:**

- Incubator settings: 7.5% CO₂ in humidified air, 37℃
- Replace the media daily (please refer to the recommended media formulation above).
- Passage every other day
- Trypsinize gently using 0.05% Trypsin/EDTA

#### **Microinjection Notes:**

• Inject 8-12 cells/C57BL/6 blastocyst.

### <u>Day 0</u>

Prepare the following culture dishes containing a monolayer of mitotically inactivated PMEF feeder cells:

- 1 x 25 cm<sup>2</sup> flask
- 1 x 75 cm<sup>2</sup> flask
- 3 x 10 cm dishes

### CHEMICON now part of Millipore

### Day 1

Thaw one vial of B6-White ES cells directly into a 25 cm<sup>2</sup> flask containing a confluent layer of inactivated PMEF cells and 5 mL of freshly prepared B6-White ES cell medium.

- Replace the PMEF cell medium with 5.5 mL of B6-White ES cell medium (see above), and allow it to equilibrate in a 37℃ incubator 1 h our before thawing the ES cells.
- Thaw one vial of B6-White ES cells by gently shaking the tube in a 37°C water bath. When the contents of the tube have thawed, spray the vial with ethanol, dry the outside of vial, and aseptically transfer the contents of the vial to the flask.
- Place the flask in a 37℃ incubator overnight.

#### <u>Day 2</u>

Passage the B6-White ES cells to a 75 cm<sup>2</sup> flask, containing a confluent layer of inactivated PMEF cells and 15 mL of pre-warmed B6-White ES cell medium.

- Examine the ES cells under the microscope. Many small phase bright ES cell colonies should be visible.
- Replace the PMEF cell medium in the prepared 75 cm<sup>2</sup> flask with 15 mL of B6-White ES cell medium, and allow the medium to equilibrate in a 37°C incubator for 1 hour before passaging the ES cells.
- Aspirate the medium from the 25 cm<sup>2</sup> flask containing the B6-White ES cells and rinse with 3 mL of PBS.
- Aspirate the PBS and add 1.5 mL of 0.5% Trypsin/EDTA, place the flask in an incubator for 5 minutes or until the ES cells are dissociated.
- Add 5 mL of B6-White ES cell medium and gently titurate the contents of the flask.
- Transfer the cell suspension to the prepared 75 cm<sup>2</sup> flask.
- Place the flask in the 37°C incubator overnight.

#### <u>Day 3</u>

Electroporation.

• Gently trypsinze the B6-White ES cells, as described above, and follow your preferred electroporation protocol.

Plate the electroporated cells on the prepared 10 cm dishes and begin selection 24 hours following electroporation.

APPEARANCE With this line the pup fur may be more pink than white. When the pups are born one can see that the pups have white blotches on their skin, but as they mature the black and white fur seem to mingle and one ends up with a "smoky" pattern with blotches of white. We suggest that all chimeras obtained are bred because germline transmittance by low coat color chimeras is common.

### CHEMICON now part of Millipore

### **REPRESENTATIVE** Germline transmission data generated using wild type B6-White ES cells:

DATA:

Wild type	Blastocysts injected	Live pups (%)	Chimeric pups (%)	Male chimeras (%)	Germline
A	40	14 (35%)	3 (5%)	2 (66%)	1 (50%)
В	78	40 (51%)	3 ( 7%)	3 (100%)	2 (70%)
С	76	32 (42%)	6 (20%)	5 (83%)	2 (40%)

### PRODUCT LICENSE

The purchase of this product conveys to the purchaser the non-transferable right to use the cells for research use only conducted by the purchaser. For commercial purchasers using this product for research, the purchase of this product conveys to the purchaser the non-transferable right to use the purchased cells for a six (6) month evaluation period after the date of purchase for research use only. After the six (6) month evaluation period, the purchaser must contact Millipore. or its designee to obtain a commercial license.

The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for any purpose. The purchaser cannot use this product or its components in whole plants, plant cells or humans.

Commercial purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research.

If the purchaser is not willing to accept the limitations of the limited use license, Millipore is willing to accept return of the products with a full refund. For information on commercial use of this product or obtaining a license to this product for purposes other than research, contact the Millipore Intellectual Property Department.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2008: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.