Instructions for Use

Wright Stain, Modified

Procedure No. WS

Intended Use

Wright Stain is intended for use in staining blood films or bone marrow films. Solutions are for "In Vitro Diagnostic Use." For professional use only. The data obtained from this qualitative procedure, used in both manual and automated staining techniques, helps demonstrate cellular elements in blood and bone marrow films of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of viral and bacterial infections, allergenic reactions, and other pathological states.

Wright Stain is a modified Romanowsky stain intended for differentially staining the cellular elements of blood. When blood films are treated as herein described, the white blood cell nucleus and cytoplasm take on characteristic blue or pink coloration. The purified dyes in the Sigma-Aldrich formulations of Wright Stain eliminate inconsistent staining and yield reproducible lot-to-lot chromogenic responses.

The procedures in this insert describe the use of Wright Stain as a manual dip stain or for use in batch stainers such as the Midas® III-Plus Automated Stainer (Cat. No. 64000H-94), Hemastainer supplied by Geometric Data, and the Fisher Stainmaster supplied by Fisher Scientific.

Reagents

Wright Stain, Modified (Cat. No. WS: WS16-500ML; WS32-1L; WS80-2.5L; WS128-4L)
Wright Stain, modified 0.3% w/v, buffered at pH 6.8 in methanol

Special Materials Required but Not Provided

- Phosphate Buffer (Cat. No. P3288-1VL or P3288-12VL)
A mixture of sodium phosphate and potassium phosphate, 0.0083 M/L, pH 7.2.
- Rinse Solution 2 (Cat. No. RS2-900ML) Ethanol solution, 18%, of wetting agent.
Contains 0.02% sodium azide as preservative.
- Methanol, Acetone Free (Cat. No. M1775-1GA)
- Microscope / Slides / Coverslips

Storage and Stability

Store Wright Stain Solutions at room temperature (18-26°C). Reagent label bears expiration date.

Store Phosphate Buffer, Rinse Solution 2 and Methanol at room temperature (18-26°C).

Store Working Phosphate Buffer at 2-8°C. Warm before use.

Deterioration

Discard Wright Stain solutions if a precipitate develops. Discard the Working Phosphate Buffer if turbidity or visible bacterial growth is present.

Preparation

Wright Stain, Modified, is supplied ready to use, although Wright stain may be diluted if placed on an automated instrument. See instructions which follow.

Prepare Working Phosphate Buffer by diluting contents of one vial Phosphate Buffer pH 7.2 at 25°C to 3.8 liters or 1 gallon with water. Mix well to dissolve.

Methanol is ready to use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Fresh whole blood films or fresh films from blood anticoagulated with EDTA must be used. Prior to preparation of films, blood must be thoroughly mixed at room temperature (18-26 °C). The films should be prepared within 1 hour of blood collection. If not stained the same day, slides should be fixed in absolute methanol and stored in a dust-free container.

Notes

- For greater cellular detail, staining time may be increased. Color (tones of blue or red) may be varied by increasing or decreasing time in deionized water.
- Rapid (15 seconds) staining is not recommended for bone marrow. For those preparations 1-3 minutes in stain and 2-6 minutes in deionized water gives satisfactory results.
- For batch staining, slide racks and dishes such as supplied by Miles Scientific for Tissue-Tek® are recommended since this system allows for vertical placement of slides.
- Color can be varied by increasing or decreasing time in deionized water. Bone marrows should

be stained for at least 90 seconds and buffered for 90 seconds to 3 minutes.

- Staining times outlined in the above procedures have given satisfactory results in our laboratories. Individual preferences may dictate adjustment of these times.
- The manual procedure times may also be used with the Hemastainer provided swing is in the off position.
- Timing may be varied to suit individual preference. (Automated procedures)
- Water rinses should be deionized water. (Automated procedures)
- If deionized water is not neutral pH, we suggest the use of Working Phosphate Buffer.
- Positive control slides should be included in each run.

Procedure

Dip Method (Rapid – Manual)

- Place approximately 50 mL Wright Stain, Modified, in a Coplin jar.
NOTE: KEEP TIGHTLY SEALED WHEN NOT IN USE. Replace when water artifacts appear in red cells or when a precipitate is evident.
- Fill another Coplin jar with deionized water.
- Place thoroughly dried blood film, feather edge DOWN, in Wright Stain, Modified, for approximately 15 seconds.
NOTE: Rapid dipping for 5-10 seconds may reduce water artifacts on films that are not thoroughly dried.
- Remove slides from stain and place in deionized water, feather edge DOWN, for approximately 30 seconds. DO NOT AGITATE SLIDE WHILE IT IS IN DEIONIZED WATER.
- Rinse briefly in running deionized water and air dry thoroughly before evaluation.

Horizontal Staining Method (Manual)

- Place thoroughly dried blood film on an appropriate staining rack.
- Flood slide with 1-2 mL Wright Stain, Modified.
- After 30 seconds, without rinsing off Wright Stain from Step 2, add an equal volume of deionized water and mix thoroughly by gently blowing on slide.
- After 1 minute, thoroughly rinse with deionized water and air dry.

Batch Staining with Hemastainer

- Set timers to provide the following times for each station:

Station 1	30 seconds
Station 2	2 minutes
Station 3	3.5 minutes
Station 4	30 seconds
Station 5	Skip and go to "air dry"

- Prepare stations for staining in the following manner:

Station 1	Absolute Methanol, 500 mL
Station 2	Wright Stain, Modified 350 mL and Absolute Methanol 150 mL
Station 3	500 mL Working Phosphate Buffer
Station 4	3.8 L of deionized water and 100 mL Working Phosphate Buffer
Station 5	Leave empty

- Turn power switch to ON.
- Set Auto-Manual switch to MANUAL.
- Set Right/Left switch to LEFT.
- Set swing switch to ON.
- Set pump switch to AUTO.
- Load basket with thoroughly air-dried blood films.
- Attach basket to hanger and tighten.
- Start process by setting Auto/Manual switch to AUTO.
- When cycle is complete and slides are thoroughly dry, set Auto/Manual switch to MANUAL. The basket will then return to starting position.

Batch Staining with Fisher Stainmaster

Set program as follows:

Event	Station	Reagent	Time (Minutes)
1	1	Absolute Methanol	0.5
2	2	Wright Stain, Modified	1.5
3	3	Working Phosphate Buffer	1.0
4	6	Deionized water	0.3
5	5	Rinse Solution 2	0.7
6	4	Deionized water	0.3
7	Dry	Air	5.0

Batch Staining with Midas® III-Plus Automated Stainer

Set program as follows:

Step	Bath	Reagent	Time (Seconds)
1	1	Absolute Methanol	30
2	2	Wright Stain, Modified	60-90
3	3	Working Phosphate Buffer	60
4	4	Running deionized water	10
5	Dry	Air	3 minutes or until dry

Unused baths may be omitted.

Performance Characteristics

Nuclei will be varying shades of purple. Cytoplasmic staining will be varying shades of blue to light pink. Fine reddish to lilac granules may be present in cytoplasm of some cell types. Basophils will demonstrate dark blue black granules in the cytoplasm. Eosinophils will demonstrate bright orange granules in the cytoplasm. Red blood cells should be pink to orange.¹

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
WS	Wright Stain	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Granules	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

WS16, WS32, WS80, WS128:

H226 Flammable liquid and vapor.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H319 Causes serious eye irritation.

H370 Causes damage to organs (Eyes, Central nervous system).

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Hematology: Principles and Procedures, Sixth Edition, Brown AB, Lea & Febiger, Philadelphia 1993 p100-102

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.SigmaAldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.SigmaAldrich.com/techservice).

Revision History

Rev. 4.0	2014
Rev. 5.0	2016
Rev. 6.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Updated automated stainer reference in Intended Use. Added Warnings and Hazards. Added CH-REP information.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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MDSS CH GmbH
MDSS GmbH