Instructions for Use

Mayer's Hematoxylin Solution

Procedure No. MHS

Intended Use

Mayer's Hematoxylin Solution is commonly used after immunohistochemistry or cytochemistry staining as a nuclear counterstain. It may also be used for standard hematoxylin and eosin (H&E) staining, but it is more commonly used where acid alcohol differentiation, or exposure to alcohol, might destroy the stained cytoplasmic component.³ Mayer's Hematoxylin Solution is formulated without alcohol, and as such will not dissolve out AEC (3-amino-9-ethylcarbazole), alkaline phosphatase/Fast Red chromogen or other soluble colored products. Mayer's Hematoxylin Solutions are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual qualitative procedure is used for the determination of chromatin in nuclei of human specimens. This data can be used as an aid to diagnosis of certain clinical conditions or pathophysiological states and should be reviewed in conjunction with other clinical diagnostic tests or information.

Hematoxylin, a common nuclear stain, is isolated from an extract of logwood.¹ The first successful biologic application of hematoxylin was described by Bohmer in 1865.¹ Mayer introduced his formulation in 1903.² Since then numerous formulations have appeared. Of these, Harris', Gill's, Mayer's and Weigert's have retained popularity. Before hematoxylin can be used as a nuclear stain, it must be oxidized to hematein and combined with a metallic ion (mordant). Most successful mordants have been salts of aluminum or iron.

Generally, hematoxylin solutions are classified as progressive or regressive based on dye concentration. Progressive stains (e.g., Mayer's hematoxylin) have a lower concentration of dye and selectively stain nuclear chromatin without staining cytoplasmic structures. The desired intensity is a function of time. If staining times are excessive, a progressive stain might act similarly to a regressive stain solution. Staining with progressive stains generally requires more time than staining with regressive stains. Regressive stains (e.g. Harris hematoxylin) color all stainable tissue components (nuclear and cytoplasmic) intensely. To arrive at the correct staining response, excess dye must be removed from the tissue section. After sufficient differentiation, a properly destained section will demonstrate nuclear staining, but will not stain cytoplasmic structures.

The final step in hematoxylin staining is the "blueing" of the tissue section. Initially tissue sections are colored either purple or a reddish purple. After exposure to alkaline solutions (warm tap water [if slightly alkaline], dilute ammonia water, Scott's tap water substitute, or lithium carbonate), the tissue section takes on the characteristic blue color of a hematoxylin stained slide.

Reagents

Mayer's Hematoxylin Solution (Cat. No. MHS: MHS16-500ML; MHS32-1L; MHS80-2.5L; MHS128-4L)

Certified hematoxylin (1.0 g/L), C.I. 75290, sodium iodate (0.2 g/L), aluminum ammonium sulfate-12 H₂O (50 g/L), chloral hydrate (50 g/L) and citric acid (1 g/L).

Special Materials Required but Not Provided

- Eosin Y Solution Counterstains:
Alcoholic (Cat. No. HT1101: HT110116-500ML; HT110132-1L; HT110180-2.5L; HT1101128-4L)
Aqueous (Cat. No. HT1102: HT110216-500ML; HT110232-1L; HT110280-2.5L; HT1102128-4L)
OR Alcoholic with Phloxine (Cat. No. HT1103: HT110316-500ML; HT110332-1L; HT110380-2.5L; HT1103128-4L)
- Reagent Alcohol (Cat. No. R8382-1GA) OR Ethanol, 100%
- Scott's Tap Water Substitute Concentrate (Cat. No. S5134-6x100ML)
- Xylene or Xylene Substitute
- Microscope, microscope slides, coverslips, and staining dishes

Storage and Stability

Store reagent at room temperature (18–26°C) protected from light. Reagent is stable until expiration date shown on the label. Do not return used solution to stock bottle.

Deterioration

Discard if staining times becomes excessive or solution turns brown.

Preparation

Filter Mayer's Hematoxylin Solution before each use. Solution is then ready to use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Standard histology texts provide necessary details for specimen collection and storage.^{4,5}

Notes

- The times given in the insert are approximate. Personal preferences will vary and the times can be adjusted to suit personal preferences. Stain solutions which are heavily used will lose their staining powers and the staining times should be lengthened or new solutions should be used.⁶
- Dilute alkaline solutions may be used in place of warm running tap water. This will shorten the time needed for the staining procedure. If using a dilute alkali solution, be sure to wash slides an additional 2–3 minutes in running tap water before proceeding to Eosin staining.
- Some tap water supplies are acidic and unsuitable for use in the "blueing" portion of this procedure. If tap water is acidic, use a dilute alkaline solution.
- Purple or red-brown nuclei are indicative of inadequate "blueing".
- If eosin staining is excessive, nuclear staining may be masked. Proper eosin staining will demonstrate a 3-tone effect. To increase differentiation of eosin, extend time in alcohols or use a first alcohol with a higher water content. The times in the alcohols may be adjusted to obtain the proper degree of Eosin staining.
- Filter working stain solution daily. Rotate alcohols and xylene/xylene substitute daily.
- Adding new stock to depleted working solutions of Mayer's Hematoxylin or Eosin is not recommended.
- Avoid excessive water carry-over into Mayer's Hematoxylin.
- Positive control slides should be included in each run.

Procedure

Procedure 1: Hematoxylin and Eosin Staining

- Prepare a 95% alcohol solution by adding 5 mL deionized water to 95 mL Reagent Alcohol or Ethanol (100%).
- Deparaffinize to water or fix and hydrate frozen sections.
- Stain in Mayer's Hematoxylin Solution for 15 minutes.
- Rinse in warm running tap water for 15 minutes.
- Place in distilled water for 30 seconds.
- If Alcoholic Eosin is to be used, place in Reagent Alcohol, 95% for 30 seconds.
- Place in Eosin Y Solution Counterstain, Alcoholic, Aqueous or Alcoholic with Phloxine for 30 – 60 seconds.
- Dehydrate and clear through 2 changes each of 95% Reagent Alcohol, Reagent Alcohol, and xylene for 2 minutes each.
- Mount with resinous mounting medium.

Procedure 2: Nuclear Counterstain for Special Stains

- Complete individual staining procedure.
- Rinse in deionized water.
- Stain in Mayer's Hematoxylin Solution 1–5 minutes.
- Rinse in running tap water or dilute alkaline solution until nuclei are blue.
- Rinse in deionized water.
- If any portion of the stain is alcohol soluble, mount in aqueous mounting media. If stain is alcohol insoluble, dehydrate in alcohol, clear in xylene or xylene substitute and mount in resinous mounting media.

Performance Characteristics

Expected Results

Nuclear chromatin should be blue. Nucleoli should be visible. Cytoplasm will display various shades of pink to pink-orange (depending upon the counterstain used) and red blood cells will be red.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
MHS	Mayer's Hematoxylin Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

MHS16, MHS32, MHS80, MHS128:

H331: Toxic if inhaled.

P261: Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.

P271: Use only outdoors or in a well-ventilated area.

P304 + P340 + P311: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

P403 + P233: Store in a well-ventilated place. Keep container tightly closed.

P405: Store locked up.

P501: Dispose of contents/ container to an approved waste disposal plant.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer

References

1. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
2. Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175
3. Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p127
4. Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
5. Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
6. Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 5.0	2016
Rev. 6.0	2022
Rev. 7.0	2022
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Added Warnings and Hazards.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271 Darmstadt,
Germany

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, S. 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

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Références

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175.
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

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Riferimenti

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17.
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175.
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p127.
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980.
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992.
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129.

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Referencias

- Conn’s Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

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Referenser

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, omarbetad av Bancroft JD och Gamble, M, Churchill Livingstone, New York, 2002, p129

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Referências

- Conn’s Biological Stains, 10th ed., RW Horobin, JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD, Gamble, M, Churchill Livingstone, New York, 2002, p 127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p 129

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- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p 127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p 129

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- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, p 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

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Referanser

- Conn's Biological Stains, 10. utg., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, s. 17
- Histopathologic Technic and Practical Histochemistry, 3. utg., RD Lillie, Editor. McGraw-Hill, New York, 1965, s. 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, s. 127
- Theory and Practice of Histotechnology, 2. utg., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4. utg., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, s. 129

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Referanslar

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, s. 17
- Histopathologic Technic and Practical Histochemistry, 3rd ed., RD Lillie, Editor. McGraw-Hill, New York, 1965, s. 175
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, s. 127
- Theory and Practice of Histotechnology, 2nd ed., Sheehan DC, Hrapchak BB, Editors, CV Mosby Co, St Louis (MO) 1980
- Laboratory Methods in Histotechnology of the Armed Forces Institute of Pathology, 4th ed., Prophet EB, Mills B, Arrington JB and Sobin LH, Editors, American Registry of Pathology, Washington DC 1992
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, s. 129

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MDSS GmbH