Instructions for Use

Giemsa Stain, Modified

Procedure No. GS-10

Intended Use

Giemsa and May Grünwald solutions are intended for use in staining blood films or bone marrow films. Solutions are for "In Vitro Diagnostic Use." For professional use only. The data obtained from this manual, qualitative procedure demonstrates cellular elements in blood and bone marrow films of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of viral and bacterial infections, allergenic reactions, and other pathological states.

Giemsa stain is a buffered thiazine-eosinate solution designed to provide coloration of blood cells similar to the original product described by Giemsa. It may be used separately or in combination with a May Grünwald Stain, also available from Sigma-Aldrich.

Reagents

Giemsa Stain, Modified (Cat. No. GS: GS500-500ML; GS1L-1L; GS80-2.5L; GS128-4L) Eosin Y, C.I. 45380, 0.15%, Azure B, C.I. 52010, 0.25% in a buffered methanol solution, pH 6.9, with stabilizer

Special Materials Required but Not Provided

- May-Grünwald Stain (Cat. No. MG: MG500-500ML; MG1L-1L; MG80-2.5L; MG128-4L). May-Grünwald Stain, 0.25% w/v, in methanol.
- Phosphate Buffer pH 7.2 at 25°C (Cat. No. P3288-1VL or P3288-12VL) A mixture of sodium phosphate and potassium phosphate, 0.0083 M/L, pH 7.2.
- Microscope / Slides / Coverslips / Staining dishes.

Storage and Stability

Store Giemsa and May-Grünwald solutions at room temperature (18-26°C). Reagent label bears expiration date.

Store Phosphate Buffer at room temperature (18-26°C).

Deterioration

Discard Giemsa and May-Grünwald solutions if a precipitate develops.

Discard the Working Phosphate Buffer if turbidity or visible bacterial growth is present.

Preparation

Giemsa and May-Grünwald solutions are supplied ready to use, although the Giemsa solution may be diluted 1:20 before use in either deionized water or in phosphate buffer solutions.

Prepare Working Phosphate Buffer by diluting contents of one vial Phosphate Buffer pH 7.2 at 25°C to 3.8 liters or 1 gallon with water. Mix well to dissolve.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Fresh whole blood films, bone marrow films, or specimens anticoagulated in EDTA should be used.

Notes

The staining procedures listed in this insert have given satisfactory results in our laboratory. Individual color preferences may necessitate slight adjustments to the times utilized.

Positive control slides should be included in each run.

Procedure

Giemsa May-Grünwald

- Dilute Giemsa Stain 1:20 with deionized water. For bluer coloration, water buffered at pH 7.2 may be used in place of deionized water.
- Place slides in May-Grünwald Stain for 5 minutes.
- Place slides in Working Phosphate Buffer or Trizma® Buffer (20-70 mmol/L). pH 7.2, for 1.5 minutes.
- Place slides in dilute Giemsa solution from Step 1 for 15-20 minutes.
- Rinse slides BRIEFLY in DEIONIZED water.
- Air dry and evaluate.

Standard Giemsa

- Fix slides in methanol 5-7 minutes.
- Air dry.
- Dilute Giemsa Stain 1:20 with deionized water. Color can be varied by diluting in buffer.
- Stain film for 15-60 minutes.
- Rinse in deionized water.
- Air dry and evaluate.

Quick Stain Giemsa

- Place air dried blood film in undiluted Giemsa Stain for 1-2 minutes.
- Place in deionized water for 2-4 minutes depending upon color preference.
- Rinse in deionized water.
- Air dry and evaluate.

Performance Characteristics

Nuclei will be varying shades of purple. Cytoplasmic staining will be varying shades of blue to light pink. Fine reddish to lilac granules may be present in cytoplasm of some cell types. Basophils will demonstrate dark blue black granules in the cytoplasm. Eosinophils will demonstrate bright orange granules in the cytoplasm. Red blood cells should be pink to orange.¹

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
GS	Giemsa Stain	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Basophils	3 of 3	3 of 3	3 of 3	3 of 3
		Eosinophils	3 of 3	3 of 3	3 of 3	3 of 3
		Red Blood Cells	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

GS128:

H225 Highly flammable liquid and vapor.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H319 Causes serious eye irritation.

H370 Causes damage to organs (Eyes, Central nervous system).

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

GS1L, GS500, GS80:

H225 Highly flammable liquid and vapor.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H370 Causes damage to organs (Eyes).

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 Wear protective gloves/ protective clothing.

P301 + P310 + P330 IF SWALLOWED: Immediately call a POISON CENTER/ doctor. Rinse mouth.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/ doctor if you feel unwell.

P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- 1. Hematology: Principles and Procedures, Sixth Edition, Brown AB, Lea & Febiger, Philadelphia 1993 p100-102

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 4.0	2014
Rev. 5.0	2016
Rev. 6.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards. Added CH-REP information.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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MDSS GmbH