


92315 Fluorescent Red 730 Reactive


Fluorescent red 730 reactive is a new fluorescent label for the near infrared area. Extremely long wavelength fluorescence is especially well suited for applications where sample autofluorescence might be critical or penetration of tissue or other matrix is required. Fluorescent red 730 reactive shows strong fluorescence.

Product Description

Formula	C ₄₂ H ₅₁ N ₃ O ₈ S
MW	757.96
Molar absorbance	185.000 l · mol ⁻¹ · cm ⁻¹ (determined in ethanol)
Abs. Max	734 nm (Ethanol)
Emission Max.	750 nm (Ethanol)
Quantity	1 mg



Spectrum



Directions for labelling of proteins with Fluorescent Red 730 reactive

1. To prepare a stock solution of the label, dissolve 1 mg of label (NHS-ester) in 50 µl absolute, amine-free DMF (final concentration: approx. 25 nmol · µl⁻¹).
2. Dissolve the desired amount of protein in bicarbonate buffer (pH 9.0, 50 mM), e.g. 1 mg of avidin in 200 µl buffer. Protein concentrations should typically be 2 mg/ml or higher. For antibodies, dialysis (e.g. two changes of buffer, one hour dialysis for each step) is recommended.
3. Transfer an appropriate volume of the label stock solution to the protein solution dropwise and under stirring. Due to the high reactivity of the NHS ester add an equimolar amount or up to an double excess of label to the protein to obtain a dye to protein ratio (D/P) between 1 and 2. Higher molar excesses of the label can lead to overlabelling of the protein causing a decrease in quantum yield of the conjugate.
4. Incubate the mixture react for one hour at room temperature.
5. Separate the obtained protein conjugate from unreacted free dye using a Sephadex column (Sephadex G25 medium; eluent PBS pH 7.2, 22 mM. Cat. no. 76847). For 1 mg of labelled protein, a column of at least 20 cm length and e.g. 6 mm width is a very good choice. First coloured band is the DY-labeled protein.

Bicarbonate buffer, pH 9.0, 50 mM

Dissolve 2.1 g of NaHCO₃ in 400 ml double distilled water. Adjust the pH to 9.0 by carefully adding small volumes of 1 M HCl or 1 M NaOH while controlling pH with a pH-meter. Add double distilled water up to a final volume of 500 ml.

Sephadex is a registered trademark of GE Healthcare

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

