

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

# **ProductInformation**

## MOPS sodium salt SigmaUltra

Product Number **M5789**Store at Room Temperature

#### **Product Description**

Molecular Formula: C<sub>7</sub>H<sub>14</sub>NO<sub>4</sub>SNa

Molecular Weight: 231.2 CAS Number: 71119-22-7

pK<sub>a</sub>: 7.2 (25 °C)

Effective buffering range: pH 6.5 - 7.9

 $\Delta pK_a/\Delta T$ : -0.015

Synonyms: 3-morpholinopropanesulfonic acid sodium salt, 3-(N-morpholino)propanesulfonic acid sodium salt

Trace elemental analyses have been performed on the SigmaUltra MOPS sodium salt. The Certificate of Analysis provides lot-specific results. SigmaUltra MOPS sodium salt is for applications which require tight control of elemental content.

The zwitterionic buffer MOPS is a structural analog to the Good buffer MES. The Good buffers were developed in the 1960's for general use in biochemistry to meet the following criteria:

- midrange pK<sub>a</sub>
- maximum water solubility and minimum solubility in all other solvents
- · minimal salt effects
- minimal change in pK<sub>a</sub> with temperature,
- · chemical and enzymatic stability,
- minimal absorption in visible or UV range
- reasonable ease of synthesis.<sup>2</sup>

The pK $_a$  of MOPS (7.2) is closer to physiological pH than that of MES (6.1), and thus MOPS may be more suitable as a physiologically relevant buffer.

MOPS buffer has been utilized in the culture of cells in such systems as *E. coli*, *Cryptococcus neoformans*, cultured human keratinocytes, and thermophilic methaogenic bacteria. <sup>3,4,5,6</sup> In protein studies, MOPS has been used in an X-ray crystallographic study of the ADP-binding site of succinyl-CoA synthetase from *E. coli*, in the characterization of the Rieske-type

ferredoxin BphF, and in an electron microscopy analysis of the engineered protein betabellin-15D. <sup>7,8,9</sup>

An investigation of the interaction of various buffers, including MOPS, with plasmid sized DNA by free solution capillary electrophoresis has been reported. A protocol describes the use of MOPS in an electrophoresis buffer for the separation of RNA in agarose gels. A procedure for preparative-scale separation of proteins by displacement chromatography that incorporates MOPS buffer has been published. 12

#### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

## **Preparation Instructions**

This product is soluble in water (231 mg/ml).

### Storage/Stability

Solutions of MOPS are not completely stable when autoclaved in the presence of glucose. <sup>13</sup> Solutions of MOPS turn yellow when autoclaved, indicating that MOPS is unstable to autoclaving.

#### References

- Ellis, K. J. and Morrison, J. F., Buffers of constant ionic strength for studying pH-dependent processes. Methods Enzymol., 87, 405-426 (1982).
- 2. Good, N. E., et al, Hydrogen ion buffers for biological research. Biochemistry, **5(2)**, 467-477 (1966).
- 3. Tucker, D. L., et al., Gene expression profiling of the pH response in *Escherichia coli*. J. Bacteriol., **184(23)**, 6551-6558 (2002).
- Petrou, M. A., and Shanson, D. C., Susceptibility of *Cryptococcus neoformans* by the NCCLS microdilution and Etest methods using five defined media. J. Antimicrob. Chemother., 46(5), 815-818 (2000).

- Sando, G. N., et al., Induction of ceramide glucosyltransferase activity in cultured human keratinocytes. Correlation with culture differentiation. J. Biol. Chem., 271(36), 22044-22051 (1996).
- Foster, M. S., et al., Improved methods for the cultivation of strictly anaerobic, extremely thermophilic methanogens. Biotechniques, 15(6), 996-998, 1000, 1002 (1993).
- 7. Joyce, M. A., et al., ADP-binding site of Escherichia coli succinyl-CoA synthetase revealed by x-ray crystallography. Biochemistry, **39(1)**, 17-25 (2000).
- 8. Couture, M. M., et al., Characterization of BphF, a Rieske-type ferredoxin with a low reduction potential. Biochemistry, **40(1)**, 84-92 (2001).

- Lim, A., et al., Engineering of betabellin-15D: a 64 residue β-sheet protein that forms long narrow multimeric fibrils. Protein Sci., 7(7), 1545-1554 (1998).
- Stellwagen, N. C., et al., DNA and buffers: are there any noninteracting, neutral pH buffers? Anal. Biochem., 287(1), 167-175 (2000).
- Molecular Cloning: A Laboratory Manual, 3rd ed., Sambrook, J. and Russell, D.W., CSHL Press (Cold Spring Harbor, NY: 2001), p. 7.32.
- Narahari, C. R., et al., Displacement chromatography of proteins using a selfsharpening pH front formed by adsorbed buffering species as the displacer. J. Chromatogr. A, 825(2), 115-126 (1998).
- 13. The Merck Index, 12th ed., Entry# 6346.

GCY/RXR 1/08