

Data Sheet

# 12Z Human Endometrial Epithelial Cell Line

Immortalized Cell Line

**SCC443****Pack Size:**  $\geq 1 \times 10^6$  viable cells/vial**Store at:** Liquid nitrogen**FOR RESEARCH USE ONLY****Not for use in diagnostic procedures. Not for Human or Animal Consumption.**

## Background

Endometriosis is a gynecologic condition that is characterized by presence of endometrial tissue outside of the uterus.<sup>1</sup> This condition affects around 6-10% of all women and could lead to pelvic pain and infertility.<sup>2</sup> The limited availability of in vitro models hinders cellular, molecular, and pathophysiologic studies on endometriosis which translates into the poor understanding of the disease and limited available treatments.<sup>3</sup>

Several endometriotic cell lines were developed by researchers from Goethe University Frankfurt in Germany.<sup>3</sup> One of these is the 12Z cell line which was derived from the endometrial peritoneum of a 37-year-old female undergoing laparoscopy. The 12Z cell line was immortalized using SV40 plasmid transfection. The cell line expresses both pan-cytokeratin and N-Cadherin.

One of the major models for the initiation of endometriosis is retrograde menstruation which then allows for endometrial tissues to be deposited in various tissues outside of the endometrium.<sup>3</sup> During this process, the cell's invasive property is thought to be essential. The 12Z immortalized human endometrial epithelial cell line is an established model for endometriosis being derived from an N-Cadherin+ endometriotic tissue, a marker for invasiveness. Moreover, this property of the cell line allows it to be used for various studies related to disease progression, pathogenesis, or therapeutics. Commonly, the 12Z cell line is used in 2D<sup>4</sup> or 3D<sup>5</sup> invasion assays and biochemical assays to study the pathophysiology of endometriosis.

## Source

12Z cell line was derived from endometrial peritoneum biopsy of a 37-year-old female. This cell line was immortalized via transfection of a plasmid containing SV40T.

## Storage and Handling

12Z cells should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting functionality.

## Short Tandem Repeat (STR) Profile


D3S1358: 14, 17	D21S11: 28, 30.2	TH01: 9, 9.3	Penta E: 11, 12
D7S820: 9, 12	D18S51: 13, 15	TPOX: 11, 12	
vWA: 17	D5S818: 12	CSF1PO: 11, 13	
FGA: 22, 23	D13S317: 12, 13	AMEL: X	
D8S1179: 10	D16S539: 13	Penta D: 9, 12	

Cancer cell lines are inherently genetically unstable. Instability may arise in the form of loss of heterozygosity of alleles at one or more genetic sites with increased passages.

## Quality Control Testing

- Each vial contains  $\geq 1 \times 10^6$  viable cells.
- Cells are tested negative for infectious diseases by a Human Essential CLEAR Panel by Charles River Animal Diagnostic Services.
- Cells are verified to be of human origin and negative for inter-species contamination from mouse, rat, Chinese hamster, Golden Syrian hamster, and Non-human Primate (NHP) as assessed by a Contamination Clear panel by Charles River Animal Diagnostic Services.
- Cells are negative for mycoplasma contamination.

## Representative Data



**Figure 1.** Bright-field image of 12Z cells two days after thaw in a T75 flask (A). Cell confluency (82%) was assessed throughout the culture using MilliCell® Digital Cell Imager (Cat. No. MDCI 10000) (B). Cell counting was performed using Scepter™ 3.0 handheld automated cell counter (PHCC360KIT) using 60 mm sensor tips (C). 12Z cells express the invasive marker CDH2 (N-Cadherin) (D), general epithelial marker, pan-cytokeratin (E), and immortalization antigen SV40T (F).

## Protocols

### Thawing Cells

1. Do not thaw the cells until the recommended medium is on hand. Cells can grow on normal tissue cultureware surfaces without any additional coating.  
Cells are thawed and expanded in 12Z Expansion Medium comprised of DMEM medium (Cat. No. SLM-120-B) with 10% FBS (Cat. No. ES-009-B), 2 mM L-Glutamine (Cat. No. TMS-002-C), and 1X penicillin/streptomycin (Cat. No. TMS-AB2-C, optional).
2. Remove the vial of frozen 12Z cells from liquid nitrogen and incubate in a 37 °C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.  
**IMPORTANT:** Do not vortex the cells.
3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
4. In a laminar flow hood, use a 1- or 2-mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
5. Using a 10 mL pipette, slowly add dropwise 9 mL of 12Z Expansion Medium (Step 1 above) to the 15 mL conical tube.

- IMPORTANT:** Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.
6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.  
**IMPORTANT:** Do not vortex the cells.
  7. Centrifuge the tube at 1100 rpm for 5 minutes to pellet the cells.
  8. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
  9. Resuspend the cells in 20 mL of 12Z Expansion Medium.
  10. Transfer the cell mixture to a T75 tissue culture flask.
  11. Incubate the cells at 37 °C in a humidified incubator with 5% CO<sub>2</sub>.

### Subculturing Cells

1. 12Z cells should be passaged at ~80-85% confluency. Do not allow the cells to grow over 85% confluency.
2. Carefully remove the medium from the T75 tissue culture flask containing the 80% confluent layer of 12Z cells. Rinse the flask with 10 mL 1X PBS. Aspirate after the rinse.
3. Apply 5 mL of Accutase™ or Trypsin/EDTA solution and incubate in a 37 °C incubator for 5 minutes.
4. Inspect the flask and ensure the complete detachment of cells by gently tapping the side of the flask with the palm of your hand.
5. Add 8 mL of 12Z Expansion Medium to the plate.
6. Gently rotate the flask to mix the cell suspension. Transfer the dissociated cells to a 15 mL conical tube.
7. Centrifuge the tube at 300 x g for 3-5 minutes to pellet the cells.
8. Discard the supernatant.
9. Add 2-5 mL of 12Z Expansion Medium to the conical tube and resuspend the cells thoroughly.  
**IMPORTANT:** Do not vortex the cells.
10. Count the number of cells using a hemocytometer or a Scepter™ 3.0 handheld automated cell counter using 60 mm sensor tips.
11. Plate the cells to the desired density. Typical split ratio is 1:5 to 1:8. The medium should be replaced every other day.

### Cryopreservation of Cells

12Z cells may be frozen in 12Z Expansion Medium and 10% DMSO using a Nalgene® slow freeze Mr. Frosty® container.

### References

1. Nat. Rev. Endocrinol 2019, 15(11): 666-682.
2. Lancet 2004, 364 (9447): 1789-1799.
3. Am. J. Pathol 2001, 159(5): 1839-1852.
- Reproductive BioMedicine Online 2022, 44(6): 976-990. DOI: <https://doi.org/10.1016/j.rbmo.2022.02.008>.
5. Sci Rep 2021, 11(1): 4115.

### Academic Use Agreement

Subject to local law

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL MILLIPORESIGMA ("SELLER") CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means 12Z Human Endometrial Epithelial Cell Line (SCC443).

"Academic Research" means any internal in vitro research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

Re-engineering or copying the Product

Making derivatives, modifications, or functional equivalents of the Product

Obtaining patents or other intellectual property rights claiming use of the Product

Using the Product in the development, testing, or manufacture of a Commercial Product

Using the Product as a component of a Commercial Product

Reselling or licensing the Product

Using the Product in clinical or therapeutic applications including producing materials for clinical trials

Administering the Product to humans

Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT [licensing@emdmillipore.com](mailto:licensing@emdmillipore.com) AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

## GMO

This product contains genetically modified organisms.

Este producto contiene organismos genéticamente modificados.

Questo prodotto contiene degli organismi geneticamente modificati.

Dieses Produkt enthält genetisch modifizierte Organismen.

Ce produit contient des organismes génétiquement modifiés.

Dit product bevat genetisch gewijzigde organismen.

Tämä tuote sisältää geneettisesti muutettuja organismeja.

Denna produkt innehåller genetiskt ändrade organismer.

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck operates  
as MilliporeSigma in the U.S. and Canada.

Merck Millicell, Scepter and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

00071644 Ver 1.0, Rev 23Mar2023, KK

