

## N-Succinyl-Leu-Leu-Val-Tyr 7-Amido-4-trifluoromethylcoumarin

Product Number **S 4939** Storage Temperature –20 °C

## **Product Description**

 $Molecular\ formula:\ C_{40}H_{50}F_3N_5O_{10}$ 

Mol. wt.: 817.9

N-Succinyl-Leu-Leu-Val-Tyr 7-amido-4-trifluoromethylcoumarin (Suc-LLVY-AFC) is a fluorogenic substrate suitable for the assay of enzymes such as chymotrypsin and proteasome protease.<sup>1</sup>

When Suc-LLVY-AFC is hydrolyzed, the free AFC produced in the reaction can be quantified by fluorometric detection (excitation 400 nm, emission 505nm) or by spectrophotometric detection at 380 nm (extinction coefficient = 12,600 at pH 7.2). When used in an enzyme assay with fluorescence detection, AFC has higher sensitivity than 4-methoxy-2-naphthylamide (MNA).<sup>2</sup>

Chymotrypsin is a serine protease with preferential cleavage at Tyr-, Trp-, Phe-, and Leu- residues. The 26S proteasome is a 2.5 Mda cellular complex that consists of approximately 31 different subunits.<sup>3</sup> It is a self-compartmentalizing protease, which digests polyubiquitinated proteins in the process of ubiquitin-mediated proteolysis.<sup>4</sup>

### **Preparation Instructions**

Prepare stock 20 mM solutions in DMSO. Also soluble in DMF.

# **ProductInformation**

### Storage/Stability

Store at -20 °C. Product is stable for at least one year, if stored as recommended.

Store stock solutions in frozen aliquots at -20 °C. Allow the material to warm to room temperature before use to ensure stability.

#### References

- Orino, E., et al., ATP-dependent reversible association of proteasomes with multiple protein components to form 26S complexes that degrade ubiquitinated proteins in human HL-60 cells. FEBS Lett., 284, 206-210 (1991).
- Johansen, H. T., et al., Colorimetric and fluorimetric microplate assays for legumain and a staining reaction for detection of the enzyme after electrophoresis. Anal. Biochem., 273, 278-283 (1999).
- 3. Voges, D., et al., The 26S proteasome: a molecular machine designed for controlled proteolysis. Annu. Rev. Biochem., **68**, 1015-1068 (1999).
- 4. Kierszenbaum, A. L., The 26S proteasome: ubiquitin-mediated proteolysis in the tunnel. Mol Reprod. Dev., **57**, 109-110 (2000).

JXU/JWM 12/01