Cell Lysis and Recombinant Protein Extraction

CelLytics

RT

CelLytic™ P Plant Cell Lysis Reagent

C 2360 This is a non-ionic detergent-based 50 mL reagent, which offers a convenient method 250 mL for efficient plant cell lysis and protein solubilization. It's a non-denaturing reagent and maintains protein immunoreactivity and biological activity. CelLytic P is efficient, rapid, and ready to use. It contains bicine buffer, which is preferable for many biological activities. Use of CelLytic P enables extraction of proteins from less than one gram to hundreds of grams of fresh or frozen leaves, employing the same short procedure. It has been tested on (but not limited to) four plant models: tobacco, tomato, spinach, and Arabidopsis. S: 23-24/25

Detection of DNA/Protein Interactions



Figure 1. Compatibility with Gel Shift Assay.

Protein extracts were prepared with the CelLytic-P reagent from spinach leaves. A double stranded 32 P labeled CREB oligonucleotide probe was incubated with 28 μ g of the whole cell extract (lanes 2-5) or without whole cell extract (lane 1, free probe). Binding reactions with the extracts were performed in the absence [-] of competitor oligonucleotide (lanes 1-2) or in the presence of x100 or x500 fold excess of unlabeled CREB binding motif oligonucleotide (specific competitor [SP], lane 3 and lane 4, respectively) or in the presence of x100 fold excess of unlabeled oligonucleotide (non specific competitor, [NS] lane 5). Binding reactions were run on a non-denaturing 6% polyacrylamide gel, dried and imaged on X-ray film. The arrows indicate the CREB-DNA complex and the free probe.

CelLytic[™] PN Extraction Kit

CELLYT-PN-1 2-8°C WET ICE	For plant leav This kit is for th and extraction leaves. Nuclei of few grams to h	es e rapid isolation of nuclei of functional nuclear proteir r nuclear proteins can be exi undreds of grams of fresh a	1 kit is from plant tracted from a ind frozen		
	leaves. The nuclear pro DNA-protein in footprinting an similar techniqu a source for ch provides a deta protein extracti tomato, spinacl 1 kit sufficient leaves Components: Nuclei Isolation Percoll, Sucrose 2.3 M, TRITON [®] X-100 Extraction Buffe Nuclei PURE Stt Filter Mesh 100 R: 20/22-37/38-4	e detection of ay, DNase-I t assay and Iso be used as A, etc. The kit tion and tobacco, n or frozen			
	Cyto Nuc				
	- 10 ¹	 RNA Polymerase II -205 -119 - 98 - 52 			
Figure 1. Detection of RNA polymerase II in Tomato Nuclear versus					

Cytoplasmic extracts, prepared with CelLytic-PN Kit. The Extracts were run on SDS-PAGE and blot-hybridized to anti-RNA Polymerase II antibody

ſ

Cell Lysis and Recombinant Protein Extraction CelLytics

Compatibility with Gel Shift Assay



Figure 2: Gel Shift analysis of protein-DNA complex formed between CREB DNA binding site and CREB protein in Tobacco Nuclear extracts.

Protein extracts were prepared with the CelLytic-PN kit from nuclei isolated at three levels of purity (high, semi-pure, crude). A double stranded ³²P labeled CREB oligonucleotide probe was incubated with 1 µg of cytoplasmic or nuclear extract. Binding reactions with the nuclear extracts were performed in the absence [–] of competitor oligonucleotide (lanes 1-5) or in the presence of x100 or x500 fold excess of unlabeled CREB binding motif oligonucleotide (specific competitor [SP], lane 6 and lane 7, respectively) or in the presence of x100 fold excess of unlabeled oligonucleotide (non specific competitor, [NS] lane 8).

Lane 1: Free probe without extract.

Lane 2: Cytoplasmic fraction present in the supernatant of "High" level of purity. Lane 3: Nuclear proteins isolated by "High" level of purity.

Lane 4: Nuclear proteins isolated by "Semi-pure" level of purity.

Lanes 5-8: Nuclear proteins isolated by "Crude" level of purity.

Binding reactions were run on a non-denaturing 6% polyacrylamide gel, dried and imaged on film. The arrows indicate the CREB-DNA complex and the free probe.

CelLytic™ B Bacterial Cell Lysis Reagent

B 3553 Standard strength

RT

A gentle, highly efficient reagent for the 500 mL extraction of proteins from bacteria (*E. coli*). A proprietary formulation of a non-ionic detergent in 20 mM Trizma[®]-HCI (pH 7.5).

Treatment of bacterial cells with CelLytic B results in rapid extraction of proteins that are suitable for affinity purification and analysis.

50 mL

Features and Benefits

- Scalable for 1 to 25 grams of bacterial cell paste
- No interference with downstream applications such as affinity chromatography, IP, and Western blotting
- Compatible with protease inhibitors, inhibitor cocktails, chaeotropes, salts, chelating agents and reducing agents

CelLytic[™] B-II Bacterial Cell Lysis Reagent

B 3678 2× strength

RT

A gentle, highly efficient reagent for the 250 mL extraction of proteins from bacteria (E. coli). A proprietary formulation of a non-ionic detergent in 20 mM Trizma[®]-HCI (pH 7.5).

Treatment of bacterial cells with CelLytic B results in rapid extraction of proteins that are suitable for affinity purification and analysis.

Features and Benefits

- Scalable for 1 to 25 grams of bacterial cell pasteNo interference with downstream applications such as
- affinity chromatography, IP, and Western blotting • Compatible with protease inhibitors, inhibitor cocktails,
- chaeotropes, salts, chelating agents and reducing agents

CelLytic™ B Plus Kit

The CelLytic[™] B Plus Kit contains all of the reagents and chemicals necessary to lyse both Gram negative and Gram positive bacteria. This kit also includes protease inhibitors to

help prevent the proteolytic breakdown of proteins. This kit is compatible with FLAG[®], HIS-Select[™], and glutathione S-transferase based affinity chromatography protein purification systems. The detergent included in the CelLytic B lysis Reagent can be removed from the protein by dialysis or ammonium sulfate precipitation, if necessary. The final purity of the recombinant product is usually higher than that obtained from traditional extraction methods. This is because the crude protein solution containing the non-ionic detergent eliminates much of the non-specific binding that occurs during

chromatography. The mild detergent does not interfere with many enzyme assays or protein assays.

Components:

CelLytic™ B Bacterial Cell Lysis Reagent Lysozyme from chicken egg white Benzonase Protease Inhibitor Cocktail for use in purification of (histidine)-

tagged proteins

CB0050 CelLytic B Plus Kit Introductory Size 1 kit

-20-0°C For bacterial lysis

WET ICE This introductory size kit can lyse up to 5 grams of fresh or frozen cell paste.

CB0500 CelLytic B Plus Kit

- -20-0°C For bacterial lysis
- WET ICE Fifty grams of fresh or frozen cell paste can be lysed with this kit.

CelLytic[™] IB Inclusion Body Solubilization Reagent

C 5236 CelLytic IB was designed to solubilize protein aggregates called inclusion bo

e 25 mL

1 kit

protein aggregates called inclusion bodies, 100 mL
 which are formed in bacteria when
 recombinant proteins are overexpressed. CelLytic IB was
 formulated to solubilize the protein of interest for immediate

analysis of protein content or refolding procedures. R: 22-36/38 S: 26-36 ĥ

ein

0

www.sigma-aldrich.com

Cell Lysis and Recombinant Protein Extraction

CelLytics

(Continuation of)

CelLytic[™] IB Inclusion Body Solubilization Reagent



Figure 1. Solubilization of Streptavidin Inclusion Bodies. Samples of solubilized inclusion body protein were assayed using BCA Reagent (BCA-1). 50 µl of each sample were incubated with BCA reagent at 60 °C for 15 min. The samples were then cooled to room temperature and assayed at 562 nm. The data above has been standardized to the protein recovery of CelLytic IB

Protease Inhibitors

Protease Inhibitor Cocktail for plant cell and tissue extracts

P 9599	DMSO solution	1 mL
–20°C	A mixture of protease inhibitors with broad	5 mL
	specificity for the inhibition of serine,	
	cysteine, aspartic, and metalloproteases, and am	inopeptidases.
	Contains 4-(2-aminoethyl)benzenesulfonyl fluori	de (AEBSF),
	bestatin, pepstatinA, E-64, leupeptin, and 1,10-p	henanthroline.
	One ml is recommended for the inhibition of pr	oteases
	extracted from 30 g of plant tissue in a total volu	ime of 100 ml.
	Tested for inhibition in extracts from Phaseolus v	<i>ulgaris</i> (kidney
	bean), Pisum sativum (pea), Triticum aestivum (v	vheat),
	Nicotiana tobaccum (tobacco), and Arabidopsis	thaliana
	(arabidopsis).	

Protease Inhibitor Cocktail for use with bacterial cell extracts

P 8465 Lyophilized powder

5 mL -20°C A mixture of protease inhibitors with broad 25 mL

specificity for the inhibition of serine, cysteine, aspartic and metallo-proteases, and aminopeptidases. Contains 4-(2-aminoethyl)benzenesulfonyl fluoride (AEBSF), pepstatin A, E-64, bestatin, and sodium EDTA. Five ml is recommended for the inhibition of proteases extracted from 20 g of Escherichia coli. Supplied with a vial of DMSO.

Protease Inhibitor Cocktail for General Use

P 2714 Lyophilized powder

-0°C

1 each A mixture of water-soluble protease inhibitors with broad specificity for the inhibition of serine, cysteine, aspartic, and metallo-proteases. Contains 4-(2aminoethyl)benzenesulfonyl fluoride (AEBSF), E-64, bestatin, leupeptin, aprotinin, and sodium EDTA.

One bottle makes 100 ml of cocktail. One ml is recommended for the inhibition of proteases equivalent to 1 mg of USP pancreatin.

Protein Extraction for Proteomics

Chloroplast Isolation Kit

CP-ISO The chloroplast isolation kit provides an 1 kit 2-8°C efficient procedure for isolating in-tact chloroplasts from plant leaves. In-tact chlorplasts serve as the ٠ WET ICE best starting material for studies of chlorplast processes such as carbon assimilation, electron flow and phosphorylation,

NEW metabolic transport, and protein targeting.

Kit Components

Name	Amount
Chloroplast Isolation Buffer 5x (CIB)	500 ml
Percoll	100 ml
Bovine Serum Albumin (BSA)	3 д
Filter Mesh 100	4 each



Ū

c o m

www.sigma-aldrich.