# **Product Information**

#### sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

## Cholesterol Esterase from bovine pancreas

Catalog Number C3766 Storage Temperature –20 °C

CAS RN 9026-00-0 EC 3.1.1.13 Synonyms: Bile salt activated lipase, sterol esterase, carboxyl ester lipase, steryl-ester acylhydrolase

#### **Product Description**

Excess cholesterol is stored intracellularly as cholesterol esters. Cholesterol esterase (CE) is a reversible enzyme that can hydrolyze or synthesize fatty acid esters of cholesterol and other sterols. Hydrolysis of water insoluble long chain fatty acid esters requires bile salt activation. Hydrolysis of water soluble esters of short chain fatty acids and lysophospholipids does not require activation by bile salts.<sup>1</sup> Cholesterol esterase catalyzes the following reaction:

Cholesterol	Cholesterol	Cholesterol
esters	esterase	+ Fatty acid

While found primarily in the pancreas and pancreatic fluid, it occurs in other tissues as well. In the bovine adrenal cortex, this reaction is one of the rate limiting steps in steroidogenesis, involving the release of cholesterol from cytoplasmic cholesterol esters.<sup>2</sup> Cholesterol esterase bound to membrane-associated heparin on brush border membranes aids in the transport of cholesterol and free fatty acid across the membrane.<sup>3</sup> This enzyme is widely used in the determination of serum cholesterol in clinical laboratories.<sup>4</sup>

Cholesterol esterse is a monomeric glycoprotein containing fucose residues.<sup>5</sup> The interaction of cholic acid with pancreatic CE results in the formation of a 400 kDa hexamer.<sup>6</sup> CE is a highly conserved protein of the  $\alpha/\beta$  hydrolase family based on the nucleotide sequence from various species.<sup>1</sup>

Molecular mass:<sup>7,8</sup> 84 kDa (SDS-PAGE)

pH optimum:<sup>9,10</sup> 7.0 (hydrolysis) 6.2 (esterification)

pH Range:<sup>7</sup> 6-8

Temperature optimum:<sup>9</sup> 37 °C

Substrates:

cholesteryl oleate<sup>12</sup> 4-nitrophenyl butyrate<sup>15</sup> triacylglycerol<sup>2</sup> 2-(diethylamino)-4*H*-3,1-benzoxazin-4-one<sup>14</sup> 2-(diethylamino)-4*H*-thieno[2,3-d][1,3]oxazin-4-one<sup>14</sup>

K<sub>M</sub>: cholesteryl oleate<sup>12</sup> 0.011 mM (in cholesteryl oleate vesicles) 0.099 mM (in liquid crystals)

4-nitrophenyl butyrate<sup>13</sup> 0.052 mM

Activators: bisphenol A diglycidyl ether<sup>11</sup> ethanol<sup>13</sup> phosphatidylcholine<sup>13</sup> methanol<sup>13</sup> phosphatidylethanolamine<sup>12</sup> *n*-butanol<sup>13</sup> cAMP-dependent protein kinase, type II<sup>8,12</sup> sodium taurocholic acid<sup>5</sup>

Inhibitors: bisphenol A methacrylate<sup>11</sup> Hg<sup>2+,8</sup> diisopropylfluorophosphate<sup>8</sup> PMSF<sup>6</sup> enolase<sup>15</sup> sodium fluoride<sup>8</sup> phosphatidic acid<sup>12</sup> phosphatidylcholine<sup>12</sup> *p*-chloromercuribenzoate<sup>6</sup> phosphatidylserine<sup>12</sup>

K<sub>i</sub>: bisphenol A methacrylate<sup>11</sup>
0.0031 mM (pH 7.0, 22 °C)
2-(diethylamino)-4*H*-3,1-benzoxazin-4-one<sup>14</sup>
0.053 mM (pH 7.0, 25 °C)

This product (Catalog Number C3766) is partially purified from bovine pancreas and is supplied as an off-white to tan lyophilized powder containing 30–65% protein (biuret), potassium phosphate, and TRITON<sup>®</sup> X-100.

Specific activity: ≥200 units/g protein

Unit definition: one unit will hydrolyze 1.0  $\mu$ mole of cholesteryl oleate to cholesterol and oleic acid per minute at pH 7.0 at 37 °C in the presence of taurocholate.

Cholesterol esterase is assayed spectrophotometrically in a 3.0 ml reaction mixture containing 287 mM potassium phosphate, pH 7.0, 0.25% (w/v) taurocholic acid, 0.25% (w/v) cholic acid, 4–6 units peroxidase, 1.4 mM cholesteryl oleate, 1.7% (v/v) polyoxyethylene 9-lauryl ether, 0.14% (w/v) NaCl, 0.083% (w/v)phenol, 0.03% (w/v) 4-aminoantipyrine, 1–1.5 units cholesterol oxidase, and 0.013–0.143 unit cholesterol esterase.

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

### **Preparation Instructions**

Cholesterol esterase is soluble in 0.4 M potassium phosphate, pH 7.0 (1 mg/ml).

### Storage/Stability

Store the product at -20 °C. When stored at -20 °C, the product retains activity for at least one year.

#### References

- Hui, D.Y., and Howles, P.N., Carboxyl ester lipase: structure-function relationship and physiological role in lipoportein metabolism and atherosclerosis. J. Lipid Res., 43, 2017-30 (2002).
- Cook, K.G., *et al.*, Hormone-sensitive cholesterol ester hydrolase of bovine adrenal cortex: identification of the enzyme protein. FEBS Lett., **132**, 10-14 (1981).
- Lopez-Candales, A., *et al.*, Cholesterol transport function of pancreatic cholesterol esterase: directed sterol uptake and esterification in enterocytes. Biochemistry, **32**, 12085-89 (1993).
- Allain, C.C., *et al.*, Enzymatic determination of total serum cholesterol. Clin. Chem., **20**, 470-75 (1974).

- Tanaka, H., *et al.*, Purification and characterization of bovine pancreatic bile salt-activated lipase. J. Biochem. (Tokyo), **125**, 883-90 (1999).
- Hyun, J., *et al.*, Cholesterol esterase a polymeric enzyme. Biochem. Biophys. Res. Comm., **44**, 819-25 (1971).
- Cook, K.G., *et al.*, Cytosolic cholesterol ester hydrolase from bovine corpus luteum. Its purification, identification, and relationship to hormone-sensitive lipase. Biochim. Biophys. Acta, **752**, 46-53 (1983).
- 8. Cook, K.G., The identity of the cholesteryl ester hydrolase of bovine corpus luteum. Biochem. Soc. Trans., **11**, 703 (1983).
- Malik, V., *et al.*, Determination of total cholesterol in serum by cholesterol esterase and cholesterol oxidase immobilized and co-immobilized on to arylamine glass. Biotechnol. Appl. Biochem., **35**, 191-97 (2002).
- Vahouny, G.V., and Treadwell, C.R., Enzymatic synthesis and hydrolysis of cholesterol esters. Methods Biochem. Anal., 16, 219-72 (1968).
- 11. Smith, R.E., *et al.*, Aromatic dental monomers affect the activity of cholesterol esterase. Biochim. Biophys. Acta, **1550**, 100-06 (2001).
- Tsujita, T., and Okuda, H., Hydrolysis of cholesteryl oleate liquid crystals by hormone-sensitive lipase. J. Biochem. (Tokyo), **113**, 264-69 (1995).
- Stout, J.S., *et a*l., Acylenzyme mechanism and solvent isotope effects for cholesterol esterasecatalyzed hydrolysis of *p*-nitrophenyl butyrate. Biochim. Biophys. Acta, **837**, 6-12 (1985).
- Pietsch, M., and Gutschow, M., Alternate substrate inhibition of cholesterol esterase by thieno[2,3d][1,3]oxazin-4-ones. J. Biol. Chem., 277, 24006-13 (2002).
- Shand, J.H., and West, D.W., Inhibition of neutral cholesteryl ester hydrolase by the glycolytic enzyme enolase. Is this a secondary function of enolase? Lipids, **30**, 763-70 (1995).
- Ikeda, I., *et al.*, Cholesterol esterase bound to intestinal brush border membranes does not accelerate incorporation of micellar cholesterol into absorptive cells. Biosci. Biotechnol. Biochem., **67**, 2381-87 (2003).

TRITON is a registered trademark of Union Carbide, Corp

KAD, RBG, JWM, MAM 12/07-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.