sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Fetuin from fetal bovine serum

BioReagent, suitable for cell culture

Catalog Number **F3385** Storage Temperature 2-8 °C

Product Description

Molecular weight: 48.4 kDa¹ CAS Number: 9014-81-7 pl: 3.3^1 ; $3.2 - 3.8^2$ Extinction Coefficient: $E^{1\%}$ = 4.1, 4.5 (278 nm)^{1,3} Synonym: Fetuin Type III.

This product is cell culture tested (0.8 g/L) and is appropriate for use in cell culture applications. Fetuin is a protein component of serum which contributes to the attachment and spreading of cells in culture medium. Researchers often add serum for this purpose. Fetuin has been used at a concentration of 500 μ g/ml to supplement serum-free F12 medium (along with insulin, transferrin and 2-mercaptoethanol) in culture of embryonal carcinoma cells.⁴

Fetuin can be used as a glycoprotein standard for carbohydrate structure in a glycoprotein. The structural elucidation of the carbohydrate portion of fetuin has been published.⁵⁻⁹ The cDNA and amino acid sequences has been reported.⁵ Fetuin contains 12 Cys residues, all involved in disulfide bridges.¹⁰

The composition of bovine fetuin (weight %) is polypeptide 74%, hexose 8.3%, hexosamines 5.5%, and sialic acid 8.7%.¹ This product usually assays at 5-8% total sialic acid. Nearly all of this sialic acid is bound. The amount of free sialic acid is typically less than 0.3%. A method for determining free sialic acid content has been published.¹¹

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

This product is soluble in water (1 mg/ml).

Storage/Stability

Solutions of fetuin are stable for approximately 2 weeks when stored at 2-8 °C and for several months when aliquoted and stored at -20 °C. This assumption is based on the observation that cells are still active in serum, which has been stored under the above conditions.

References

- Spiro, R. G., Studies on fetuin, a glycoprotein of fetal serum. I. Isolation, chemical composition and physiocochemical properties. *J. Biol. Chem.*, 235, 2860-2869 (1960).
- 2. Kithier, K., and Poulik, M. D., Comparative studies of bovine-fetoprotein and fetuin. *Biochim. Biophys. Acta*, **278(3)**, 505-516 (1972).
- 3. Murray, A. C., et al., Circular dichroism studies on native fetuin and some of its derivatives. *Biochim. Biophys. Acta*, **175(2)**, 331-338 (1969).
- Bottenstein, J. et al., The Growth of Cells in Serum-Free Hormone-Supplemented Media. *Methods Enzymol.*, 58, 94-109 (1979).
- Dziegielewska, K. M., et al., The complete cDNA and amino acid sequence of bovine fetuin. Its homology with alpha 2HS glycoprotein and relation to other members of the cystatin superfamily. *J. Biol. Chem.*, **265(8)**, 4354-4357 (1990).
- Yet, M. G., The covalent structure of individual N-linked glycopeptides from ovomucoid and asialofetuin. *J. Biol. Chem.*, **263(1)**. 111-117 (1988).
- Takasaki, S. and Kobata, A., Asparagine-linked sugar chains of fetuin: Occurrence of tetrasialyl triantennary sugar chains containing the gal α (1→3)GlcNAc sequence. *Biochemistry*, **25(19)**, 5709-5715 (1986).
- Spiro, R. G., and Bhoyroo, V. D., Structure of the o-glycosidically linked carbohydrate units of fetuin. *J. Biol. Chem.*, **249(18)**, 5704-5717 (1974).

- Rafelson, M. E., Jr., et al., Neuraminidase (sialidase) from influenza virus. Methods Enzymol., VIII, 677-680 (1966).
- Chin, C. C., and Wold, F., The use of tributylphosphine and 4-(aminosulfonyl)-7-fluoro-2,1,3-benzoxadiazole in the study of protein sulfhydryls and disulfides. *Anal. Biochem.*, **214(1)**, 128-134 (1993).
- 11. Aminoff, D., Methods for the quantitative estimation of n-acetylneuraminic acid and their application to hydrolysates of sialomucoids. *Biochem. J.*, **81**, 384-392 (1961).

PHC 03/13-1