

ProductInformation

CHELATING RESIN Sigma Prod. No. C7901 Chelex 100

PHYSICAL PROPERTIES:

Appearance: White moist powder (swollen beads) Analytical grade resin Capacity of Sodium Form: 0.7 meq/ml

STRUCTURE:

A styrene-divinylbenzene copolymer containing paired iminodiacetate ions which act as chelating groups in binding polyvalent metal ions. It is considered a weakly acidic resin.

USAGE:

Product is autoclavable in sodium form. It has a maximum operating temperature of 75°C and is not soluble in water.

Binding is a function of pH. Absorption is very low below pH 2 and increases sharply from pH 2 to 4. It reaches a maximum above pH 4. Optimum binding for many divalent cations is at pH 6.5 or higher.¹ Its selectivity for divalent over monovalent ions is approximately 5000 to 1, and it has a very strong attraction for transition metals, even in highly concentrated salt solutions. Actual selectivity values for any particular system depend on the pH, ionic strength, and presence of other complex-forming species. Metals can be removed using either the batch or column technique, although the column technique is more efficient. With 50-100 mesh, rapid flow rates are obtained and large volumes of solution can be processed in very little time.

Resin can be regenerated by washing in the following sequence: 2 bed volumes of 1 N HCl, 5 bed volumes water, 2 bed volumes of 1 N NaOH, 5 bed volumes water.

APPLICATIONS:

Removal of metals from enzyme solutions² Removal of metals from cell suspensions³ Calcium removal from dinucleotides⁴ Calcium removal from calmodulin and buffer⁵ Removal of calcium from erythrocyte lysates⁶ Reducing calcium and magnesium concentrations in tissue culture medium⁷

CHELATING RESIN Sigma Prod. No. C7901 Chelex 100

APPLICATIONS: (continued)

Effect of pH and ionic strength on chelating properties⁸ Stability of metal complexes⁹ Trace metal studies¹ As a medium for extraction of DNA from forensic-type samples¹⁰ Removal of metal ion from guinea pig complement¹¹

REFERENCES:

- 1. Figura, P. and McDuffie, *B., Anal. Chem.*, vol. 49, 1950 (1977).
- 2. Dunn, M.F., et al., *Biochemistry*, vol. 19, 718 (1980).
- 3. Borson, W.F., et al., *Biochemistry*, vol. 14, 2275 (1975).
- 4. Reinhardt, C.G. and Krugh, T.R., *Biochemistry*, vol. 17, 4845 (1978).
- 5. Crouch, T.H. and Klee, C.B., *Biochemistry*, vol. 19, 3692 (1980).
- 6. Raymond, F.A., and Weinshilboum, R.M., *Clin. Chem. Acta*, vol. 58, 185 (1975).
- 7. Brennan, J.K., et al, *In Vitro*, vol. 11, 354 (1975).
- 8. Luttrell, G.H., More, C. and Kenner, C.T., *Anal. Chem.*, vol. 43, 1370 (1971).
- 9. Chaberek, S. and Martell, A.E., J. Am. Chem. Soc., vol. 74, 5052 (1952).
- 10. Walsh, P.S., Metzger, D.A. and Higuchi, R., *BioTechniques*, vol. 10, 506 (1991).
- 11. Amiraian, K., McKinney, J.A. and Duchna, I., *Immunology*, vol. 26, 1135 (1974).

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.