Instructions for Use

Acid Fast Stain

Procedure No. HT80

Intended Use

Acid Fast Stain is for the histological demonstration of acid fast bacteria in tissue sections and smears. Acid Fast Stain reagents are for "In Vitro Diagnostic Use". For professional use only. This manual, qualitative procedure identifies acid fast bacteria in tissue sections and smears of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of acid fast bacterial infections.

Carbol-fuchsin acid fast stain is used to demonstrate acid fast bacteria in tissue sections and on smears. This procedure uses pararosaniline dye in a phenol solution. Phenol appears to combine the pararosaniline dye within the acid fast bacilli.¹ Dimethyl sulfoxide (DMSO) is added to the solution to accelerate the infiltration of dye. When stain is applied, all cells including hard to stain acid fast, are colored red. The next step is to decolorize all tissue components except the acid fast bacilli. This procedure incorporates the decolorizer with the counterstain in the malachite green solution.² It eliminates problems of over decolorization and inconsistent counterstaining. Included is an acid fast stain technique for tissue paraffin sections using rapid staining in a microwave oven.³⁻⁵

Reagents

Carbol-Fuchsin Solution (Cat. No. HT8018-250ML; HT80116-500ML)

Pararosaniline (certified), C.I. 42500, 0.85%; phenol, 2.5%; glycerol, 5%; DMSO, 5%; and ethyl alcohol, 5%, in deionized water.

Malachite Green Solution (Cat. No. HT8028-250ML; HT80216-500ML)

Malachite green oxalate (certified), C.I. 42000, 1.5%; acetic acid, 10%; and glycerol, 17%, in deionized water. Malachite Green Solution contains both counterstain and differentiator.

Special Materials Required but Not Provided

- Positive control slides, such as Acid Fast TISSUE-TROLTM control slides, Catalog No. TTR001, should be included in each run.
- Xylene or xylene substitute

Optional Reagents/Instruments

- Absolute Alcohol
- Coplin jars
- Microwave oven

Storage and Stability

Store Carbol-Fuchsin Solution and Malachite Green Solution at room temperature (18–26°C). Reagents are stable until expiration dates shown on the labels.

Preparation

Carbol-Fuchsin Solution and Malachite Green Solution are ready to use. Malachite Green Solution should be used once and properly discarded.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Tissue

Any well-fixed paraffin section cut 5–6 microns or frozen section.

Smears

Smear a direct or concentrated specimen on clean unused slide. Air dry and fix smears by flaming or a 30-second immersion in absolute alcohol.

Notes

It is recommended that gloves be worn when performing these procedures.

If methylene blue is used as a counterstain, do not over stain with methylene blue. Overstaining with methylene blue can mask any organism present. If the slides are overstained, the slides should be placed in acid alcohol to remove the methylene blue counterstain. Wash in water and then repeat the counterstaining step.⁶

DO NOT dehydrate through graded alcohols after staining. Alcohol dehydration may result in over-differentiation of tissue sections.

Procedure

Tissue Paraffin Sections – Manual

- Deparaffinize slides and hydrate to deionized water.
- Stain in Carbol-Fuchsin Solution for **30–45 minutes**.
- Wash in deionized water to remove excess stain.
- Differentiate and counterstain in Malachite Green Solution for **2 minutes**.
- Rinse in tap water and let air dry, or blot completely dry.
- Dip in xylene and coverslip.

Microwave Procedure

- Deparaffinize slides and hydrate to deionized water.
- Place slides in **40 mL** Carbol-Fuchsin Solution contained in a plastic Coplin jar. Loosely cover jar with lid or use vented lids.
- Microwave on **400 watts** for **30 seconds**. Stir with beral pipette or applicator stick. Let incubate for **10 minutes**.
- Wash in deionized water to remove excess stain.
- Place slides in ^{40 mL} Malachite Green Solution for **1 minute at room temperature**.
- Rinse in tap water and let air dry, or blot completely dry.
- Dip in xylene and coverslip.

Smears and Frozen Sections – Manual

- Place slide in Carbol-Fuchsin Solution for **3 minutes**.
- Wash in tap water.
- Counterstain and differentiate in Malachite Green Solution for **1–2 minutes**.
- Rinse in tap water and let air dry, or blot completely dry.
- Dip in xylene and coverslip.

Performance Characteristics

Acid fast bacilli	Red
Background	Green

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT801	Carbol-Fuchsin Solution	Acid Fast Bacilli	3 of 3	3 of 3	3 of 3	3 of 3
		Background	3 of 3	3 of 3	3 of 3	3 of 3
HT802	Malachite Green Solution	Acid Fast Bacilli	3 of 3	3 of 3	3 of 3	3 of 3
		Background	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT80116, HT8018:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

H412 Harmful to aquatic life with long lasting effects.

P202 Do not handle until all safety precautions have been read and understood.

P264 Wash skin thoroughly after handling.

P273 Avoid release to the environment.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

HT80216:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H412 Harmful to aquatic life with long lasting effects.

P273 Avoid release to the environment.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

HT8028:

- H226 Flammable liquid and vapor.
- H315 Causes skin irritation.
- H319 Causes serious eye irritation.
- H412 Harmful to aquatic life with long lasting effects.
- P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P233 Keep container tightly closed.
- P240 Ground and bond container and receiving equipment.
- P273 Avoid release to the environment.
- P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
2. Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Bejamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Jornal of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
6. Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 4.0	2014
Rev. 5.0	2016
Rev. 6.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Updated reference 2. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards. Added CH-REP information.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Bejamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Jorنال of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986
- Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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Références

- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Benjamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Jorنال of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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Bibliografia

- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980.
- Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Benjamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Journal of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554.
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986.
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986.
- Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183.

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Referencias

1. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
2. Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Bejamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Jorنال of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
6. Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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Referencer

1. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
2. Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Benjamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Jorنال of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
6. Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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Referências

- Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979
- Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 151
- Luna, Lee G., Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 106

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- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Bejamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Journal of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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Reference

1. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
2. Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Benjamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Jorنال of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986
6. Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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Referanslar

1. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
2. Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Bejamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Journal of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
6. Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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Referenties

- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Glenn E. Pottz, M.S., James H. Rampey, M.S., Furmandean Benjamin, M.T., A Method for Rapid Staining of Acid-Fast Bacteria in Smears and Sections of Tissue, American Journal of Clinical Pathology, Volume 42, Issue 5, November 1964, Pages 552-554
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation, J Pathol 148: 183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
- Histotechnology: A Self-Instructional Text, 2nd Edition by Carson FL, ASCP Press, Chicago (IL), 1997, p 183

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MDSS GmbH