


3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

# **ProductInformation**

# Ionomycin from Streptomyces conglobatus

Catalog Number **I9657** Storage Temperature -20 °C

CAS RN: 56092-81-0



Molecular formula: C<sub>41</sub>H<sub>72</sub>O<sub>9</sub> Molecular weight: 709.00

### **Product Description**

Ionomycin is a polyether antibiotic produced by *Streptomyces conglobatus*, ATCC 31005. 1 It is highly specific for divalent cations.

Ion selectivity is as follows: Ca<sup>2+</sup> > Mg<sup>2+</sup> >> Sr<sup>2+</sup> = Ba<sup>2+</sup> Binding of Sr<sup>2+</sup> and Ba<sup>2+</sup> is insignificant and binding to monovalent cations or rubidium is negligible. La<sup>2+</sup> is also bound to some extent. Complexation with a cation is always in a 1:1 stoichiometry and pH-dependent. Essentially no binding of Ca<sup>2+</sup> occurs below pH 7.0 and maximum binding takes place at pH 9.5.<sup>2</sup>

lonomycin has significant advantages for use in studies of Ca<sup>2+</sup> transport across biological membranes. <sup>3</sup> It is also used to equilibrate intracellular and extracellular calcium ion levels for *in situ* calibrations of fluorescent indicators. <sup>4</sup> The resultant calcium flux leads to several downstream effects, such as up-regulation of CD7 in T cells (signal of activation)<sup>5</sup>, or the hydrolysis of phosphoinositides and activation of Protein Kinase C in T cells. <sup>6</sup> It was found to have antiproliferative effects on the human bladder cancer cell line both *in vitro* and in *vivo*. <sup>7</sup> Ionomycin can serve as an inducer of apoptosis <sup>8</sup> which was suggested to act by activation of a latent, calcium-responsive endonuclease <sup>9</sup>.

Purity: ≥98% (HPLC)

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices

#### **Preparation Instructions**

The product, a lyophilized semi-solid, is soluble in methanol and in DMSO. It is routinely tested in acetonitrile at 10 mg/ml, yielding a clear solution.

### Storage/Stability

Store desiccated and protected from light at -20 °C. Under these conditions the product is stable for 2 years. Stock solutions in DMSO, stored at -20 °C, and protected from light, are stable for several months.

## References

- 1. Liu, W-C. et al., Ionomycin, a new polyether antibiotic., *J. Antibiot.*, **31**, 815-19 (1978).
- Liu, C-M., and Hermann, Characterization of ionomycin as a calcium ionophore., T. E., *J. Biol. Chem.*, 253, 5892-94 (1978).
- Beeler, T.J.,. et al., The effect of ionomycin on calcium fluxes in sarcoplasmic reticulum vesicles and liposomes., *J. Biol. Chem.*, 254, 6229-31 (1979).
- 4. Kao, J.P., Practical aspects of measuring [Ca<sup>2+</sup>] with fluorescent indicators., *Methods Cell Biol.*, **40**, 155-81 (1994)
- Ware, R.E., et al., Induction of T cell CD7 gene transcription by nonmitogenic ionomycin-induced transmembrane calcium flux., *J. Immunol.*, 147, 2787-94 (1991)
- 6. Chatila, T., et al., Mechanisms of T cell activation by the calcium ionophore ionomycin, *.J. Immunol.*, **143**, 1283-89 (1989)

- 7. Miyake, H., et al., Calcium ionophore, ionomycin inhibits growth of human bladder cancer cells both *in vitro* and *in vivo* with alteration of Bcl-2 and Bax expression levels., *J. Urol.*, **162**, 916-21 (1999)
- 8. Miyake, H., et al., Stress protein GRP78 prevents apoptosis induced by calcium ionophore, ionomycin, but not by glycosylation inhibitor, tunicamycin, in human prostate cancer cells., *J. Cell Biochem.*, **77**, 396-408 (2000).
- 9. Aagaard-Tillery, K.M. and Jelinek, D.F., Differential activation of a calcium-dependent endonuclease in human B lymphocytes. Role in ionomycin-induced apoptosis., *J. Immunol.*, **155**, 3297-307 (1995).

NDH,PHC 03/06-1