


92544 Abberior® FLIP 565, maleimide

Description

This photoswitchable marker belongs to the photosensitive spiroamide compounds which, upon photoswitching with UV light at 360–375 nm, form transiently fluorescent species. In the dark, the colored and fluorescent form of Abberior FLIP 565 returns to the nonfluorescent state via thermal relaxation. The compound is performing very well in PALM, STORM and GSDIM microscopy. Below is shown an image taken with Abberior FLIP 565 with a Nikon N-STORM microscope at the Nikon Imaging Application Center in Hamburg. Due to the very short lifetime of the emitting form the compound is less suited for conventional fluorescence microscopy.



Key Features

- Ideal for PALM, STORM, GSDIM microscopy
- Compatible with 2-photon excitation (for optical sectioning)
- Multiple switching cycles

Applications

Abberior FLIP 565 can also undergo two-photon photoactivation with intense IR light (~ 760 nm). Consequently, this label is dedicated for the recording images with optical sectioning, i.e. activating a thin layer of ~ 500 nm.



Chemical Data : Abberior® FLIP 565

Structure:	on request
Formula:	on request
Molecular weight:	907.9 g/mol (NHS ester), 932.9 (maleimide)
Solubility:	water, acetonitrile, methanol, DMSO, DMF
Polarity:	moderately hydrophilic
Charge:	–1 (when caged or conjugated)
Purity:	> 90 %

Photophysical Data: Abberior® FLIP 565

Absorption Maximum (off-state), λ_{\max} :	314 nm (PBS, pH 7.4),
Extinction Coefficient, $\epsilon(\lambda_{\max})$:	47,000 M ⁻¹ cm ⁻¹ (PBS, pH 7.4)
Correction Factor, $CF_{260} = \epsilon_{260}/\epsilon_{\max}$:	tbd
Correction Factor, $CF_{280} = \epsilon_{280}/\epsilon_{\max}$:	tbd
Fluorescence Maximum, λ_{fl} :	580 nm (PBS, pH 7.4)
Photoactivation Wavelength:	310-380 (one-photon activation) 650-800 (two-photon activation)
Fluorescence Quantum Yield, η :	0.38 (PBS, pH 7.4),
Fluorescence Lifetime, τ :	tbd

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.



Literature

V. N. Belov et.al. "Rhodamines NN: A Novel Class of Caged Fluorescent Dyes", *Angew. Chem. Int. Ed.* 49, 3520–3523 (2010)

G. Y. Mitronova et.al. "New Fluorinated Rhodamines for Optical Microscopy and Nanoscopy", *Chem. Eur. J.* **16**, 4477–4488 (2010)

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

