

Quick Start Guide

GenElute™ HP Plasmid Miniprep Kit

NA0150, NA0160

Reagents to Prepare

- Spin RNase Solution tube briefly and add 13 μ L (10 preps), 78 μ L (70 preps) or 500 μ L (350 preps) of RNase A Solution to the Resuspension Solution. Store at 2-8°C.
- Dilute the Wash Solution 2 with 10 mL (10 prep), 48 mL (70 prep) or 300 mL (350 prep) of 95-100% ethanol. After each use, tightly cap the diluted Wash Solution to prevent ethanol evaporation.

Protocol

All spins at $\geq 12,000 \times g$.

Harvest & Lyse Bacteria

1. Pellet cells from 1-5 mL overnight culture by spinning for 1 min. Discard supernatant.
2. Resuspend cells in 200 μ L Resuspension Solution. Pipet or vortex to mix.
3. Add 200 μ L Lysis Solution. Gently invert 6-8 times to mix. Allow to clear for ≤ 5 min.

Prepare Cleared Lysate

4. Add 350 μ L Neutralization Solution. Gently invert 4-6 times.
5. Pellet debris 10 min.

Bind DNA to Column

6. Transfer cleared lysate into column in a microcentrifuge tube.
7. Spin 1 min. Discard flowthrough.

Wash to Remove Contaminants

8. Add 500 μ L Wash Solution 1 to column. Spin 1 min. Discard flowthrough.
9. Add 750 μ L Wash Solution 2 to column. Spin 1 min. Discard flowthrough.
10. Spin 1 min to dry column.

Elute Purified DNA

11. Transfer column to new collection tube.
12. Add 100 μ L Elution Solution to column. Spin 1 min.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/TechService.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/Terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/Offices.