

Product Information


MJ 33

Product Number **M 3315**

Storage Temperature 2-8 °C

CAS #: 199106-13-3

Synonyms: 1-Hexadecyl-3-(trifluoroethyl)-*sn*-glycero-2-phosphomethanol, lithium



Product Description

Molecular Formula: C₂₂H₄₃F₃O₆PLi

Molecular Weight: 498.5

Appearance: white solid

Purity: 90% by NMR

MJ 33 is an active site-directed specific inhibitor of Type I (pancreatic) and bee venom phospholipase A₂ (PLA₂). PLA₂ designates a class of enzymes with an approximate molecular weight of 14 kDa that hydrolyzes the *sn*-2 ester of glycerophospholipids to produce a fatty acid and a lysophospholipid. Based on amino acid sequences, the known PLA₂s can be divided into ten distinct groups with specific functions and tissue distribution. Since naturally-occurring phospholipids are membrane constituents that are essentially insoluble in water, in order to access the substrate PLA₂s must bind to the lipid-aqueous interface as functionally active monomers. Subsequently PLA₂ binds to the substrate molecule to form a complex that then becomes activated on an ionic interface.¹

Many compounds have been identified that are nonselective inhibitors of PLA₂ isozymes and, thus, cannot be used to study the role of specific PLA₂ isozymes in complex cellular processes. Determination of crystal structures of the PLA₂ enzyme made it possible to identify anionic binding sites on PLA₂. Over 100 amphiphilic phosphoesters were tested for their ability to inhibit the active site of phospholipase A₂.² MJ 33 is a competitive, active site-directed inhibitor, that selectively complexes with Type 1B PLA_{2A} (bee venom, pancreatic, and acidic lung PLA₂). In crystallographic studies, MJ 33 and five sulfate or phosphate anions are bound between the two subunits

of the PLA₂ dimer. The *sn*-2-phosphate of MJ 33 binds to the active site of the A subunit of PLA_{2A} and the alkyl chain extends into the active site slot of the B subunit across the subunit interface.¹

MJ 33 at 3 mol% produced 95% inhibition of pancreatic PLA₂. Mol % is defined as ([MJ 33]/([MJ 33] + [PG])) × 100, where PG is the concentration of phosphatidyl-glycerol. Bee venom, pancreatic, and acidic lung PLA₂s share common characteristics. All are soluble enzymes that are active at acidic pH (4.5–5.0) and are calcium-independent. MJ 33 has poor affinity for the Type II human synovial PLA₂ and has only moderate affinity toward lysosomal PLA₂ isolated from macrophages (IC₅₀ = 15 mol%).^{3,4,5} Thus, MJ 33 is a valuable tool for determining the role of Type 1B PLA₂ in cellular processes.

Preparation Instructions

MJ 33 is soluble in water at 10 mg/ml with warming to 60 °C.

Storage/Stability

Store at 2 to 8 °C under nitrogen, in a desiccator and protected from light.

References

1. Pan, Y. H., et al., Five coplanar anion-binding sites on one face of phospholipase A₂: relationship to interface binding. *Biochemistry*, **40**, 609-617 (2001).
2. Jain, M. K., et al., Active-site-directed specific competitive inhibitors of phospholipase A₂: novel transition-state analogues. *Biochemistry*, **30**, 10256–10268 (1991).
3. Fisher, A. B., Dodia, C., Lysosomal-type PLA₂ and turnover of alveolar DPPC. *Am. J. Physiol. Lung Cell Mol. Physiol.*, **280**, L748-L754 (2001).

4. Akiba, S., et al., Characterization of acidic Ca^{2+} -independent phospholipase A_2 of bovine lung. .Comp. Biochem. Physiol. B. Biochem. Mol. Biol., **120**, 393-404 (1998).
5. Shinozaki, K., Waite, M., A novel phosphatidyl-glycerol-selective phospholipase A_2 from macrophages. Biochemistry., **38**, 1669-1675 (1999)

AH 7/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.