

Product Information

MSMLS™ Plate 6 (Lipophilic)

Supplied by IROA Technologies, LLC.

Catalog Number **MSMLS06**

Storage Temperature $-20\text{ }^{\circ}\text{C}$

Product Description

MSMLS06 (Mass Spectrometry Metabolite Library Plate 6) is a collection of high quality small biochemical molecules that span a broad range of primary metabolism. These are high purity (>95%) compounds supplied in an economical, ready-to-use format.

The library is most commonly used to provide retention times and spectra for key metabolic compounds, help optimize analytical mass spectrometry protocols, and qualify and quantify mass spectrometry sensitivity and limit of detection.

MSMLS™ Plate 6 comes with MLSDiscovery™, a software tool to support the extraction, manipulation, and storage of the data generated when using the MSMLS06. For further information on the software, to download, and for manual and video links please visit: www.sigmaaldrich.com/catalog/product/sigma/msmls06

Components

MSMLS Plate 6 compounds are conveniently provided at 5 µg per well, enough for multiple injections, suitable for manual and automated workflow.

The library is intended to be used for mass spectrometry metabolomics applications and provides a broad representation of primary metabolism.

Occasionally the plate map will change due to the availability of compounds. The compounds of each row have distinct molecular masses and can be multiplexed, but users should refer to the plate map for their specific lot before proceeding.

The plate map contains descriptors and represents information gathered from multiple databases. Please note this data may contain some errors. It is recommended for users to carefully review plate composition and description information prior to use. To help build a better database, please report any observed discrepancies.

MSMLS06 includes:

- 1 polypropylene plate in 96 well format
Polypropylene deepwell (1.2 mL, total volume per well) plate (MasterBlock®, Greiner Number 780215) in combination with seals (VIEWseal™, Greiner Number 676070)
- 5 µg (dried weight) of each metabolite

MLSDiscovery software includes:

- Plate map
- Alpha-numeric assigned position
- Descriptors:
 - Name
 - Parent CID
 - KEGG ID where available or ChemSpider ID
 - Molecular formula
 - Molecular weight
 - CAS Registry number
 - ChEBI
 - HMDB ID/YMDB ID
 - PubChem Compound and Substance ID
 - Metlin ID

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

The following are suggestions and dependent on user chromatography and instrumentation.

This plate contain primarily lipid-like compounds with the exception of the water soluble sugar compounds. It is recommended to solubilize the lipid-like compounds using a chloroform:methanol (1:1) mixture. For sugar compounds it is recommended to solubilize with ultrapure water.

Pool compounds for multiplexing. Again, be sure to check the plate map to ensure one can adequately separate the compounds using the chromatographic system prior to pooling.

Storage/Stability

Store the plate at $-20\text{ }^{\circ}\text{C}$.

Once the metabolites are dissolved, the plate should be resealed and kept at $-20\text{ }^{\circ}\text{C}$ or $-80\text{ }^{\circ}\text{C}$ for long-term storage and protected from light. Avoid repeated freeze/thaw cycles.

Procedure

The compounds of MSMLS Plate 6 can either be used as standards and injected individually or mixed in such a way that the entire library may be examined with reasonable efficiency. Mixing compounds by row mixtures allows multiple compounds to be analyzed per injection. Be sure to check the plate map to ensure one can adequately separate the compounds using the chromatographic system prior to pooling.

The following are only suggestions and depend on user chromatography and instrumentation.

1. Individual Injections

As standards, each well represents a single compound. The entire library may be examined in great detail with several injections, one for each of the unique metabolites (Total volumes for each well of $250\text{ }\mu\text{L}$ – 1 mL may be considered).

2. Simple multiplex injections

If the rows of the plate are pooled, then the entire collection may be analyzed with 8 injections of simple mixtures. Keep the total volume for each well to $150\text{ }\mu\text{L}$ or less to prevent loss due to dilution and use 5 – $10\text{ }\mu\text{L}$ of each well for the pooled sample. Then, inject 2 , 4 , or $6\text{ }\mu\text{L}$ of the pooled material as needed.

References

1. Wishart, D.S. et al., HMDB: the Human Metabolome Database. *Nucleic Acids Res.*, 2007 Jan; 35 (Database issue):D521-6. 17202168.
2. Wishart, D.S. et al., HMDB: a knowledge base for the human metabolome. *Nucleic Acids Res.*, 2009 37 (Database issue):D603-610. 18953024.
3. Wishart, D.S. et al., HMDB 3.0 — The Human Metabolome Database in 2013. *Nucleic Acids Res.*, 2013. Jan 1;41(D1):D801-7. 23161693
4. Hastings, J. et al., The ChEBI reference database and ontology for biologically relevant chemistry: enhancements for 2013. *Nucleic Acids Res.*, 2013.
5. Kanehisa, M., and Goto, S., "KEGG: Kyoto Encyclopedia of Genes and Genomes". *Nucleic Acids Res.*, 2000 28 (1): 27– 30. doi:10.1093/nar/28.1.27. PMC 102409.PMID 10592173.
6. Tautenhahn, R. et al., An accelerated workflow for untargeted metabolomics using the METLIN database. *Nature Biotechnology*, 2012 30: 826–828. doi:10.1038/nbt.2348.
7. Smith, C.A. et al., METLIN: a metabolite mass spectral database. *The Drug Monit.*, 2005 27 (6): 747–51. doi:10.1097/01.ftd.0000179845.53213.39. PMID 16404815.
8. Bolton, E. et al., PubChem: Integrated Platform of Small Molecules and Biological Activities. Chapter 12 IN Annual Reports in Computational Chemistry, Volume4, American Chemical Society, Washington, DC. 2008 Apr [free author manuscript]

MSMLS and MLSDiscovery are trademarks of IROA Technologies, LLC.

MasterBlock is a registered trademark and VIEWseal is a trademark of Greiner Bio-One GmbH.

SS,AA,MAM 04/19-1