

ProductInformation

Bisindolylmaleimide X Hydrochloride

Product Number **B 3931**


Storage Temperature -20°C

Synonym: 2-(8-Aminomethyl-6,7,8,9-tetrahydropyrido[1,2-a]indol-3-yl)-3-(1-methylindol-3-yl)maleimide, HCl Salt;
3-(8-Aminomethyl-6,7,8,9-tetrahydropyrido[1,2-a]indol-10-yl)-4-(1-methylindol-3-yl)-1H-pyrrole-2,5-dione, HCl Salt; Ro 31-8425

Product Description

Molecular formula: $\text{C}_{26}\text{H}_{24}\text{N}_4\text{O}_2\cdot\text{HCl}$

Mol. wt.: 461.0



Bisindolylmaleimides are potent, selective inhibitors of protein kinase C (PKC). They are structurally similar to the naturally occurring molecule, staurosporine, but they are more selective for PKC over other protein kinases. Bisindolylmaleimides are used to selectively probe for PKC-mediated pathways for transduction of hormone, cytokine, and growth factor signals.

Bisindolylmaleimides inhibit PKC by interacting with the catalytic subunit. Inhibition is competitive with ATP. Studies of structure-activity relationships of analogs indicate that cationic substituents at the indole nitrogen increase the potency as an inhibitor of PKC.

Bisindolylmaleimide X inhibits PKC within intact platelets and T cells. It also inhibits Fas-mediated apoptosis and T cell-mediated autoimmune diseases.

The selectivity of bisindolylmaleimide X for rat brain PKC over two other protein kinases is shown in the table below.

Enzyme	IC ₅₀
Protein Kinase C	15 nM
cAMP-Dependent Protein Kinase	2800 nM
Phosphorylase Kinase	1300 nM

Preparation Instructions

Prepare stock solutions in DMSO.

Storage/Stability

Store product at -20°C . Protect from light.

References

- Davis, P.D., et al., Inhibitors of protein kinase C. 2. Substituted bisindolylmaleimides with improved potency and selectivity. *J. Med. Chem.*, **35**, 994-1001 (1992).
- Wilkinson, S.E. et al., Isozyme specificity of bisindolylmaleimides, selective inhibitors of protein kinase C. *Biochem. J.*, **294**, 335-337 (1993).
- Zhou, T., et al., Bisindolylmaleimide VIII facilitates Fas-mediated apoptosis and inhibits T cell-mediated autoimmune diseases. *Nat. Med.*, **5**, 42-48 (1999).
- Bit, R. A., et al., Inhibitors of protein kinase C. 3. Potent and highly selective bisindolylmaleimides by conformational restriction. *J. Med. Chem.*, **36**, 21-29 (1993).

JWM 12/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.