

ECMatrix™-411 E8 Laminin Substrate

Stem Cell Reagent

Cat. # CC162-1050UG

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.
NOT FOR HUMAN OR ANIMAL CONSUMPTION.

Pack size: 1050 µg

Store at 2-8°C



Data Sheet

page 1 of 2

Description

Endothelial cells line the interior surface of blood and lymphatic vessels, forming an interface between circulating blood or lymph in the lumen and the rest of tissue or organ. Laminin-411 is known to bind to the integrin $\alpha 6 \beta 1$ which is located on the cell surface of pluripotent stem cells and is known to induce endothelial cell differentiation. Professor Kiyotoshi Sekiguchi's group (Matrixome, Inc.) have produced a recombinant E8 fragment of laminin-411 at large-scale while retaining the full integrin binding activity. ECMatrix™-411 E8 Laminin Substrates are recombinant Laminin-411-E8 fragments which bind to integrin $\alpha 6 \beta 1$ and increases the differentiation of pluripotent stem cells (ES/iPSCs) into acetyl-LDL⁺/CD31⁺ endothelial cells.

Storage and Handling

ECMatrix™-411 E8 Laminin Substrates should be stored at 2-8°C. Avoid multiple freeze-thaw cycles and protect from light.

Presentation

1) 6 X 175 µg ECMatrix™-411 E8 Laminin Substrate (0.5 mg/mL in PBS). Expressed in CHO-S cells.

Quality Control Testing

- Purity (SDS-Page): > 95%
- Endotoxin Test: ≤ 750 EU/mg
- Mycoplasma Test: Negative
- Sterility Test: Negative
- Integrin Binding Assay (kDa): ≤ 10 nM

Protocol

Precoating Method

1. Dilute the 0.5 mg/mL stock solution with sterile PBS to achieve a 2.5 µg/mL working solution.
2. Coat dishes with ECMatrix™-411 at 0.25 µg/cm² (for example, for one well of a 6-well plate add 1 mL of the 2.5 µg/mL working solution).
3. Incubate for 1 hour at 37°C, 3 hours at room temperature or overnight at 4°C.
4. Before use, remove remaining fluid from the coated surface (do not rinse).
5. Detach cells into small clumps using Accutase.
6. Plate the cells at desired density.

Note: Do not allow the plates to dry, briefly spin down all liquids in the tube before use, avoid repeated freeze-thaw cycles.

References


1. Takayama K, et al. Laminin 411 and 511 promote the cholangiocyte differentiation of human induced pluripotent stem cells. *Biochem Biophys Res Commun.* 2016 May 20;474(1):91-96.
2. Ohta R, et al. Laminin-guided highly efficient endothelial commitment from human pluripotent stem cells. *Sci Rep.* 2016 Nov 2;6:35680.

SPECIES LEGEND: H Human Ca Canine M Mouse R Rat Rb Rabbit B Bovine P Porcine WR Most Common Vertebrates

Please visit www.millipore.com for additional product information and references.

Submit your published journal article, and earn credit toward future purchases. Visit www.millipore.com/publicationrewards to learn more!

**KhES-1 Human ES Cell
Derived Endothelial Cells**



**253G4 Human ES Cell
Derived Endothelial Cells**




Figure 1. Differentiation of human embryonic stem cells into endothelial cells. Human ES cell lines (KhES-1, 253G4) differentiated into endothelial cells grown on ECMatrix™-411 E8 Laminin Substrates express acetyl-LDL (red) and CD31 (green). Cells counterstained with Dapi (blue).

ECMatrix™ is a trademark of Merck KGaA
Matrigel is a registered trademark of Corning Incorporated

■ antibodies ■ Multiplex products ■ biotools ■ cell culture ■ enzymes ■ kits ■ proteins/peptides ■ siRNA/cDNA products

Please visit www.millipore.com for additional product information, test data and references

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EMD Millipore®, the M mark, Upstate®, Chemicon®, Linco® and all other registered trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2019 Merck KGaA, Darmstadt, Germany. All rights reserved.



We Buy 100% Certified
Renewable Energy