NovaSeal™ Battery Powered Crimping Tool

User and Maintenance Guide

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Millipore Corporation. Millipore Corporation assumes no responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall Millipore Corporation be liable for incidental or consequential damages in connection with or arising from the use of this manual.

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Symbols

The following symbols are used throughout this guide:

Safety Warnings

Pinching may occur.

Read this manual

Introduction

The NovaSeal Battery Powered Crimping Tool is a hand-held electrohydraulic crimping tool which is a component of the NovaSeal Crimping Solution. The Crimping Tool allows a Mobius® disposable assembly to be cut and crimped at a specific location while maintaining sterility and integrity of the fluid paths.


Figure 1: NovaSeal Battery Powered Crimping Tool

Catalogue Numbers

Catalogue Number	Description
NM0029	NovaSeal Battery Powered Crimping Tool 220-240 V European plug
NM0030	NovaSeal Battery Powered Crimping Tool 100-120 V US and Japanese plug
NM0038	NovaSeal Battery Powered Crimping Tool 220-240 V UK plug
NM0021	Die kit $\frac{1}{4} \times \frac{1}{2}$ in. (ID x OD) tubing
NM0022	Die kit $\frac{3}{8} \times \frac{5}{8}$ in. (ID x OD) tubing
NM0023	Die kit $\frac{1}{2} \times \frac{3}{4}$ in. (ID x OD) tubing
NM0024	Battery, 9.6 V Ni-Cd
NM0039	$\frac{1}{4} \times \frac{1}{2}$ in. (ID x OD) Crimped pinch pipe cover
NM0041	$\frac{3}{8} \times \frac{5}{8}$ in. (ID x OD) Crimped pinch pipe cover
NM0040	$\frac{1}{2} \times \frac{3}{4}$ in. (ID x OD) Crimped pinch pipe cover

Operator and Equipment Safety

Any attempt to use the tool in a manner not specified by Millipore Corporation may result in damage to the tool, voiding of product warranty, and possible operator injury.

All operators of the tool must read and understand this User and Maintenance Guide before using this tool. Failure to follow operating instructions could result in operator injury, improper crimp or damage to the tool.

Prior to operation, the tool must be fully assembled as shown in this guide.

Use appropriate personal protective equipment and eye protection when operating the tool.

DO NOT expose the tool to liquid spray or immerse the tool in any liquid.

DO NOT autoclave or steam the tool.

DO NOT use without a properly installed die.

DO NOT place the battery where it may come in contact with conductive materials.

Electric engine may cause sparks. Do not use near flammable liquids or gasses.

Do not insert fingers into tool, pinching may occur.

Unpacking

The NovaSeal Battery Powered Crimping Tool is shipped in two containers, the Crimping Tool Box and the Die Box. If any components are missing or damaged, contact your local Millipore representative.


Figure 2: Crimping Tool Box contents (User and Maintenance Guide and Certificate of Quality not shown)




Figure 3: Die Box contents (User and Maintenance Guide not shown).

Charging the Battery

1. Inspect the battery charger before using to insure there are no foreign objects in it. Plug the battery charger into an appropriate power receptacle.
2. Remove the battery from the tool.

3. Place the battery into the battery charger. Allow the battery to charge for approximately 60 minutes. The green LED will illuminate when the battery is fully charged. Two or three charging cycles may be required before a new battery is fully charged.

4. The battery should be recharged after 50 crimping cycles.

Battery charger LED indicators

Indicator	Explanation
 Flashing green	Ready for battery.
 Steady red	Charging battery.
 Steady green	Battery fully charged.
 Flashing red	Battery too hot or too cold.
 Flashing green and red	Defective battery.
 No light	No power to battery charger or battery is polarized incorrectly.

Assembling the Crimping Tool



DO NOT install the battery on the Crimping Tool until the Crimping Tool is fully assembled.


Pinching may occur.

1. Push the locking pin toward the tool (A), turn it counter-clockwise (B), then pull the locking pin out and hold (C).
3. Push the locking pin in until it locks the jaws. A “click” will be heard. Once the jaws are installed they cannot be removed.

2. Insert the jaws into the Crimping Tool.


Installing and Removing the Die


 **DO NOT** install the battery on the Crimping Tool until the Crimping Tool is fully assembled and the die is installed.

The die should be changed after approximately 1000 crimps or when worn, whichever occurs first. Millipore recommends performing a test crimp each time a die is installed to ensure proper crimping.


1. Inspect the cutting edge on the die and ensure that it is not worn and that it is free of foreign particles.

2. Verify that the die and fitting sizes match. No more than 1 mm on each side of the fitting should extend beyond the die.

3. Place the die on a stable, flat surface. The lowest edge of the die should face away from the tool. Gently squeeze the jaws to open them just enough to accommodate the die. Slide the jaws onto the die. Pull the die locks outward so that the jaws slide down over the die. Release the locks once the die is in place to lock the die in position.


4. To remove the die, place the die on a stable, flat surface. Pull the die locks outward while raising the tool so that the jaws slide up over the die.

Installing and Removing the Battery


DO NOT install the battery on the Crimping Tool until the Crimping Tool is fully assembled and the die is installed.

1. Slide the battery into the bottom of the tool until it locks into position. The LED will flash two times to indicate the battery is installed properly.

2. The Crimping Tool is now ready to use.

To remove the battery, press the battery lock on the Crimping Tool and slide the battery out of the tool.

Operating the Crimping Tool

It is recommended that a crimping test be performed using these test pieces after the tool is fully assembled or after the die is changed.

Once properly assembled, the crimping head may be rotated 360°.

Use appropriate personal protective equipment and eye protection when operating the tool.

DO NOT use in temperatures exceeding 50 °C (122 °F).

Pinch pipes are plated with nickel (Ni).

DO NOT use on systems where tubing pressure will exceed 0.5 bar (7.25 psi).

DO NOT use lubricants or liquid between the pinch pipe and tubing.

After 50 crimping cycles, allow the tool to cool down for 15 minutes before resuming use.


Do not insert fingers into tool, pinching may occur.


Notes:

- A complete crimping cycle is required for a proper crimp.
- The crimping cycle may be stopped at any time by releasing the trigger.
- Inspect the cutting edge on the die and ensure that it is not worn and that it is free of foreign particles.
- Inspect pinch pipes for damage before each crimp. Do not use deformed pinch pipes.


1. Open the crimping jaws by squeezing the jaws closest to the handle.


2. The die and fitting sizes must match. No more than 1 mm on each side of the fitting should extend beyond the die. Inspect the die to be sure it is free of foreign objects. Place the tubing fitting into the die and ensure that it is properly seated in the die.

3. Close the jaws by slowly releasing them. Holding the crimper with both hands, press and hold the trigger for ten seconds. If the LED illuminates, refer to the Troubleshooting section in this guide.


4. Release the trigger.
5. Open the jaws by squeezing the jaws closest to the handle. Verify that the tubing pinch pipe has been cut completely. Remove the tubing from the tool by pulling the tubing straight out.

The crimped edge may be sharp.

6. Inspect the crimped end of each pinch pipe. The edges should be sharp, the crimps should be dimpled and the crimps on each tubing piece should be of equal length.
7. Install optional end covers if desired.

Do not attempt to remove the pinch pipe from the tubing once the crimping cycle is complete as the integrity of the crimp will be compromised.


Maintenance

The die should be cleaned with a soft bristle brush to remove debris.

There are no user-serviceable parts inside the Crimping Tool. Service should be performed by Millipore Service Technicians only.

The Crimping Tool should be serviced after 10,000 crimping cycles or one year, whichever comes first.

The crimping Tool may be wiped down with a soft cloth dampened with Ethanol.

The user must dispose of the tool and battery in accordance with local regulations. In the European Union (EU), this unit is subject to the scope of the WEEE (2002/96/EG) directive which regulates the collection and recycling of electric and electronic units and the RoHS (2002/95/EEC) directive which bans new electrical and electronic equipment containing more than 0.1 weight percentage lead, mercury, hexavalent chromium, polybrominated diphenyl ethers (PBDE) and 0.01 weight percentage cadmium per homogenous material. The battery must be disposed of according to the EEC battery guideline or local regulations.

Troubleshooting

In most instances, the suggested actions/resolutions should resolve the problem. If the problem persists or is not included in this list, contact Millipore Technical Service.

Problem	Possible Causes	Remedy
Crimp pressure is low	Service required	Return to Millipore
Tool leaks oil	Service required	Return to Millipore
Rough edge on cut	Worn die	Replace die
Tool LED remains lit for 20 seconds after crimp	Battery discharged	Recharge battery
Tool LED flashes for 20 seconds after crimp	Service required	Return to Millipore
Tool remains lit for 20 seconds then flashes intermittently	Battery discharged, service required	Return to Millipore
Battery charger flashes green and red	Defective battery.	Replace battery.

Specifications

Parameter	Value
Weight	2.7 kg (6 lbs)
Compression Time	5 to 10 seconds
Thrust Force	15 kN
Battery Voltage	9.6 V
Battery Capacity	2.0 Ah
Battery Charging Time	approximately 60 minutes
Cycles per charge	50
Operating Temperature Range	5 to 50 °C
Sound Level	70.6 dB (A) at a distance of 1 m
Vibrations	less than 2.5 m/s ²

Standard Warranty

Millipore Corporation (“Millipore”) warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore’s published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon. In the event of a breach of the foregoing warranty, Millipore’s sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER INDIRECT DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.

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Millipore は Millipore Corporation の登録商標です。NovaSeal および Mobius は Millipore Corporation の商標です。© 2007, 2009 Millipore Corporation. All rights reserved.
Millipore and Mobius are registered trademarks of Millipore Corporation. NovaSeal and the M mark are trademarks of Millipore Corporation. © 2007, 2009 Millipore Corporation. All rights reserved.