



## BAIRD PARKER SULPHA SUPPLEMENT

Product Number **B2052**  
Storage Temperature 2-8°C

### Product Description

Baird Parker Sulpha Supplement is an antibacterial supplement that can be used for the suppression of the *Proteus* species growing on Baird Parker Agar Base.

### Components

| <u>Item</u>     |       |
|-----------------|-------|
| Sulphamethazine | 50 mg |

### Precautions and Disclaimer

For laboratory use only. Not for drug, household or other uses.

## Product Information

### Preparation Instructions

Rehydrate the contents of 1 vial with 5 mls sterile filtered 0.01 N sodium hydroxide. Mix well and aseptically add to 1000 mls of sterile, molten Baird Parker Agar Base B1051 containing 50 mls of Egg Yolk Tellurite Emulsion E7774. Mix well and pour into sterile petri plates. plates. Allow the tubed medium to solidify in a slanted position.

### Storage/Stability

Store the supplement at 2-8°C.

### References

1. Baird-Parker, A.C., (1962). J. Appl. Bact. 25,12.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.