Product Information

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Glycerol 3-phosphate Oxidase from *Aerococcus viridans*

Catalog Number **G9888** Storage Temperature –20 °C

CAS RN 9046-28-0 EC 1.1.3.21 Synonyms: GPO, α-Glycerophosphate oxidase, sn-glycerol-3-phosphate:oxygen 2-oxidoreductase

Product Description

Many bacteria and yeast can utilize glycerol as a carbon source. After uptake by the cell, glycerol is phosphorylated to α -glycerol-3-phosphate, which in turn is oxidized to enter the glycolytic pathway. α -Glycero-phosphate oxidase (GPO) catalyzes the oxidation of α -glycerol-3-phosphate to dihydroxyacetone phosphate by the following reaction:

GPO

Glycerol-3-PO₄ + O₂ \rightarrow dihydroxyacetone-PO₄ + H₂O₂

GPO has been used for sensitive metabolite assays of starch and lipid synthesis, pyrophosphate, ATP, ADP, and most glycolytic intermediates in *Arabidopsis* seeds.¹ GPO is part of the dihydroxyacetone phosphate:glycerol-3-phosphate cycle in the bloodstream form of *Trypanosoma brucei*.²

```
Molecular mass:<sup>3</sup> 63 kDa (SDS-PAGE)
Cofactor:<sup>4</sup> FAD
Optimal pH:<sup>3</sup> 7.5–8.0
Optimal temperature:<sup>2</sup> 37 °C
K_{M}:<sup>4</sup> 2.3 mM
```

Activator:² Glycerol Inhibitors:⁴ benzylformic acid, glyoxylic acid, methylglyoxal

This product is purified from *Aerococcus viridans*. It is supplied as a lyophilized powder.

Protein: ≥60% (Lowry), balance primarily sucrose

Specific activity: ≥70 units/mg solid

Unit definition: One unit will oxidize 1.0 μ mole of L-glycerol-3-phosphate to dihydroxyacetone phosphate with the formation of hydrogen peroxide per minute at pH 8.1 at 37 °C.

GPO is assayed spectrophotometrically in a 1.05 ml reaction mixture containing 60 mM Tris HCl, 5 units peroxidase, 0.06% (v/v) Triton[™] X-100, 0.01% (w/v) 4-aminoantipyrine, 0.02% (w/v) phenol, 95 mM DL-α-glycerophosphate, 0.01% bovine serum albumin, and 0.004–0.009 units GPO, at pH 8.1 at 37 °C.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

GPO is soluble (0.45 unit/ml) in cold 20 mM Tris HCl, pH 7.5 at 37 °C, containing 0.2% (w/v) bovine serum albumin. Dissolve immediately before use. A study of solution storage of GPO has been reported.⁴

Storage/Stability

Store the product at -20 °C. When stored at -20 °C, GPO should retain activity for two years.

References

- 1. Gibon, Y. et al., Plant J., 30(2), 221-235 (2002).
- 2. Opperdoes, P. et al., Eur. J. Biochem., **76(1)**, 29-39 (1977).
- Streitenberger, S.A. et. al., Appl. Microbiol. Biotechnol., 57(3), 329-333 (2001).
- 4. Macková, M. et al., Lett. Appl. Microbiol., **30(3)**, 188-191 (2000).

Triton is a trademark of The Dow Chemical Company or an affiliated company of Dow.

GRO, JWM, GCY, MAM 09/16-1

©2016 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.