


# LC-MS Contaminants

Avoid, identify, minimize.



## Table of Contents

Introduction: the consequences of contamination	4
Sample components that can interfere with LC-MS results	5
General system care, maintenance and laboratory practice	6
Sample preparation is crucial for minimizing contamination	7
Proper mobile phase preparation to minimize contaminants	17
Effects of column choice on LC-MS performance	24
Useful LC-MS resources	28
Appendix I. Common mass spectrometry contaminants and their sources	29
Appendix II. Monoisotopic ion masses of commonly observed repeating units in LC-MS	35

# How to identify and avoid contaminants in LC-MS

(Liquid Chromatography-Mass Spectrometry)

You thought your LC-MS analysis would be straightforward, but there are peaks you didn't expect.

In this technical bulletin, you will learn some tips on identifying LC-MS contaminants and avoiding contamination. Our 110 years of separation expertise, combined with our precision-manufactured products, give you the greatest chance of obtaining reproducible, clean data.

# Introduction: the consequences of contamination

Of the key parameters determining LC-MS success (sensitivity, resolution, reproducibility, efficiency, selectivity, speed, column lifetime), sensitivity, reproducibility and column lifetime are negatively affected by the presence of contaminants.

Parameters negatively affected by contaminants	What it is	How to minimize contaminants
Sensitivity	Lowest level of analyte detectable above background; sensitivity is reduced by ion suppression	Wash columns to mitigate column bleeding
		Effective sample preparation
		Use certified mobile phase
Reproducibility	Can mean column-to-column reproducibility or run-to-run reproducibility	Prepare sample with devices that do not introduce extractable impurities
		Run sufficient controls to verify run-to-run stability within a batch
		Remove sample components that interfere with separation, ionization and fragmentation
		Use high-quality HPLC columns
Column Lifetime	Number of injections on a column without change in selectivity and efficiency	Use a robust column with high matrix tolerance: e.g., monolithic columns have high tolerance
		Eliminate sample contaminants that adsorb strongly or ionize easily. Avoid polymers that, when fragmented, result in multiple peaks of varying m/z. Use guard column and proper sample preparation (such as centrifugation, filtration, and extraction) to remove particles and extend column lifetime.
		Fully elute/clean column after each sample

**Table 1.**  
LC-MS performance parameters.

# Sample components that can interfere with LC-MS results

Some biological matrices, such as plasma, contain high amounts of phospholipids. If not removed prior to chromatography, separating phospholipids from analytes of interest can require long chromatography run times and high concentrations of organic solvents. Furthermore, phospholipids can build up on analytical column, and unexpectedly elute in future runs. Drug formulation agents, such as polysorbitans and polyethylene glycol, can also interfere and cause ionization suppression.

Besides sample-derived contaminants, additional sources of contamination are sampling devices, solvent impurities, containers, sample preparation devices, volatile organics introduced as a result of handling personal care products, and even columns themselves.

Plasticizers from labware can interfere with LC-MS, resulting in the need to lengthen the chromatography run in order to resolve these peaks from analyte peaks.

A list of common contaminants, their molecular weight, and possible sources can be found in Table 2.

## Sample components that can interfere with LC-MS results include:

- Metabolites
- Detergents
- Salts/Buffer components
- Degradation products
- Counterions
- Matrix

Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
74.06059	[M+H] <sup>+</sup>	C <sub>3</sub> H <sub>7</sub> NO	Dimethyl formamide	Solvent
102.12827	[M+H] <sup>+</sup>	C <sub>6</sub> H <sub>15</sub> N	TEA	Triethylamine, buffer
107.0782	[A <sub>2</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
123.09222	[M+H] <sup>+</sup>	C <sub>7</sub> H <sub>10</sub> N <sub>2</sub>	DMAP	Dimethylaminopyridine, solvent
153.13917	[M+H] <sup>+</sup>	C <sub>9</sub> H <sub>16</sub> N <sub>2</sub>	DBU	1,8-Diazabicyclo[5.4.0]undec-7-ene
214.09018	[M+H] <sup>+</sup>	C <sub>10</sub> H <sub>15</sub> NO <sub>2</sub> S	n-BBS	n-butyl benzenesulfonamide, plasticizer
242.28477	M <sup>+</sup>	C <sub>16</sub> H <sub>36</sub> N	TBA	Tetrabutylammonium, buffer
279.15964	[M+H] <sup>+</sup>	C <sub>16</sub> H <sub>22</sub> O <sub>4</sub>	Dibutylphthalate	Plasticizer, phthalate ester
371.1018	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>5</sub>	Polysiloxane	Polysiloxane, followed by m/z 388
371.31614	[M+H] <sup>+</sup>	C <sub>22</sub> H <sub>42</sub> O <sub>4</sub>	DEHA	Bis(2-ethylhexyl) adipate, plasticizer
391.28484	[M+H] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticizer
445.12060	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>6</sub>	Polysiloxane	Polysiloxane, followed by m/z 462
447.2934	[M+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
515.41341	[M+H] <sup>+</sup>	C <sub>30</sub> H <sub>58</sub> O <sub>4</sub> S	DDTDP	Didodecyl 3,3'-thiodipropionate, antioxidant
519.13940	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>7</sub>	Polysiloxane	Polysiloxane, followed by m/z 536
593.15820	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>8</sub>	Polysiloxane	Polysiloxane, followed by m/z 610

**Table 2.**  
List of selected contaminants observed in mass spectra (ESI, positive mode, ion mass ≤1000 Da). Refer to Appendix I for a more complete list, or one of the databases listed in the resources section on page 28.

# General system care, maintenance and laboratory practice


**In addition to good laboratory practices, such as wearing powder-free, nitrile gloves and monitoring laboratory air (which can contain siloxanes and phthalates), follow these tips for minimizing contamination in LC-MS.**

- Flush HPLC system with organic eluent (preferably isopropanol or methanol; acetonitrile [ACN] can polymerize and block valves if system is stopped for several weeks) regularly to prevent microbial contamination. The interval of flushing depends on the eluents and buffers used and should be between two and four weeks.
- Pump debris is collected in the pump outlet filter. Some of these components can leach and be detected by MS. Replace the filter every 1–2 months or after changing from ACN to methanol (or vice versa) for lower baseline noise and general system protection.
- Filter frits attached to the inlets of the mobile phase tubing to protect the LC system from particulate matter should be made out of stainless steel. Cleaning of glass frits is time-consuming (buffer residue is hard to remove); in addition, silica and alkali are dissolved from the glass filter and form adducts  $[M+X]^+$ .

# Sample preparation is crucial for minimizing contamination

**Without sample preparation, samples contain components that are incompatible with HPLC/UHPLC/MS analyses:**

- Undissolved particles/precipitates in a sample clog and reduce the life of the chromatography column.
- Sample matrices may contain many impurities, making chromatograms challenging to interpret; for example, sample matrix contains components that either elute at the same point in the LC-MS chromatogram as the analyte (potentially causing ionization suppression) or affect analyte signal intensity.
- Particles held up on the column can leach contaminants into the mobile phase (in the current sample and subsequent samples), thereby increasing background.



**Figure 1.** Without sample preparation, the presence of contaminants in the sample results in more peaks in the chromatogram (as represented in the lefthand schematic), making analysis challenging.

Select a sample preparation method that brings the sample into a solution that is free of particles. Additional points of consideration include concentrating the analyte and reducing sample complexity. For example, a plasma sample might benefit from solid phase extraction, which removes contaminants (proteins, lipids) and also concentrates the sample, whereas fruit and vegetable juices, with their high particle load, might benefit simply from dilution and filtration.

Depending on the method chosen, sample preparation may be used, for example, to selectively enrich analytes, increase analyte concentration, or remove impurities that cause ionization suppression.

## How can you tell if your LC-MS analysis is suffering from ionization suppression?

### Consider performing the following steps to test for ionization suppression:

1. First, assess the detector response to a calibration standard under conditions of zero ionization suppression.
2. Spike an identical concentration of this standard into prepared sample matrix. Assess the detector response again to determine the effect of ionization suppression.
3. Assess the detector response when the spiked sample prepared in step 2 is processed using the sample preparation method(s) being considered.
4. Finally, add additional calibration standard to determine if the expected increase in signal is observed.

### To mitigate the interfering effects of ionization suppression, consider performing these steps:

- Dilute sample or reduce volume injected.
- Reduce ESI flow rate to the nL/min range—this will generate smaller, highly charged droplets that can resist the effects of nonvolatile salts, in case those have not been removed from the sample. Note that you should never use nonvolatile salts in the mobile phase.
- Choose a sample preparation method that removes contaminants causing ionization suppression. Using solid phase extraction instead of protein precipitation, for example, can reduce ionization suppression by phospholipids. Phospholipid removal is discussed further below.
- Change the strength of the mobile phase or the slope of the gradient, so that the analytes of interest may elute further from the solvent front and from the end of the gradient. In these regions of the chromatogram, ionization suppression is most likely to occur.

## Types of sample preparation commonly used for LC-MS

### Ten of the most popular sample preparation procedures currently in use (as ranked by percentage of survey respondents who reported using each method):

- |                   |  |                                    |
|-------------------|--|------------------------------------|
| 1. Filtration     | 5. pH adjustment   | 8. Sonication                      |
| 2. Centrifugation | 6. Vortexing   | 9. Solid phase extraction (SPE)    |
| 3. Dilution       | 7. Concentration (e.g., by ultrafiltration, precipitation) | 10. Liquid-liquid extraction (LLE) |
| 4. Evaporation    |  |                                    |

According to a 2013 survey of users performed by LC/GC magazine, filtration was the most commonly used sample preparation method. For complex samples containing components that contribute to high background and/or interfere with analyte ionization and fragmentation, filtration alone cannot provide the sample necessary for analysis, but forms an integral part of an overall sample preparation strategy, which involves other sample preparation techniques, like extraction, centrifugation and depletion.

## Sample preparation tips

### Choose an appropriate membrane filter to remove particles from your sample.

The presence of particles in a sample can reduce the signal-to-noise ratio, reduce column lifetime, and increase backpressure in the LC system, potentially causing system failure. Filtration through a microporous membrane is a simple and effective method for

removing particles from a sample. However, particle retention ability is different between different membranes and between different suppliers. As Figure 2 suggests, PTFE filters with polypropylene housing consistently deliver high particle retention.




Figure 2.

Particle retention ability differs between different membranes and between different suppliers. In principle, 100% of filtration membranes should retain particles. To test this hypothesis, microporous membranes from two different vendors (pore size 0.2  $\mu\text{m}$ ) were tested for latex particle retention following filtration of a suspension of 0.3  $\mu\text{m}$  latex particles in water. PTFE=polytetrafluoroethylene; PP=polypropylene; PVDF=polyvinylidene fluoride; PES=polyethylenesulfone.

### Minimize extractables (contaminants) from sample prep device

Extractable impurities can generate interfering peaks in a chromatogram or mass spectrum, making it difficult or impossible to identify or quantify analytes of interest. Therefore, it is important to use a sample preparation device that leaches minimal impurities into the sample.

Though a number of syringe filters are certified as “low-extractable” for use in high performance LC (HPLC), most of those filters are certified using HPLC

coupled to detection of ultraviolet (UV) absorbance. Though this method provides information about the levels of UV-absorbing extractables coming from a filter, this information does not necessarily correlate with data obtained from an MS detector.

Select vendors now validate syringe filters using mass spectrometric analysis of extractables, which provides valuable guidance in choosing an appropriate membrane for filtering your sample.

#### TIP

In general, hydrophilic PTFE syringe filters provide the cleanest samples (with the lowest levels of extractable impurities). Presence of polymeric extractable impurities (such as from polypropylene syringe filters) complicate analysis of small molecular analytes.

	Hydrophilic PTFE	Polypropylene (Vendor A)	Polypropylene (Vendor B)	Nylon (Vendor A)	Nylon (Vendor B)
Reproducibility					
Range of Compatibility with Organic Solvents	Broad	Broad	Broad	Moderate	Moderate
Extractable Level	Low	High	Medium	High	High
Nature of Extractables	MW 100–400 Da	Polymeric	Variable	Polymeric-Variable	Polymeric-Variable

Table 3.

Overall mass spectral signal intensity for five different types of HPLC-certified syringe filters when tested using eight different solvents. The range of chemical compatibility with solvents may indicate the general level of extractables leached by a particular membrane. Millex® LCR filters, which contain hydrophilic PTFE and have broad compatibility with solvents, show the lowest level of signal intensity (and therefore background noise). On the other hand, polypropylene syringe filters from vendor A as well as nylon syringe filters from vendors A and B all show very high levels of extractables, impacting background signal.

Consider these parameters for evaluating the suitability of a membrane filter for LC-MS:

Solvent compatibility of device

When selecting a filter, determine if constituents in the liquid being filtered will chemically attack the filter. If the filter undergoes chemical degradation, its performance will be compromised, and it may release foulants into the sample stream.

Some solvents may be incapable of dissolving the filter, but could be absorbed into the polymer matrix, causing it to swell over time, altering the effective pore size of the filter and changing its performance.

Lot-to-lot reproducibility of extractables level

This parameter reflects the consistency with which filters are manufactured. Since there are very few MS-certified filters, this parameter helps select the right filter for MS applications and indicates the degree of variation in levels of extractables when different lots of syringe filters are used.

Intensity of signal contribution from extractables: Total Ion Current (TIC) chromatograms

LC-MS-certified membrane filters should be supplied with a total ion current chromatogram that shows the intensity of all peaks detected under a specified set of experimental conditions, normalized to an internal standard. The TIC chromatograms can enable comparisons of extractable profiles between membranes and different filter vendors.

Type of extractables: low molecular weight, discrete peaks vs. polymeric peaks

Any type of extractables can confound downstream analysis, but the discrete peaks from low molecular weight extractables are typically less problematic than peaks from polymeric extractables, which typically show peaks separated by a common mass difference ranging over a wide m/z range. (See Appendix II for a table of mass differences of repeating units derived from common contaminating extractables.) Polymeric extractables are also difficult to remove from the sample or mass spectrometer, even after extensive cleaning of the mass spectrometer.

Adsorption of analyte to device

Because the internal surface area of polymeric microporous membranes is 100–600 times as great as the frontal surface area, there is a vast internal surface area available for nonspecific binding.

Choosing a membrane filter with low nonspecific analyte binding ensures that the overall molecular composition of the filtrate is minimally altered upon passing through the device.

Common extractable contaminants

- Polyethylene glycol (PEG)
- Metal ions (e.g., lithium, sodium, potassium, copper, platinum, iron)
- Phthalates (present in many plastics)
- Slip agents (amides)




Figure 3.

Few extractable impurities from Millex® LCR syringe filter (left) containing 0.45 µm pore hydrophilic PTFE membrane as detected by MS. In contrast, a syringe filter containing 0.45 µm pore polypropylene membrane (non-MilliporeSigma, right) shows significant extractables. Presence of polymeric extractable impurities (from polypropylene syringe filters) complicate analysis of small molecular analytes. Millex® LCR filters showed a highest peak intensity of about 8 x 10<sup>5</sup> for extractable masses, whereas non-MilliporeSigma polypropylene syringe filters showed extractable levels about 30 times higher (2.5 x 10<sup>6</sup>). Such high signal intensity, which can be comparable to the signal from the analyte of interest, can make analyte quantitation very challenging.

TIP

Prerinse the filter with sample/solvent to reduce the extractables.




Figure 4.

**Nylon syringe filters are a common source of extractables.** In this experiment, the extractable peaks seen after filtration through nylon were greatly reduced when the first milliliter of filtrate was discarded and the second milliliter was analyzed.

Another potential source of contamination is the syringe used to filter and/or inject the sample. Table 4 shows the level of zinc contamination from various types of syringes.



Syringe Used	Time of Contact	Zn Contamination (ppb)
Plastic with air gap	No contact	< 10
Plastic with black piston seal	15 min.	96
Plastic with black piston seal	30 min.	171
Glass with metal Luer fitting & PTFE piston	30 min.	470

Table 4.

Level of zinc contamination with respect to syringe used.

TIP

Use a plastic syringe with an air gap between the sample and the piston. Any surface that comes in contact with the sample has the potential to introduce extractables as well as contribute to analyte binding.


Ultrafiltration separates free from protein-bound analytes

Centrifugal ultrafilters, particularly devices with regenerated cellulose membrane that have defined nominal molecular weight cutoffs, are ideal for separating free from protein-bound microsolute in serum, plasma, and other biological samples,

as illustrated in Figure 5. This sample preparation method has been cited in LC-MS analyses for deproteinizing samples to reduce complexity or matrix interference. The method has also been used for LC-MS analyses of binding studies in new drug investigations.

TIP

If filtering many samples at a time, increase throughput while maintaining consistency by using 96-well or 384-well filter plates and a microplate-compatible vacuum manifold.



When filtration isn't enough


For more complex sample matrices, use more specific sample preparation methods, such as solvent evaporation, protein precipitation, liquid/liquid extraction, QuEChERS, and SPE to transform samples into forms suited for LC-MS.

Some useful tools for these procedures include:

- Separatory funnel
- LiChrolut® product range – for SPE
- EXTrelut® pre-packed columns for extraction of lipophilic compounds from aqueous solutions – for sorbent-supported LLE workflows
- Supel™-Select Polymeric SPE – HLB and ion-exchange phases for a wide range of applications and pH conditions
- Solvents, acids, bases, salts – for protein precipitation


For samples with high salt load (e.g., food, body fluids or tissue) a desalting (sample preparation step) using LiChrolut® cartridges is recommended.

Sample characteristics determine your SPE procedure



Protein precipitation, followed by filtration, is often an effective, simple way to reduce the complexity of the sample matrix (Figure 6). For this process, it can be advantageous to use a filter plate, which enables


precipitation and filtration in a single device, eliminating the need for sample transfer and thereby improving analyte recovery.



**Figure 6.**  
**Analyte recovery after protein precipitation.** Three different drugs (chlorpheniramine, imipramine and niflumic acid) were spiked into plasma at various concentrations. Protein precipitation was carried out using 1:4 water:acetonitrile as the precipitating solvent. The samples were filtered through various multiwell filter plates with polypropylene, PTFE, and washed PTFE (washed with solvent). Drug recovery in the filtrate was determined using LC-MS/MS analysis of the filtrate.

**TIP**

Prepare samples for nano LC-MS using ZipTip® pipette tips (Figure 7). This sample preparation microdevice is a 10 µL pipette tip with a 0.6 or 0.2 µL bed of chromatography media fixed at its end with no dead volume. It is ideal for concentrating and purifying samples for sensitive analyses such as nano LC-MS or MALDI-ToF MS.



**Figure 7.**  
ZipTip® pipette tips increase sensitivity of mass spectrometric analysis. Plasma sample from rat dosed with 10 mg/kg prazosin was injected into an LTQ/Orbitrap mass spectrometer by nanoelectrospray (a) before and (b) after preparation using a C18 ZipTip® pipette tip. Adapted from Erve JCL et al, Rapid Commun. Mass Spectrom. 2008; 22: 3015–3026.

**Phospholipids: a concern or LC-MS analysis of small molecules in biological matrices**


**Phospholipids are present as a major component of all cell membranes.**

They are therefore present in all biological sample matrices including serum, plasma and whole blood and can be a problem in LC-MS analysis of small molecules because they often co-elute and ionize along with the analytes of interest. This co-elution results in ion suppression (an erroneous decrease) of the MS signal that can cause variability and impact LC-MS result accuracy. Even if phospholipids do not co-elute with the analyte of interest, they can accumulate on your analytical column.

**Phospholipid removal techniques:**

To overcome the problem of phospholipid-induced ion suppression, some analysts try traditional SPE. Traditional SPE often requires time-consuming and complex method development, but still only removes nominal amounts of phospholipids. A variety of products designed specifically for the removal of both proteins and phospholipids are now commercially available, including HybridSPE® plates and cartridges (Figure 8).





**Figure 8.**  
Removing phospholipids can improve the signal-to-noise ratio in LC-MS.


Guide to Sample Preparation Tools

Description	Catalog No.
<b>Millex® Syringe Filters</b>	
Millex®-LCR Filter, 0.45 µm, Hydrophilic PTFE, 25 mm, non-sterile, 50/pk	<b>SLCR025NS</b>
Millex®-LCR Filter, 0.45 µm, Hydrophilic PTFE, 25 mm, non-sterile, 250/pk	<b>SLCR025NB</b>
Millex®-LCR Filter, 0.45 µm, Hydrophilic PTFE, 25 mm, non-sterile, 1000/pk	<b>SLCR025NK</b>
Millex®-LCR Filter, 0.45 µm, PTFE, 13 mm, non-sterile, 100/pk	<b>SLCR013NL</b>
Millex®-LCR Filter, 0.45 µm, PTFE, 13 mm, non-sterile, 1000/pk	<b>SLCR013NK</b>
<b>ZipTip® Pipette Tips</b>	
ZipTip® with 0.6 µL C 4 resin, 8/pk	<b>ZTC04S008</b>
ZipTip® with 0.6 µL C 4 resin, 96/pk	<b>ZTC04S096</b>
ZipTip® with 0.6 µL C 4 resin, 960/pk	<b>ZTC04S960</b>
ZipTip® with 0.2 µL C 18 resin, 8/pk	<b>ZTC18M008</b>
ZipTip® with 0.2 µL C 18 resin, 96/pk	<b>ZTC18M096</b>
ZipTip® with 0.2 µL C 18 resin, 960/pk	<b>ZTC18M960</b>
ZipTip® with 0.6 µL C 18 resin, 8/pk	<b>ZTC18S008</b>
ZipTip® with 0.6 µL C 18 resin, 96/pk	<b>ZTC18S096</b>
ZipTip® with 0.6 µL C 18 resin, 960/pk	<b>ZTC18S960</b>
ZipTip® with 0.6 µL strong cation resin, 8/pk	<b>ZTSCXS008</b>
ZipTip® with 0.6 µL strong cation resin, 96/pk	<b>ZTSCXS096</b>

Description	Catalog No.
<b>Samplicity® Filtration Systems</b>	
Millex Samplicity® Filters, 0.45 µm Hydrophilic PVDF, 96/pk	<b>SAMPHV001</b>
Millex Samplicity® Filters, 0.45 µm Hydrophilic PVDF, 384/pk	<b>SAMPHV004</b>
Millex Samplicity® Filters, 0.45 µm Hydrophilic PTFE, 96/pk	<b>SAMPLCR01</b>
Millex Samplicity® Filters, 0.45 µm Hydrophilic PTFE, 384/pk	<b>SAMPLCR04</b>
Millex Samplicity® Filters, 0.20 µm Hydrophilic PTFE, 96/pk	<b>SAMPLG001</b>
Millex Samplicity® Filters, 0.20 µm Hydrophilic PTFE, 384/pk	<b>SAMPLG004</b>
Samplicity® Filtration System, Bold Blue	<b>SAMPSYSBL</b>
<b>Ultrafree®-MC and -CL Centrifugal Microfiltration Units</b>	
Ultrafree®-MC Filter, 0.22 µm Hydrophilic PTFE, 25/pk	<b>UFC30LG25</b>
Ultrafree®-MC Filter, 0.45 µm Hydrophilic PTFE, 25/pk	<b>UFC30LH25</b>
Ultrafree®-CL Filter, 0.22 µm Hydrophilic PTFE, 25/pk	<b>UFC40LG25</b>
Ultrafree®-CL Filter, 0.45 µm Hydrophilic PTFE, 25/pk	<b>UFC40LH25</b>
Centrifree® Ultrafiltration Device with Ultracel® Membrane	<b>4014</b>
<b>MultiScreen® Filter Plates</b>	
MultiScreen® Solvinert 96-well Plate, 0.45 µm Hydrophilic PTFE, 50/pk	<b>MSRLN0450</b>
MultiScreen® Deep Well Solvinert 96-well Plate, 0.45 µm Hydrophilic PTFE, 10/pk	<b>MDRLN0410</b>
<b>Solid Phase Extraction</b>	
EXtrelut® NT 20 pre-packed columns for extraction of lipophilic compounds from aqueous solutions (20 mL sample)	<b>115096</b>
LiChrolut® RP-18 E (40 – 63 µm) 500 mg 3 mL standard PP-tubes 50 extraction tubes per package	<b>119849</b>

Proper mobile phase preparation to minimize contaminants

For LC-MS, use the highest quality of pure solvents and reagents and avoid further contamination by careful handling. Any impurity could cause signal suppression and/or adduct formation with target molecules and therefore decrease sensitivity (signal-to-noise ratio) and/or increase complexity of the mass spectrum.



Chromatographic conditions	
<b>System</b>	Bruker Esquire 3000+ ion trap MS
<b>Detection</b>	Pos. ESI-MS, m/z range 50 – 2000
<b>Flow rate</b>	1.0 mL/min
<b>Mobile phase</b>	A: Water (Cat. No. <a href="#">115333</a> ) B: Various ACN qualities as indicated in the graphics
<b>Gradient</b>	0 min 98% A, 30 min 0% A, 38 min 98% A
<b>Temperature</b>	25°C

**Figure 9.**  
Combined TICs of the blank runs of four different acetonitrile qualities. All solvents were delivered to the MS source via an LC system.

## Hypergrade and gradient grade solvents minimize contaminant peaks

Figure 9 illustrates the influence of LiChrosolv® acetonitrile quality on the background noise intensity in mass spectra. MilliporeSigma solvents labeled “hypergrade for LC-MS LiChrosolv®” are dedicated for use with MS systems and deliver minimized

contaminant peaks, ion suppression, adduct formation and background noise and therefore maximize sensitivity. Gradient grade solvent quality (labeled “gradient grade for liquid chromatography LiChrosolv®”) are suitable for LC-UV gradient runs.


## Use ultrapure water (bottled or freshly delivered)

Ultrapure water for LC-MS applications can be either bottled or freshly delivered from a water purification system. The choice is mainly determined by daily consumption. Demineralized tap water is not recommended for use in combination with LC-MS setups because of possible system contamination. The quality of LiChrosolv® bottled water for chromatography and fresh Milli-Q® ultrapure water produced from laboratory water purification systems is consistently high and generally independent of the regularity of use.

As a result, the now-improved water quality is perfectly suitable for the production of mobile phase, buffers, blanks, standards preparation, sample dilution, glassware rinsing or extraction used in these critical applications.

### TIP

If sharing a Milli-Q® ultrapure water system with molecular biology researchers requiring a BioPak® application-specific polisher at the point of use, we recommend adding a separate point-of-delivery (POD) unit fitted with an LC-Pak® polisher, because water passing through a BioPak® unit may not be suitable for LC-MS applications. The LC-Pak® polisher is optimized for mobile phase, buffers, blanks, standards preparation, sample dilution, glassware rinsing or extraction used in LC-MS.



**Figure 10.**

Exploiting the full potential of Milli-Q® systems via proper handling. Top: MS spectra of two samples of water delivered at different points of time; bottom: TICs of the same samples and one additional sample. All analyses were performed via direct injection of the preconcentrated solvents into the MS operated in positive ESI mode.

## Solvent storage

### Tips for maintaining your solvent purity:

- Store all eluents (water and organic) in surface-treated amber glass bottles (original packaging of all MilliporeSigma LC-MS grade solvents) or in borosilicate glass (if solvents have to be decanted).
- Select a solvent storage system that is appropriate for usage volume and withdrawal frequency (Table 5).
- Do not use standard glass bottles; silica and alkali dissolve and form adducts  $[M+X]^+$  with analytes.
- Use MilliporeSigma HPLC bottle caps/adapters with tube connections and membrane filter mounted directly on the original brown glass bottle. This protects both solvents and environment.
- Avoid decanting; it is a possible source of contamination.
- Avoid improvised repairs for fixing solvent tubing; this may cause leakages and/or release of contaminants to the eluents.
- Do not use plastic devices (bottles, funnels, etc.) to handle or store solvents, buffers, etc. Solvents extract additives (anti-static agents, stabilizers, plasticizers) from plastic, a source of contaminant ghost peaks and increased background noise.

Solvent storage systems	Storage container volume
Bottle top adapters for directly connecting solvent bottles to LC system	1 L, 2.5 L, 4 L (for infrequent withdrawal)
Stainless steel barrel with adapter for direct withdrawal	10 L, 30 L (for frequent withdrawal)
Stainless steel barrel directly connected to LC system	
Central storage of stainless steel barrels with adapters to supply solvent to multiple different laboratories	

**Table 5.**  
Select solvent storage options based on volume of usage and frequency of withdrawal to minimize contamination.

## Using water as a mobile phase?

### Keep in mind these additional considerations:

- Keep 5% organic solvent in your eluent if chromatographic conditions allow. This avoids microbial contamination of bottle, tubing and LC system.
- Keep 5% of aqueous eluent in the organic mobile phase to avoid buffer precipitation in the system, e.g., in valves, and subsequent tedious cleaning procedures.

## Solvent container cleaning: avoid the dishwasher


Dishwashers are standard laboratory equipment, but they are operated using chemicals such as strong bases and surfactants. Strong bases can lead to dissolution of silica and alkali from glassware and cause the formation of adducts  $[M+X]^+$  with analytes, while traces of surfactants remain on the glass surface after the cleaning process and decrease MS sensitivity by increasing background noise.

The easiest way to avoid dishwashing is “cleaning” of all equipment via simple evaporation of both solvents and additives. All chemicals dedicated to the application in LC-MS are volatile; therefore, this procedure is straightforward as long as chemicals are highly pure and microbial growth can be eliminated. In case of equipment contamination, flushing with LiChrosolv® or Milli-Q® water or organic hypergrade solvents has to be performed to achieve sustainable cleaning.

## Buffers and additives

### When working with buffer to adjust pH of eluents, keep in mind:

- Use volatile salts (such as ammonium formate, ammonium acetate, or trimethylamine). Nonvolatile salts (e.g., phosphates, borates, sulfates or citrates) precipitate in and block the MS source, requiring tedious cleaning procedures.
- Total ionic strength of the eluent should not exceed 20 mM. Adjust buffer concentration in the aqueous solvent accordingly. Buffers for LC-MS should be prepared using the purest salt and acid/base quality available. If possible, avoid working with an ammonium bicarbonate buffer. The salt is often highly contaminated — see comparison with ammonium acetate (Figure 11).



**Figure 11.**  
**Comparison of MS spectra (top) and TIC chromatograms (bottom) of the two buffer systems, ammonium bicarbonate and ammonium acetate.** Both mixtures were prepared using Milli-Q® ultrapure water and the same acetic acid source and were analyzed via direct injection into the mass spectrometer operated in positive ESI mode. Note the high MS background signals observed when utilizing ammonium bicarbonate as a buffer.

Buffer pH is generally adjusted via a titration with the respective acid or base and monitored with a pH electrode. The unavoidable contamination of the buffer solution with alkali ions from the pH electrode can be decreased by using a miniaturized system available from several suppliers. Unlike standard equipment with a diameter of approximately 10 mm, the diameter of miniaturized electrodes is only 3 mm.

Buffers not only adjust the pH and ionize a target molecule [M], they can also form adducts [M+buffer], e.g., with ammonium, alkali, halogens, formate or acetate. This leads to the detection of additional peaks in the MS spectrum. Even a complete suppression of the analyte signal is possible when the vapor pressure of the resulting adduct (mainly alkali) is decreased significantly. As a result of this phenomenon and in order to keep the ESI source clean, volatile buffers are recommended.

TIP

Avoid using TFA. Trifluoroacetic acid (TFA) is widely used as an ion pairing reagent to improve the liquid chromatographic separation of peptides or proteins when using standard UV for detection. However, TFA can cause strong ion suppression in mass spectrometry (mainly in negative ESI mode) and also contaminates the LC-MS system. Formic acid (0.1%) is commonly used instead as a mobile phase modifier that is compatible with LC-MS.

Guide to Mobile Phase Preparation Reagents


Description	Catalog No.
Milli-Q® Advantage Water Purification System	Z00Q0V0WW*
Milli-Q® Integral Water Purification System	ZRXQ010WW*
LC-Pak® Application-Specific Polisher	LCPAK0001
LiChrosolv® Solvents	
Acetonitrile hypergrade for LC-MS LiChrosolv®	100029
Methanol hypergrade for LC-MS LiChrosolv®	106035
Ethanol gradient grade for liquid chromatography LiChrosolv®	
2-Propanol gradient grade for liquid chromatography LiChrosolv®	101040
Toluene for liquid chromatography LiChrosolv®	108327
Water for chromatography LiChrosolv® (LC-MS)	115333
Suprapur® Inorganic Acids and Bases	
Acetic acid (glacial) 100% Suprapur®	100066
Ammonia solution 25% Suprapur®	105428
Formic acid 98–100% Suprapur®	111670
Hydrochloric acid 30% Suprapur®	100318
Other Reagents	
Acetic acid (glacial) 100% anhydrous for analysis EMSURE® ACS, ISO, Reag. Ph Eur	100063
Ammonia solution 28–30% for analysis EMSURE® ACS, Reag. Ph Eur	105423
Ammonium acetate for analysis EMSURE® ACS, Reag. Ph Eur	101116
Dichloromethane for organic trace analysis UniSolv®	106454
Formic acid 98–100% for analysis EMSURE® ACS, Reag. Ph Eur	100264
n-Hexane for organic trace analysis UniSolv®	104369
n-Pentane for organic trace analysis UniSolv®	107288
Petroleum benzine boiling range 40–60°C for organic trace analysis UniSolv®	116740
2-Propanol for analysis EMSURE® ACS, ISO, Reag. Ph Eur	109634

\*Contact your local representative for detailed ordering information.

For even more LC-MS solvents and reagents, visit: [SigmaAldrich.com/lc-ms](https://sigmaaldrich.com/lc-ms)

TIP

Avoid equilibrating columns with more than 10 column volumes of mobile phase (or one blank gradient run with subsequent equilibration). Contaminants in solvents and additives can accumulate on a stationary phase. Figure 12 shows this effect for plasticizers dissolved in the eluent on a reversed phase column after equilibration for 0, 15 and 60 minutes. While these compounds would become eluted as very broad peaks under isocratic conditions (and cause an increased background noise), they elute as distinct, intensive peaks under gradient conditions and can interfere with analyte signals. Instead, run samples immediately after two or three blank runs to ensure that the system is stable prior to sample analysis.



Chromatographic conditions	
System	Bruker Esquire 6000plus
Detection	Pos. ESI-MS, TICs
Flow rate	0.4 mL/min
Mobile phase	A: Water (Cat. No. <a href="#">115333</a> ) / acetonitrile (Cat. No. <a href="#">100029</a> ) 95/5 (v/v) + 0.1 % formic acid (Cat. No. <a href="#">100264</a> ) B: Acetonitrile + 0.1% formic acid
Gradient	0 min 100% A, 3 min to 95% B, 5 min at 95% B, then back to 100% A
Temperature	25°C
Sample	Plasticizers (*) were added by the immersion of plastic tubing in aqueous solvent B

**Figure 12.** Accumulation of contaminants on an HPLC column for various periods of time and elution via a gradient profile.

# Effects of column choice on LC-MS performance

For samples available in small amounts (such as plasma or serum) or where analyte concentrations are low, a setup consisting of both highly sensitive separation and MS detection techniques is necessary for proper identification of the target molecules.

Suboptimal column choice and misuse of the column could decrease the signal-to-noise ratio and increase background noise.

## Tips for proper column choice:

- 1. Refer to comprehensive column selection guides** for full guidance on selecting a column with optimal stationary phase and dimensions and to match method specifications. Column selection guides can be found at:
  - [MerckMillipore.com/chromatography](https://www.merckmillipore.com/chromatography)
  - [SigmaAldrich.com/choose-column](https://www.sigmaaldrich.com/choose-column)
  - The technical resource, “Solutions for Mass Spectrometry,” which is available at: [MerckMillipore.com/mass-spectrometry](https://www.merckmillipore.com/mass-spectrometry)
- 2. Pick column in accordance with eluent pH.** Using an eluent pH that is too high (e.g., >8) can dissolve the backbone of silica-based HPLC columns. Using a pH that is too low (<2) can strip the stationary phase (C18, etc.). Both options can lead to additional signals in your spectrum, increased background noise and/or signal suppression. Both scenarios may decrease column lifetime.
- 3. Consider using polymeric columns for highly alkaline samples.** Polymeric columns are more stable at high pH than silica columns, where the silica may dissolve. However, polymeric columns possess smaller phase ratios and therefore, lower resolution. In addition, they are prone to swelling in organic solvents, leading to changed chromatographic characteristics. Furthermore, due to micropores in the stationary phase, column performance may be lower as compared to silica-based columns.
- 4. Highly endcapped stationary phases** are another good option for sample analysis at high pH. Endcapping or a cross-linked C18 modification leads to more pH-stable columns.

- 5. Column diameter influences the sensitivity of the analysis.** The sensitivity increases with decreasing column internal diameter (or increasing mass of the injected sample). For example, when changing from a 4.6 mm i.d. column to a 0.1 mm i.d. capillary column, sensitivity theoretically increases by a factor of approximately 2000 (Table 6). Hence, a combination of capillary chromatography coupled to mass spectrometry may be the best combination for high sensitivity analysis. However, it is important to note that extra-column effects may impact signal-to-noise ratio when column diameter is decreased—for example, the system dwell volume, dead volume in the system, the ability of the system pump to deliver accurate gradient and the volume of the detector cell can all result in peak broadening and loss of sensitivity.


Column i.d. (mm)	Typical flow rate (µL/min)	Relative sensitivity
4.6	1000 – 6000	1
2.0	200 – 800	5.3
0.2	0.5 – 20	530
0.1	0.4 – 3	2100
0.05	0.1 – 0.8	8500

**Table 6.**  
Effect of decreasing column diameter on flow rate and relative sensitivity.

- 6. If the sample and analyte allow for HILIC chromatography,** consider using HILIC instead of reversed phase columns. In HILIC chromatography, analysis is performed under highly organic conditions (e.g., 80% acetonitrile and 20% aqueous buffer) and polar compounds elute later than in reverse phase chromatography. Under these conditions, the eluent is vaporized more easily in the MS source, resulting in better sensitivity (signal-to-noise ratio).

## Tips for column usage

1. Wash columns after each use with appropriate strong eluent to remove all adsorbed compounds. (See “Column bleeding” and Figure 13.)
2. Use the column at proper operating temperature (refer to each column User Guide) in order to avoid loss of stationary phase or dissolution of column backbone, both of which contribute to the appearance of additional signals on the spectrum that may interfere with analysis.
3. You can use a guard column (usually either 5 mm or 10 mm long) directly in front of the main column to protect the column against contamination (such as particles, sample matrix) or harsh eluent pH. Guard columns should be changed frequently in order to keep system backpressure low and in order to maximize the lifetime of the analytical column.




**Figure 13.**  
Decreasing column diameter may improve sensitivity. Typical fast separation of four compounds in less than two minutes using a Chromolith® 4.6 mm i.d. column at a flow rate of 4 mL/min (left). The same separation was achieved on a Chromolith® 3 mm i.d. column (right). Both chromatograms exhibit excellent column efficiency and peak resolution, however the 3 mm i.d. column demonstrates improved sensitivity at just 1.7 mL/min, thus saving 57% of solvents.



Column bleeding

The stationary phase of every HPLC column (except for normal phase systems) is made out of covalently bound organic entities altering its physical properties. Depending on the quality of both phase modification and a subsequent washing step, these entities (e.g., octadecyl, cyano, phenyl) can be stripped off the column during a chromatographic run and cause weak to severe interfering signals.

This unwanted phenomenon is referred to as “column bleeding” and leads to a decreased sensitivity in MS. It can be avoided by flushing the column prior to analysis using isopropanol and 0.1% formic acid as a solvent at half optimum flow for one hour. This process removes unbound or weakly bound organic entities, minimizes column bleeding and hence increases sensitivity by decreasing background noise (Figure 14).



**Figure 14.** Column washing can compensate for column bleeding. Total ion current (TIC) of a competitor column after 20 gradient runs: left– unwashed column; right– column washed with 0.1% formic acid in isopropanol at 0.2 mL/min for 48 min.

Overview of common column types for LC-MS samples

1. Reversed phase, for hydrophobic molecules

- a. Purospher® STAR RP-18 endcapped (standard)
  - i. Small particles deliver high separation efficiency/peak capacity.
  - ii. Suitable for cleaner samples and MS, after removal of matrix components.
- b. Use C8 instead of C18 column to separate lipophilic analytes from matrix.
- c. Fused-core columns; for example, Ascentis® Express columns, which feature narrower particle size distribution, more consistent packed bed, and shorter diffusion path compared to traditional, fully porous particles. The result is increased resolution, added sensitivity and faster runs.

2. HILIC (ZIC®-HILIC/cHILIC [small molecules])

- a. SeQuant® ZIC®-HILIC/cHILIC/pHILIC bonded zwitterionic stationary phases
  - i. For polar hydrophilic molecules, i.e., the majority of the endogenous molecules.
  - ii. Combines perfectly with ESI-MS detection due to the applied solvents and additives.
  - iii. Significant increase in sensitivity in comparison with reversed phase chromatography.
- b. Strongly retained polar analytes can be removed from HILIC columns by changing to a more polar eluent.

Monolithic silica columns have high matrix tolerance

The analysis of samples with high matrix load requires tedious and time-consuming sample preparation steps. For cost-effective investigations, sample handling has to be kept as short as possible and combined with robust LC columns displaying a high matrix tolerance and long lifetime.

The 50-2 mm monolithic silica column is well-suited for fast gradient run liquid chromatography, and the applied low flow rates make it the perfect choice for MS detection. Analysis of dirty samples, such as food or tissue, can be performed on this robust column type without the need for a guard column or tedious and complex sample preparation procedures.

The advantages of Chromolith® monolithic silica columns are:

- Exceptional robustness or lifetime—described as number of injections—enabling cost savings.
- High matrix tolerance decreases tedious sample preparation steps, speeds up all processes and allows for fast and simple HPLC analyses.
- Low backpressure, fast analytical speed and high reproducibility on standard HPLC systems without the need for upgrading to high-priced UHPLC solutions.

Consult a column selection guide today!

[MerckMillipore.com/mass-spectrometry](https://www.merckmillipore.com/mass-spectrometry)



# Useful LC-MS resources

## Solutions for Mass Spectrometry

Since the development of an aluminum oxide for adsorption chromatography in 1904, MilliporeSigma has offered chromatography products and solutions of exceptional, application-oriented quality. Our columns, plates and mobile phases are used in mass spectrometry labs throughout the world.

This compendium of selected articles is intended to help and support users of chromatographic separation methods in mass spectrometric laboratories worldwide, by providing an overview of different actual topics in this field of analytical science.

Download at:  
[MerckMillipore.com/mass-spectrometry](https://www.merckmillipore.com/mass-spectrometry)

## John Dolan’s LC Troubleshooting Bible (2014 edition)

For over 30 years, John Dolan has been writing the LC Troubleshooting column of LC/GC Magazine, answering questions from chromatographers around the world. Browse the columns all in one place (1984 – 2013) at the LC Troubleshooting Bible website.

Questions addressed cover all aspects of chromatography, ranging from “Importance of Sample Matrix” to “Gradient Elution Surprises.”

Visit: [lcsources.com/tsbible](https://lcsources.com/tsbible)

## The ChromBook™ App

The MilliporeSigma ChromBook™ App is more than a mobile device-friendly catalog of sorbents, analytical and preparative columns, sample preparation and thin layer chromatography products; it also includes tips for LC method development and a column selection guide based on analyte properties.

Download at:  
[MerckMillipore.com/mass-spectrometry](https://www.merckmillipore.com/mass-spectrometry)

## chromacademy.com

Described as the “world’s largest e-Learning website for analytical scientists,” CHROMacademy is a well-organized, all-in-one-place resource for basic LC/LC-MS training, troubleshooting, learning about advanced topics and even obtaining technical support from consultants. It includes video guides and quick tips, such as in “My LC-MS Isn’t Behaving!”

Many features are available by subscription only, but for analysts doing a lot of LC, it’s worth it.

### Training modules within CHROMacademy include:

- Sample Prep
- HPLC
- MS

## Sigma-Aldrich® LC-MS Resource Guide

A veteran of separation science, Sigma-Aldrich® has published this practical guide specifically for LC-MS users analyzing biomolecules, food and beverage, environmental samples and pharmaceutical products. It covers many advanced technologies, including enantiomer separation, phospholipid removal, solid phase microextraction (SPME) and molecularly imprinted polymers (MIP).

Download at:  
[SigmaAldrich.com/lc-ms](https://www.sigmaaldrich.com/lc-ms)

## Virtual Conference: 5-Star LC-MS Practices

Held in 2014, the conference addressed the need for the use of high-quality tools in LC-MS, including five major types of materials: ultrapure water, solvents, reagents, filters, and liquid chromatography columns. During the conference, MilliporeSigma experts from around the world presented a series of webinars focusing on LC-MS best practices for their respective areas. Virtual conference participants learned why the use of “5-star” materials is critical for their LC-MS work, and how they can maximize their research using a variety of tips and tricks explained in the webinars.

View on demand at:  
[MerckMillipore.com/mass-spectrometry](https://www.merckmillipore.com/mass-spectrometry)

# Appendix I.

## Common mass spectrometry contaminants and their sources

This list of potential interfering or contaminant ions in mass spectrometry (ESI positive mode, mass ≤1000 Da) is adapted from an excerpt of “Interferences and contaminants encountered in modern mass spectrometry” Bernd O. Keller, Jie Sui, Alex B. Young and Randy M. Whittall Analytica Chimica Acta 627, Issue 1, 3 October 2008, Pages 71-81. However, we have updated the masses listed in the previous publication by calculating the singly charged monoisotopic ion mass of each listed ion based on its molecular formula.

Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
33.0340	[M+H] <sup>+</sup>	CH <sub>3</sub> OH	Methanol	Acetonitrile, solvent
42.0344	[M+H] <sup>+</sup>	CH <sub>3</sub> CN	ACN	Acetonitrile, solvent
59.0609	[M+NH <sub>4</sub> ] <sup>+</sup>	CH <sub>3</sub> CN	ACN	Acetonitrile, solvent
63.0446	[A <sub>1</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
64.0163	[M+Na] <sup>+</sup>	CH <sub>3</sub> CN	ACN	Acetonitrile, solvent
65.0603	[M <sub>2</sub> +H] <sup>+</sup>	CH <sub>3</sub> OH	Methanol	Methanol, solvent
74.0606	[M+H] <sup>+</sup>	C <sub>3</sub> H <sub>7</sub> NO	Dimethyl formamide	solvent
74.0606	[A <sub>1</sub> B <sub>1</sub> +H] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (CH <sub>3</sub> OH) <sub>m</sub>	Acetonitrile/Methanol	ESI solvents
77.0603	[A <sub>1</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
79.0218	[M+H] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
83.0609	[M <sub>2</sub> +H] <sup>+</sup>	CH <sub>3</sub> CN	Acetonitrile	ESI solvents
85.0265	[A <sub>1</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
85.0594	[M+H] <sup>+</sup>	C <sub>2</sub> D <sub>6</sub> OS	d6-DMSO	d6-Dimethylsulfoxide, solvent
88.0399	[A <sub>1</sub> B <sub>1</sub> +H] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (HCOOH) <sub>m</sub>	Acetonitrile/Formic Acid	ESI solvents
96.0425	[A <sub>1</sub> B <sub>1</sub> +Na] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (CH <sub>3</sub> OH) <sub>m</sub>	Acetonitrile/Methanol	ESI solvents
99.0422	[A <sub>1</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
101.0841	[M+H] <sup>+</sup>	C <sub>5</sub> H <sub>10</sub> NO	NMP	N-methyl 2-pyrrolidone; solvent, floor stripper
101.0005	[A <sub>1</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
101.0037	[M+Na] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
101.0814	[A <sub>2</sub> B <sub>2</sub> +H] <sup>+</sup>	[MeOH] <sub>n</sub> [H <sub>2</sub> O] <sub>m</sub>	Methanol/Water	ESI solvents
102.0555	[A <sub>1</sub> B <sub>1</sub> +H] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (CH <sub>3</sub> COOH) <sub>m</sub>	Acetonitrile/Acetic Acid	ESI solvents
102.1283	[M+H] <sup>+</sup>	C <sub>6</sub> H <sub>15</sub> N	TEA	Triethylamine, buffer
103.9561	[M+ <sub>63</sub> Cu] <sup>+</sup>	C <sub>2</sub> H <sub>3</sub> N	ACN	Acetonitrile, solvent
104.9928	[M+Na] <sup>+</sup>	C <sub>2</sub> H <sub>3</sub> O <sub>2</sub> Na	Sodium acetate	ESI solvents
105.0429	[M <sub>2</sub> +Na] <sup>+</sup>	C <sub>2</sub> H <sub>3</sub> N	ACN	Acetonitrile, solvent
41.0265	[M+ <sub>65</sub> Cu] <sup>+</sup>	C <sub>2</sub> H <sub>3</sub> N	ACN	Acetonitrile, solvent
107.0708	[A <sub>2</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
115.0161	[A <sub>1</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
115.0871	[A <sub>1</sub> B <sub>1</sub> +H] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (C <sub>3</sub> H <sub>7</sub> NO) <sub>m</sub>	Acetonitrile/Dimethylformamide	solvent
120.0483	[M+CH <sub>3</sub> CN+H] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
122.0817	[M+H] <sup>+</sup>	C <sub>4</sub> H <sub>11</sub> NO <sub>3</sub>	TRIS	TRIS, buffer
123.0633	[A <sub>2</sub> B <sub>2</sub> +Na] <sup>+</sup>	[CH <sub>3</sub> OH] <sub>n</sub> [H <sub>2</sub> O] <sub>m</sub>	Methanol/Water	ESI solvents
123.0922	[M+H] <sup>+</sup>	C <sub>7</sub> H <sub>10</sub> N <sub>2</sub>	DMAP	Dimethylaminopyridine, solvent
124.0374	[A <sub>1</sub> B <sub>1</sub> +Na] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (CH <sub>3</sub> COOH) <sub>m</sub>	Acetonitrile/Acetic Acid	ESI solvents
129.0528	[A <sub>2</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
130.1596	[M+H] <sup>+</sup>	C <sub>8</sub> H <sub>19</sub> N	DIPEA	Diisopropylethylamine, solvent
132.9054	M <sup>+</sup>	Cs	Cs-133	Cesium, from Cesium Iodide used as calibrant
133.1076	[A <sub>3</sub> B <sub>2</sub> +H] <sup>+</sup>	[CH <sub>3</sub> OH] <sub>n</sub> [H <sub>2</sub> O] <sub>m</sub>	Methanol/Water	ESI solvents
135.1021	[A <sub>2</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
137.0749	[M+CH <sub>3</sub> CN+NH <sub>4</sub> ] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
142.0303	[M+CH <sub>3</sub> CN+Na] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
144.1752	[M+H] <sup>+</sup>	C <sub>9</sub> H <sub>21</sub> N	TPA	Tripropylamine, solvent

Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
144.9827	[M <sub>2</sub> + <sub>63</sub> Cu] <sup>+</sup>	CH <sub>3</sub> CN	ACN	Acetonitrile, solvent, together with m/z 147
145.0267	[A <sub>2</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
146.0694	[M <sub>3</sub> +Na] <sup>+</sup>	CH <sub>3</sub> CN	ACN	Acetonitrile, solvent
146.9809	[M <sub>2</sub> + <sub>65</sub> Cu] <sup>+</sup>	CH <sub>3</sub> CN	ACN	Acetonitrile, solvent, together with m/z 145
147.1134	[A <sub>2</sub> B <sub>2</sub> +H] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (CH <sub>3</sub> OH) <sub>m</sub>	Acetonitrile/Methanol	ESI solvents
149.0239	[f+H] <sup>+</sup>	C <sub>8</sub> H <sub>4</sub> O <sub>3</sub>	Pthalic Anhydride	fragment ion originating from phthalate esters
150.1283	[M+H] <sup>+</sup>	C <sub>10</sub> H <sub>15</sub> N	Phenyldiethylamine	solvent
151.0970	[A <sub>3</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
153.1392	[M+H] <sup>+</sup>	C <sub>9</sub> H <sub>16</sub> N <sub>2</sub>	DBU	1,8-Diazabicyclo[5.4.0]undec-7-ene
155.0895	[A <sub>3</sub> B <sub>2</sub> +Na] <sup>+</sup>	[CH <sub>3</sub> OH] <sub>n</sub> [H <sub>2</sub> O] <sub>m</sub>	Methanol/Water	ESI solvents
157.0357	[M <sub>2</sub> +H] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
157.0841	[A <sub>2</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
158.9646	[M+Na] <sup>+</sup>	C <sub>2</sub> F <sub>3</sub> O <sub>2</sub> Na	NaTFA	Sodium trifluoroacetate, salt
163.0395	[M-CH <sub>3</sub> OH+H] <sup>+</sup>	C <sub>10</sub> H <sub>10</sub> O <sub>4</sub>	Dimethyl phthalate	Phthalate esters, plasticizer
163.1334	[M+H] <sup>+</sup>	C <sub>8</sub> H <sub>18</sub> O <sub>3</sub>	DGBE	Diethylene glycol monobutyl ether, cpd. In scintillation cocktail
169.0953	[A <sub>2</sub> B <sub>2</sub> +Na] <sup>+</sup>	(CH <sub>3</sub> CN) <sub>n</sub> (CH <sub>3</sub> OH) <sub>m</sub>	Acetonitrile/Methanol	ESI solvents
170.1188	[M <sub>2</sub> +H] <sup>+</sup>	C <sub>2</sub> D <sub>6</sub> OS	d6-DMSO	d6-Dimethylsulfoxide, solvent
171.0058	[f+Na] <sup>+</sup>	C <sub>8</sub> H <sub>4</sub> O <sub>3</sub>	Phthalic anhydride	from phthalate esters, plasticizer
173.0580	[A <sub>2</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
173.0790	[A <sub>3</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
179.0176	[M <sub>2</sub> +Na] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
181.1229	[M+H] <sup>+</sup>	C <sub>11</sub> H <sub>16</sub> O <sub>2</sub>	BHA	Butylated hydroxyanisole, antioxidant additives
183.0810	[M+H] <sup>+</sup>	C <sub>13</sub> H <sub>10</sub> O	DPK	Diphenyl ketone
183.1444	[A <sub>4</sub> B <sub>2</sub> +H] <sup>+</sup>	[CH <sub>3</sub> OH] <sub>n</sub> [H <sub>2</sub> O] <sub>m</sub>	Methanol/Water	ESI solvents
185.1154	[M+Na] <sup>+</sup>	C <sub>8</sub> H <sub>18</sub> O <sub>3</sub>	GE	glycol ether
186.2222	[M+H] <sup>+</sup>	C <sub>12</sub> H <sub>27</sub> N	TBA	Tributylamine, solvent
189.0529	[A <sub>3</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
193.1440	[A <sub>3</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
195.0657	[M+H] <sup>+</sup>	C <sub>10</sub> H <sub>10</sub> O <sub>4</sub>	Dimethyl phthalate	Phthalate esters, plasticizer
195.1232	[A <sub>4</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
203.1048	[M+Na] <sup>+</sup>	C <sub>11</sub> H <sub>16</sub> O <sub>2</sub>	BHA	Butylated hydroxyanisole, antioxidant additives
205.1263	[A <sub>4</sub> B <sub>3</sub> +Na] <sup>+</sup>	[CH <sub>3</sub> OH] <sub>n</sub> [H <sub>2</sub> O] <sub>m</sub>	Methanol/Water	ESI solvents
214.0902	[M+H] <sup>+</sup>	C <sub>10</sub> H <sub>15</sub> NO <sub>2</sub> S	n-BBS	n-butyl benzenesulfonamide, plasticizer
215.1259	[A <sub>3</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
217.1052	[A <sub>4</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
221.1905	[M+H] <sup>+</sup>	C <sub>15</sub> H <sub>24</sub> O	BTH	Butylated hydroxytoluene, Antioxidant
225.1967	[M+H] <sup>+</sup>	C <sub>13</sub> H <sub>24</sub> N <sub>2</sub> O	DCU	N,N'-Dicyclohexylurea
231.0999	[A <sub>3</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
231.1167	[M+NH <sub>4</sub> ] <sup>+</sup>	C <sub>10</sub> H <sub>15</sub> NO <sub>2</sub> S	n-BBS	n-butyl benzenesulfonamide, plasticizer
233.0791	[A <sub>4</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
236.0721	[M+Na] <sup>+</sup>	C <sub>10</sub> H <sub>15</sub> NO <sub>2</sub> S	n-BBS	n-butyl benzenesulfonamide, plasticizer
239.1495	[A <sub>5</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
239.2254	[(M.H <sub>35</sub> Cl) <sub>2</sub> -Cl] <sup>+</sup>	C <sub>6</sub> H <sub>15</sub> N	TEA.HCl	Triethylamine-hydrochloride, buffer
241.2224	[(M.H <sub>37</sub> Cl) <sub>2</sub> -Cl] <sup>+</sup>	C <sub>6</sub> H <sub>15</sub> N	TEA.HCl	Triethylamine-hydrochloride, buffer
242.2848	M <sup>+</sup>	C <sub>16</sub> H <sub>36</sub> N	TBA	Tetrabutylammonium, buffer
243.1174	M <sup>+</sup>	C <sub>19</sub> H <sub>15</sub>	Trityl cation	Trityl cation, [Ph3C] <sup>+</sup>
243.1725	[M+Na] <sup>+</sup>	C <sub>15</sub> H <sub>24</sub> O	BTH	Butylated hytroxytoluene, Antioxidant additives
251.1858	[A <sub>4</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
251.2011	[AB <sub>1</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
257.0316	[M <sub>3</sub> +Na] <sup>+</sup>	C <sub>2</sub> H <sub>6</sub> OS	DMSO	Dimethylsulfoxide, solvent
261.1314	[A <sub>5</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
265.2168	[AB <sub>1</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
267.1725	[M+H] <sup>+</sup>	C <sub>12</sub> H <sub>27</sub> O <sub>4</sub> P	TBP	Tributylphosphate
273.1279	M <sup>+</sup>	C <sub>20</sub> H <sub>17</sub> O	MMT	Monomethoxytrityl cation
273.1678	[A <sub>4</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
273.1830	[AB <sub>1</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
277.1053	[A <sub>5</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
279.0939	[M+H] <sup>+</sup>	C <sub>18</sub> H <sub>15</sub> OP	TPO	Triphenylphosphine oxide
279.1596	[M+H] <sup>+</sup>	C <sub>16</sub> H <sub>22</sub> O <sub>4</sub>	Dibutylphthalate	Plasticiser, phtalate ester
279.2300	[AB <sub>1</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
281.0517	[M+H-CH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>4</sub>	Polysiloxane	Polysiloxane, (neutral methane loss from m/z 297)
282.2797	[M+H] <sup>+</sup>	C <sub>18</sub> H <sub>35</sub> NO	Oleamide	Slip agent in polyethylene films
283.1757	[A <sub>6</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
284.2953	[M+H] <sup>+</sup>	C <sub>18</sub> H <sub>37</sub> NO	Stearamide	Slip agent in polyethylene films

Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
287.1987	[AB <sub>1</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
288.2539	[M+H] <sup>+</sup>	C <sub>16</sub> H <sub>33</sub> NO <sub>3</sub>	n,n-DDA	n,n-bis(2-hydroxyethyl) dodecanamide, anti-static agent
289.1417	[A <sub>4</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
293.2457	[AB <sub>1</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
295.2273	[AB <sub>2</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
297.0830	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>4</sub>	Polysiloxane	Polysiloxane, followed by m/z
301.1416	[M+Na] <sup>+</sup>	C <sub>16</sub> H <sub>22</sub> O <sub>4</sub>	Dibutylphthalate	Dibutylphthalate, plasticizer
304.2616	[M+Na] <sup>+</sup>	C <sub>18</sub> H <sub>35</sub> NO	Oleamide	Slip agent in polyethylene films
305.1576	[A <sub>6</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
306.2773	[M+Na] <sup>+</sup>	C <sub>18</sub> H <sub>37</sub> NO	Stearamide	Slip agent in polyethylene films
309.2277	[A <sub>5</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
309.2430	[AB <sub>2</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
315.2535	[M+H] <sup>+</sup>	C <sub>18</sub> H <sub>34</sub> O <sub>4</sub>	DBS	Dibutyl sebacate, plasticizer
317.1155	[M+K] <sup>+</sup>	C <sub>16</sub> H <sub>22</sub> O <sub>4</sub>	Dibutylphthalate	Dibutylphthalate, plasticizer
317.2093	[AB <sub>2</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
321.1316	[A <sub>6</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
323.2562	[AB <sub>2</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
325.2590	[M <sub>2</sub> +H] <sup>+</sup>	C <sub>8</sub> H <sub>18</sub> O <sub>3</sub>	DGBE	Diethylene glycol monobutyl ether, cpd. In scintillation cocktail
327.0786	[M+H] <sup>+</sup>	C <sub>18</sub> H <sub>15</sub> O <sub>4</sub> P	TPP	Triphenyl phosphate, flame retardant in plastics
327.2019	[A <sub>7</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
331.2097	[A <sub>5</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
331.2249	[AB <sub>2</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
337.1190	[M+H] <sup>+</sup> ; ( <sub>120</sub> Sn)	C <sub>13</sub> H <sub>28</sub> O <sub>2</sub> S <sub>n</sub>	Tributyl tin formate	Tributyl tin formate, catalyst
337.2719	[AB <sub>2</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
338.3423	[M+H] <sup>+</sup>	C <sub>22</sub> H <sub>43</sub> NO	Erucamide	Erucamide, (Cis-13-docosenoic amide)
339.2535	[AB <sub>3</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
347.1836	[A <sub>5</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
349.1838	[A <sub>7</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
353.2692	[AB <sub>3</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
355.0705	[M+H-CH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>5</sub>	Polysiloxane	Polysiloxane, (neutral methane loss from m/z 371)
355.3688	[M-Cl] <sup>+</sup>	C <sub>22</sub> H <sub>47</sub> N <sub>2</sub> OCl	PATC	Palmitamidopropyl-trimonium chloride, personal care products additive
360.3242	[M+Na] <sup>+</sup>	C <sub>22</sub> H <sub>43</sub> NO	Erucamide	Erucamide, (Cis-13-docosenoic amide)
361.2355	[AB <sub>3</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
365.1578	[A <sub>7</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
367.2696	[A <sub>6</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
367.2824	[AB <sub>3</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
368.4256	[M-Cl] <sup>+</sup>	C <sub>25</sub> H <sub>54</sub> NCI	BTAC-228	Behentrimonium chloride, personal care product additive
371.1018	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>5</sub>	Polysiloxane	Polysiloxane, followed by m/z 388
371.2281	[A <sub>8</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
371.3161	[M+H] <sup>+</sup>	C <sub>22</sub> H <sub>42</sub> O <sub>4</sub>	DEHA	Bis(2-ethylhexyl) adipate, plasticizer
371.3161	[M+H] <sup>+</sup>	C <sub>22</sub> H <sub>42</sub> O <sub>4</sub>	DOA	Diocetyl adipate, plasticizer
375.2511	[AB <sub>3</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
381.2981	[AB <sub>3</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
383.2797	[AB <sub>4</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
388.1284	[M+NH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>5</sub>	Polysiloxane	Polysiloxane, (see m/z 371)
389.2515	[A <sub>6</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
391.2848	[M+H] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticiser
393.2101	[A <sub>8</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
397.2954	[AB <sub>4</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
405.2255	[A <sub>6</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
405.2617	[AB <sub>4</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
409.1840	[A <sub>8</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
411.3086	[AB <sub>4</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
413.2668	[M+Na] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticiser
415.2543	[A <sub>9</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
419.2773	[AB <sub>4</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
425.3114	[A <sub>7</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
425.3243	[AB <sub>4</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
427.3060	[AB <sub>5</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
429.0893	[M+H-CH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>6</sub>	Polysiloxane	Polysiloxane, (neutral methane loss from m/z 445)
429.2407	[M+K] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticiser
437.2363	[A <sub>9</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether



Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
441.3216	[AB <sub>5</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
445.1206	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>6</sub>	Polysiloxane	Polysiloxane, followed by m/z 462
447.2934	[A <sub>7</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
449.2879	[AB <sub>5</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
449.3856	[M <sub>7</sub> +H] <sup>+</sup>	C <sub>13</sub> H <sub>24</sub> N <sub>2</sub> O	DCU	N,N'-Dicyclohexylurea
453.2102	[A <sub>9</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
453.3441	[M+H] <sup>+</sup>	C <sub>24</sub> H <sub>44</sub> N <sub>4</sub> O <sub>4</sub>	nylon	Cyclic oligomer of polyamide 66, (adipic acid/hexylmethylenediamine condensation)
454.2933	[M+CH <sub>3</sub> CN+Na] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticiser
455.3349	[AB <sub>5</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
459.2805	[A <sub>10</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
462.1471	[M+NH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>6</sub>	Polysiloxane	Polysiloxane (see m/z 445)
463.2673	[A <sub>7</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
463.3036	[AB <sub>5</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
469.3505	[AB <sub>5</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
471.3322	[AB <sub>6</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
481.2625	[A <sub>10</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
483.3533	[A <sub>8</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
485.3478	[AB <sub>6</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
493.3141	[AB <sub>6</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
494.5665	[M-Cl] <sup>+</sup>	C <sub>34</sub> H <sub>72</sub> NCI	DPDMA	Dipalmitylidimethylammonium chloride, catalyst, personal care products additive
497.2364	[A <sub>10</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
499.3611	[AB <sub>6</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
484.3376	[M+H-CH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>7</sub>	Polysiloxane	Polysiloxane, (neutral methane loss from m/z 519)
503.3068	[A <sub>11</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
505.3353	[A <sub>8</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
507.3298	[AB <sub>6</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
513.3767	[AB <sub>6</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
515.3584	[AB <sub>7</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
515.4134	[M+H] <sup>+</sup>	C <sub>30</sub> H <sub>58</sub> O <sub>4</sub> S	DDTDP	Didodecyl 3,3'-thiodipropionate, antioxidant
519.1394	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>7</sub>	Polysiloxane	Polysiloxane, followed by m/z 536
521.3092	[A <sub>8</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
522.5978	[M-Cl] <sup>+</sup>	C <sub>36</sub> H <sub>76</sub> NCI	SPDMA	Stearyl-palmitylidimethylammonium chloride, catalyst, personal care product additive
525.2887	[A <sub>11</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
529.3740	[AB <sub>7</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
531.4083	[M+H] <sup>+</sup>	C <sub>30</sub> H <sub>58</sub> O <sub>5</sub> S	DDTDP	Didodecyl 3,3'-thiodipropionate oxidized to sulfoxide, antioxidant
531.4777	[M+H] <sup>+</sup>	C <sub>35</sub> H <sub>62</sub> O <sub>3</sub>	Irganox	Irganox 1076, antioxidant in synthetic polymers, antioxidant
536.1659	[M+NH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>7</sub>	Polysiloxane	Polysiloxane (see m/z 519)
537.3403	[AB <sub>7</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
537.8796	[M <sub>6</sub> -6H+ <sub>3</sub> Fe+O] <sup>+</sup>	C <sub>2</sub> H <sub>4</sub> O <sub>2</sub>	Acetic acid-Fe-O-complex	during ESI with metal tips and acetic acid
541.2626	[A <sub>11</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
541.3952	[A <sub>9</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
543.3873	[AB <sub>7</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
547.3330	[A <sub>12</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
547.4032	[M+H] <sup>+</sup>	C <sub>30</sub> H <sub>58</sub> O <sub>6</sub> S	DDTDP	Didodecyl 3,3'-thiodipropionate oxidized to sulfone, antioxidant
550.6291	[M-Cl] <sup>+</sup>	C <sub>38</sub> H <sub>80</sub> NCI	DSDMA	Distearylidimethylammonium chloride, catalyst, personal care products additive
551.3560	[AB <sub>7</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
553.3903	[M+Na] <sup>+</sup>	C <sub>30</sub> H <sub>58</sub> O <sub>5</sub> S	DDTDP	Didodecyl 3,3'-thiodipropionate oxidized to sulfoxide, antioxidant
553.4597	[M+Na] <sup>+</sup>	C <sub>35</sub> H <sub>62</sub> O <sub>3</sub>	Irganox	Irganox 1076, antioxidant in synthetic polymers, antioxidant
555.8902	[M <sub>6</sub> -6H+H <sub>2</sub> O+3Fe+O] <sup>+</sup>	C <sub>2</sub> H <sub>4</sub> O <sub>2</sub>	Acetic acid-Fe-O-complex	during ESI with metal tips and acetic acid
557.4029	[AB <sub>7</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
559.3846	[AB <sub>8</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
563.3771	[A <sub>9</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
569.3149	[A <sub>12</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
573.4003	[AB <sub>8</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
577.1269	[M+H-CH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>8</sub>	Polysiloxane	Polysiloxane, (neutral methane loss from m/z 593)
579.3511	[A <sub>9</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
581.3666	[AB <sub>8</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents

Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
585.2888	[A <sub>12</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
587.4135	[AB <sub>8</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
591.3592	[A <sub>13</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
593.1582	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>8</sub>	Polysiloxane	Polysiloxane, followed by m/z 610
595.3822	[AB <sub>8</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
597.9007	[M <sub>7</sub> -6H+3Fe+O] <sup>+</sup>	C <sub>2</sub> H <sub>4</sub> O <sub>2</sub>	Acetic acid-Fe-O-complex	during ESI with metal tips and acetic acid
599.4370	[A <sub>10</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
601.4292	[AB <sub>8</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
603.4108	[AB <sub>9</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
610.1847	[M+NH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>8</sub>	Polysiloxane	Polysiloxane (see m/z 593)
613.3411	[A <sub>13</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
615.4043	[M+H] <sup>+</sup>	C <sub>32</sub> H <sub>58</sub> N <sub>2</sub> O <sub>2</sub> S	CHAPS	3-[(3-Cholamidopropyl)dimethylammonio]-1-propanesulfonate
617.4265	[AB <sub>9</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
621.4190	[A <sub>10</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
621.9735	[M <sub>6</sub> -6H+3Fe+O] <sup>+</sup>	C <sub>3</sub> H <sub>6</sub> O <sub>2</sub>	Propionic acid Fe-O complex	during ESI with metal tips and acetic acid
625.3928	[AB <sub>9</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
629.3151	[A <sub>13</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
631.4397	[AB <sub>9</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
635.3854	[A <sub>14</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
637.3929	[A <sub>10</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
639.4084	[AB <sub>9</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
645.4554	[AB <sub>9</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
647.4370	[AB <sub>10</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
651.1456	[M+H-CH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>9</sub>	Polysiloxane	Polysiloxane, (neutral methane loss from m/z 667)
657.3673	[A <sub>14</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
657.4789	[A <sub>11</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
661.4527	[AB <sub>10</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
667.1769	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>9</sub>	Polysiloxane	Polysiloxane, followed by m/z 684
669.4190	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
673.3413	[A <sub>14</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
675.4659	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
679.4116	[A <sub>15</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
679.4608	[A <sub>11</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
679.5122	[M+H] <sup>+</sup>	C <sub>36</sub> H <sub>66</sub> N <sub>6</sub> O <sub>6</sub>	nylon	Cyclic oligomer of polyamide 66, (adipic acid/hexylmethylenediamine condensation)
683.4346	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
684.2035	[M+NH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>9</sub>	Polysiloxane	Polysiloxane (see m/z 667)
689.4816	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton, reduced	101R Detergents
691.4633	[AB <sub>11</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
695.4348	[A <sub>11</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
701.3936	[A <sub>15</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
705.4789	[AB <sub>11</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
713.4452	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
715.5208	[A <sub>12</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
717.3675	[A <sub>15</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
719.4921	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
723.4378	[A <sub>16</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
725.1644	[M+H-CH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>10</sub>	Polysiloxane	Polysiloxane, (neutral methane loss from m/z 741)
727.4608	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
733.5078	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
735.4895	[AB <sub>12</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
737.5027	[A <sub>12</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
741.1957	[M+H] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>10</sub>	Polysiloxane	Polysiloxane, followed by m/z 758
745.4198	[A <sub>16</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
749.5051	[AB <sub>12</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
753.4766	[A <sub>12</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
757.4714	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
758.2223	[M+NH <sub>4</sub> ] <sup>+</sup>	[C <sub>2</sub> H <sub>6</sub> SiO] <sub>10</sub>	Polysiloxane	Polysiloxane (see m/z 741)
761.3937	[A <sub>16</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
763.5184	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
767.4640	[A <sub>17</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
771.4871	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
773.5626	[A <sub>13</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether

Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
777.5340	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
779.5157	[AB <sub>13</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
789.4460	[A <sub>17</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
793.5313	[AB <sub>13</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
795.5446	[A <sub>13</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
798.5884	[M <sub>2</sub> +NH <sub>4</sub> ] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticiser
801.4976	[AB <sub>13</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
803.5438	[M <sub>2</sub> +Na] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticiser
805.4199	[A <sub>17</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
807.5446	[AB <sub>13</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
809.4875	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>18</sub> H <sub>34</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 20
811.4903	[A <sub>18</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
811.5185	[A <sub>13</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
815.5133	[AB <sub>13</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
819.5177	[M <sub>2</sub> +K] <sup>+</sup>	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>	Diisooctyl phthalate	Diisooctyl phthalate, plasticiser
821.5602	[AB <sub>13</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
823.5419	[AB <sub>14</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
831.6045	[A <sub>14</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
833.4722	[A <sub>18</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
837.5575	[AB <sub>14</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
845.5238	[AB <sub>14</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
849.4461	[A <sub>18</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
851.5708	[AB <sub>14</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
853.5137	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>18</sub> H <sub>34</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 20
853.5864	[A <sub>14</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
855.5165	[A <sub>19</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
859.5395	[AB <sub>14</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
865.5501	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>22</sub> H <sub>42</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 40
865.5864	[AB <sub>14</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
867.5681	[AB <sub>15</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
869.5604	[A <sub>14</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
877.4984	[A <sub>19</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
881.5838	[AB <sub>15</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
889.5501	[AB <sub>15</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
889.6464	[A <sub>15</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
891.5657	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>24</sub> H <sub>44</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 80
893.4724	[A <sub>19</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
893.5814	[AB <sub>10</sub> +Na] <sup>+</sup>	[C <sub>24</sub> H <sub>46</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 60
897.5399	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>18</sub> H <sub>34</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 20
899.5427	[A <sub>20</sub> B+H] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
903.5657	[AB <sub>15</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
905.6803	[M+H] <sup>+</sup>	C <sub>48</sub> H <sub>88</sub> N <sub>8</sub> O <sub>8</sub>	nylon	Cyclic oligomer of polyamide 66, (adipic acid/hexamethylene diamine condensation)
909.5763	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>22</sub> H <sub>42</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 40
909.6127	[AB <sub>15</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
911.5943	[AB <sub>16</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
911.6283	[A <sub>15</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
921.5246	[A <sub>20</sub> B+Na] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
925.6100	[AB <sub>16</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
927.6022	[A <sub>15</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
933.5763	[AB <sub>16</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
935.5919	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>24</sub> H <sub>44</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 80
937.4986	[A <sub>20</sub> B+K] <sup>+</sup>	[C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub> H <sub>2</sub> O	PEG	Polyethylene glycol, ubiquitous polyether
937.6076	[AB <sub>11</sub> +Na] <sup>+</sup>	[C <sub>24</sub> H <sub>46</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 60
939.6232	[AB <sub>16</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
941.5661	[AB <sub>13</sub> +Na] <sup>+</sup>	[C <sub>18</sub> H <sub>34</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 20
947.5919	[AB <sub>16</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
947.6882	[A <sub>16</sub> B+H] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
953.6025	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>22</sub> H <sub>42</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 40
953.6389	[AB <sub>16</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
955.6205	[AB <sub>17</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
969.6362	[AB <sub>17</sub> +H] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
969.6702	[A <sub>16</sub> B+Na] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
977.6025	[AB <sub>17</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents
979.6181	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>24</sub> H <sub>44</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 80

Mono-isotopic ion mass (singly charged)	Ion type	Formula for M or subunit or sequence	Compound ID or species	Possible origin and other comments
981.6338	[AB <sub>12</sub> +Na] <sup>+</sup>	[C <sub>24</sub> H <sub>46</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 60
983.6494	[AB <sub>17</sub> +Na] <sup>+</sup>	[C <sub>14</sub> H <sub>28</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	X-100R, X-114R, X-405R, or X-45R Detergents
985.5923	[AB <sub>14</sub> +Na] <sup>+</sup>	[C <sub>18</sub> H <sub>34</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 20
985.6441	[A <sub>16</sub> B+K] <sup>+</sup>	[C <sub>3</sub> H <sub>6</sub> O] <sub>n</sub> H <sub>2</sub> O	PPG	Polypropylene glycol, ubiquitous polyether
991.6181	[AB <sub>17</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>24</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	101 Detergents
997.6287	[AB <sub>13</sub> +Na] <sup>+</sup>	[C <sub>22</sub> H <sub>42</sub> O <sub>6</sub> ][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Tween®	Tween® 40
997.6651	[AB <sub>17</sub> +Na] <sup>+</sup>	[C <sub>15</sub> H <sub>30</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®, reduced	101R Detergents
999.6470	[AB <sub>18</sub> +H] <sup>+</sup>	[C <sub>14</sub> H <sub>22</sub> O][C <sub>2</sub> H <sub>4</sub> O] <sub>n</sub>	Triton®	X-100, X-114, X-405, or X-45 Detergents

# Appendix II.

## Monoisotopic ion masses of commonly observed repeating units in LC-MS

### Positive ion

Mass difference	Origin
14.0157	-[CH <sub>2</sub> ]-, alkane chains, waxes, fatty acids, methylation
15.9949	O, oxidation
18.0106	H <sub>2</sub> O, water clusters
28.0313	-[C <sub>2</sub> H <sub>4</sub> ]-, natural alkane chains such as fatty acids
32.0262	CH <sub>3</sub> OH, methanol clusters
41.0266	CH <sub>3</sub> CN, acetonitrile clusters
42.0470	-[C <sub>3</sub> H <sub>6</sub> ]-, propyl repeating units, propylation
44.0262	-[C <sub>2</sub> H <sub>4</sub> O]-; polyethylene glycol, PEG, and related components such as Triton®- and Tween®-containing buffers
49.9968	-[CF <sub>2</sub> ]-, from perfluoro compounds
53.0032	NH <sub>4</sub> Cl salt adducts/clusters
56.0626	-[C <sub>4</sub> H <sub>8</sub> ]-, butyl repeating units, butylation
57.9586	NaCl, sodium chloride clusters
58.0419	-[C <sub>3</sub> H <sub>6</sub> O]-; polypropylene glycol and related compounds, PPG, and related compounds
63.0320	CHOONH <sub>4</sub> , ammonium formate adducts/clusters
67.9874	NaHCO <sub>2</sub> , sodium formate clusters
67.9874	CHOONa, sodium formate adducts/clusters
72.0395	-OH replacement with -OSi(CH <sub>3</sub> ) <sub>3</sub> ,(=[C <sub>3</sub> H <sub>8</sub> Si]), trimethylsiloxane, endcapping reagent
73.9326	KCl adducts/clusters
74.0188	-[O-Si(CH <sub>3</sub> ) <sub>2</sub> ]-, polysiloxane, silicone rubber polymer (typical series at m/z's 355, 429, 503, 593, 667, 741, 815)
78.0139	C <sub>2</sub> H <sub>6</sub> OS, DMSO adducts/clusters, dimethylsulfoxide solvent
82.0031	NaCH <sub>3</sub> CO <sub>2</sub> , sodium acetate clusters
84.0516	C <sub>2</sub> D <sub>6</sub> OS, deuterated DMSO adducts/clusters, NMR solvent
135.9748	NaCF <sub>3</sub> CO <sub>2</sub> , sodium trifluoroacetate clusters
162.0528	-[C <sub>6</sub> H10O <sub>5</sub> ]-, polysaccharides residues
226.1681	-[C1 <sub>2</sub> H <sub>22</sub> N <sub>2</sub> O <sub>2</sub> ]-, cyclic oligomers from polyamide 66 (series observed with m/z 453, 679, 905)
259.8099	CsI, cesium iodide clusters, used as calibration

## To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645  
Germany: 069 86798021  
Italy: 848 845 645  
Spain: 901 516 645 Option 1  
Switzerland: 0848 645 645  
United Kingdom: 0870 900 4645

For other countries across Europe,  
please call: +44 (0) 115 943 0840

For other countries across Europe  
and the world, please visit:

**[MerckMillipore.com/offices](http://MerckMillipore.com/offices)**

For Technical Service visit:

**[MerckMillipore.com/techservice](http://MerckMillipore.com/techservice)**  
**[MerckGroup.com/life-science](http://MerckGroup.com/life-science)**

Merck KGaA  
Frankfurter Strasse 250  
64293 Darmstadt  
Germany

