3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**



## **Protein Array Procedure**

### **Product Description**

The binding specificity of the purified Prestige Antibodies® is determined on protein arrays, using PrEST protein fragments, to ensure high specificity toward its antigen and low background binding.

#### **Preparation Instructions**

PBS -8.1 mM Na<sub>2</sub>HPO<sub>4</sub>, 1.5 mM KH<sub>2</sub>PO<sub>4</sub>, 137 mM NaCl, and 2.7 mM KCl, pH 7.4

PBST - PBS with 0.1% (v/v) TWEEN® 20

Blocking Buffer – PBST with 3% bovine serum albumin

### **Procedure**

All washes are performed at room temperature (RT) on a shaker.

Antibody dilution with PBST:

- for antibody concentrations ≤0.04 mg/ml, 1:500
- for antibody concentrations ≥0.05 mg/ml, 1:3,000 <u>Note</u>: The specified working dilutions of the antibodies are to be considered as guidelines only. Optimal dilutions must be determined by the user.
- The PrEST protein (the antigen) is diluted 1:10 in PBS containing BSA (100 μg/ml) and printed with 383 other PrEST antigens on epoxy coated microarray slides.
- 2. The slides are incubated overnight at 37 °C.
- Unreacted epoxide groups are blocked by incubating the printed slide in Blocking Buffer for 60 minutes.

- 4.  $60~\mu l$  of the primary antibody solution are applied to each well and the slide is incubated for 60 minutes. The arrays are separated by a silicon mask.
- The slide is washed twice, each time for 5 minutes with PBST followed by a third wash with PBS for 5 minutes.
- 6. The secondary antibody solution is prepared by diluting the goat anti-rabbit IgG Alexa Fluor® 647 conjugate 60,000-fold in PBST.
  Note: The secondary antibody is fluorescently labeled, and thus, light sensitive. The slide should be kept in the dark during the rest of the procedure.
- 7. The secondary antibody solution is added to the slide, which is incubated for 60 minutes in the dark.
- 8. The slide is washed twice, each time for 5 minutes with PBST followed by a third wash with PBS for 5 minutes in the dark.
- 9. The slide is spin dried for 10–15 seconds using a table centrifuge prior to scanning.

Prestige Antibodies is a registered trademark of Sigma-Aldrich® Biotechnology LP and Sigma-Aldrich Co. Alexa Fluor is a registered trademark of Molecular Probes.

TWEEN is a registered trademark of Croda International PLC.

AKN, KB, LPG, MAM 02/11-1