

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Hexadecyltrimethylammonium bromide

Catalog Number **H5882**Store at Room Temperature

CAS Number: 57-09-0 Synonym: CTAB

Product Description

Molecular Formula: C₁₉H₄₂BrN Molecular Weight: 364.5 Melting Point: 237-243 °C¹

Critical Micelle Concentration (CMC): 0.92 to 1.0 mM

(water)2,3

Aggregation number: 61 (water, 25 °C),

169 (13 mM KBr)^{2,4}

pH: 6.0-7.5 (0.1 M H₂O, 0 °C).⁵

Hexadecyltrimethylammonium bromide is a bactericidal, cationic detergent. Its activity is neutralized by soaps and anionic detergents (i.e., SDS). Trimethylammonium bromide compounds form insoluble complexes with SDS.³ It is active at alkaline pH against both Grampositive and Gram-negative organisms.⁶

CTAB has been used for the isolation of plant high molecular weight DNA (by a rapid method)⁷ as well as plant DNA for use in PCR analysis.^{8,9} It has been used to precipitate nucleic acids.^{7,8,10,11}

CTAB provides enhancement of Concanavalin A mediated agglutination. ¹² It has been used for the determination of protein molecular weights in electrophoretic systems ^{13,14} and for determination of critical micelle concentration of detergents. ¹⁵ CTAB has also been used as a titrant for potentiometric titration of perchlorate ¹⁶ and as a phase-transfer catalyst in reduction of arenes and heterocyclic compounds. ¹⁷

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

The product is soluble in water (100 mg/ml). It is freely soluble in alcohol; sparingly soluble in acetone. It is practically insoluble in ether and benzene.¹

Storage/Stability

The product is stable in acid solution.¹

References

- 1. The Merck Index, 11th ed., Entry# 2018.
- 2. Neugebauer, J. M., Detergents: An Overview. Meth. Enzymol., **182**, 239-253 (1990).
- Protein Purification Applications: A Practical Approach, Harris, E. L. V., and Angal, S., eds., IRL Press at Oxford University Press (New York, NY: 1990), p. 71.
- 4. Helenius, A., and Simons, K., Solubilization of Membranes by Detergents. Biochim. Biophys. Acta., **415(1)**, 29-79 (1975).
- Biochemika 1989 Handbook, p. 557.
- Data for Biochemical Research, 3rd ed., Dawson,
 R. M. C., et al., Oxford University Press (New York,
 NY: 1986), p. 287.
- Murray, M. G., and Thompson, W. F., Rapid Isolation of High Molecular Weight Plant DNA. Nucleic Acids Res., 8(19), 4321-4325 (1980).
- Short Protocols in Molecular Biology, 3rd ed., Ausubel, F. M., et al., Wiley (New York, NY: 1995), pp. 2-10.
- 9. Dempster, E. L., et al., Rapid DNA Extraction from Ferns for PCR-based Analyses. Biotechniques, **27(1)**, 66-68 (1999).
- Gustincich, S., et al., A Fast Method for Highquality Genomic DNA Extraction from Whole Human Blood. Biotechniques, 11(3), 298-302 (1991).
- Morimoto H., et al., Spectrophotometric Analysis of RNA and DNA Using Cetyltrimethylammonium bromide. Anal. Biochem., 62(2), 436-448 (1974).

- 12. Isomaa, B., and Hagerstrand, H., ConA-mediated Agglutination in Rat T-lymphocytes is Enhanced by the Membrane-perturbing Agent CTAB. Cell Mol. Biol., **32(5)**, 627-630 (1986).
- Panyim, S., et al., A Simplified Gel Electrophoretic System and Its Validity for Molecular Weight Determinations of Protein-cetyltrimethylammonium Complexes. Anal. Biochem., 81(2), 320-327 (1977).
- Akins, R. E., et al., Cetyltrimethylammonium Bromide Discontinuous Gel Electrophoresis: Mrbased Separation of Proteins with Retention of Enzymatic Activity. Anal. Biochem., 202(1), 172-178 (1992).
- Samsonoff, C., et al., The Use of Coomassie Brilliant Blue for Critical Micelle Concentration Determinatio of Detergents. J. Colloid Interface Sci., 109(2), 325-329 (1986).
- 16. Talanta, 26, 1061 (1979).
- Januszkiewicz, K. R., and Alper, H., Exceedingly Mild, Selective and Stereospecific Phase-Transfer-Catalyzed Hydrogenation of Arenes. Organometallics, 2, 1055-1057 (1983).

FEB,RXR,MAM 01/16-1