SigmaAldrich.com

Sigma-Aldrich.

Product Information

Peroxidase from horseradish

Sigma Type VI-A, essentially salt-free, lyophilized powder, 950-2000 units/mg solid (using ABTS), \geq 250 units/mg solid (using pyrogallol)

P6782

Product Description

EC Number: 1.11.1.7 CAS Registry Number: 9003-99-0 Synonym: Hydrogen peroxide oxidoreductase; HRP

Horseradish peroxidase (HRP) is isolated from the roots of horseradish (*Amoracia rusticana*) and belongs to the ferroprotoporphyrin group of peroxidases. HRP readily combines with hydrogen peroxide (H_2O_2) . The resultant [HRP-H₂O₂] complex can oxidize a wide variety of hydrogen donors:

Donor + $H_2O_2 \rightarrow Oxidized Donor + 2 H_2O$

HRP will oxidize various substrates (see Table 1):

- Chromogenic
- Chemiluminescent (such as luminol or isoluminol)
- Fluorogenic (such as tyramine, homovanillic acid, or 4-hydroxyphenyl acetic acid)

HRP is a single chain polypeptide that contains four disulfide bridges. HRP is a glycoprotein that contains 18% carbohydrate. The carbohydrate composition consists of galactose, arabinose, xylose, fucose, mannose, mannosamine, and galactosamine, depending upon the specific isozyme.¹

HRP is a widely used label for immunoglobulins in many different immunochemistry applications, including immunoblotting, immunohistochemistry, and ELISA. HRP can be conjugated to antibodies by several different methods, including glutaraldehyde, periodate oxidation, through disulfide bonds, and also via amino and thiol directed cross-linkers. HRP is the most desired label for antibodies, since it is the smallest and most stable of the three most popular enzyme labels (peroxidase, β -galactosidase, alkaline phosphatase) and its glycosylation leads to lower non-specific binding.² A review of glutaraldehyde and periodate conjugation methods has been published.³

Peroxidase is also utilized for the determination of glucose⁴ and peroxides⁵ in solution. Several publications⁶⁻¹⁸ and dissertations¹⁹⁻⁴⁰ have cited use of P6782 in their research protocols.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Reagent

This product is supplied as an essentially salt-free, lyophilized powder. Please note that Cat. No. P6782 is sold by weight (as MG, or milligrams), rather by enzyme units (such as KU, or kilounits).

Specific Activity:

- \geq 250 units/mg solid (pyrogallol as substrate)
- 950-2,000 units/mg solid (ABTS as substrate)

Note: Using 2,2'-Azino-bis(3-ethylbenzthiazoline-6-sulfonic acid) (ABTS) tablets (Cat. No. A9941) as the substrate, ~4× the activity is observed when compared to pyrogallol (purpurogallin unit).

Unit definition (purpurogallin): One unit will form 1.0 mg of purpurogallin from pyrogallol in 20 seconds at pH 6.0 at 20 °C. This unit is equivalent to ~18 μ M units per minute at 25 °C.

RZ (Reinheitszahl): 2.5 – 4.0

RZ is the absorbance ratio A_{403}/A_{275} determined at 0.5-1.0 mg/mL in deionized water. RZ is a measure of hemin content, **not** enzymatic activity. Even preparations with high RZ values may have low enzymatic activity.

Total molecular mass:⁴¹ ~44 kDa (~44,000 Da)

- Polypeptide chain: 33,890 Da
- Hemin plus Ca²⁺: ~700 Da
- Carbohydrate: 9,400 Da

Extinction coefficient:⁴² $E^{mM} = 100 (403 \text{ nm})$ Optimal pH range:⁴³ 6.0-6.5 (activity at pH 7.5 is 84% of the maximum)



The enzyme is most stable in the pH range of 5.0-9.0. Isoelectric point:¹ isozymes range from 3.0-9.0 (at least seven isozymes)

Inhibitors:⁴⁴ sodium azide; cyanide; L-cystine; dichromate; ethylenethiourea; hydroxylamine; sulfide; vanadate; *p*-aminobenzoic acid; Cd²⁺, Co²⁺, Cu²⁺, Fe³⁺, Mn²⁺, Ni²⁺, Pb²⁺ ions

Preparation Instructions

This product is soluble in water or 0.1 M phosphate buffer, pH 6.0.

Storage/Stability

Store the product at 2-8 °C. The enzyme remains active for at least five years. The product may also be stored at -20 °C, if desired. The lyophilized powder will retain at least 80% of its activity after two weeks at 37 °C.

A 1 mg/mL solution in 0.1 M phosphate buffer, pH 6.0, remains active for at least two weeks at room temperature. A solution retains activity after 5 freeze-thaw cycles.

References

- Shannon, L.M. et al., J. Biol. Chem., 241(9), 2166-2172 (1966).
- 2. Deshpande, S.S., *Enzyme Immunoassays, From Concept to Product Development*. Chapman and Hall (New York, NY), pp. 169-171 (1996).
- Harlow, E. and Lane, D., *Antibodies: A Laboratory Manual*. Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY), pp. 346-348 (1988).
- Bergmeyer, H.-U. *et al.*, *Methods of Enzymatic Analysis*, 2nd ed. (Bergmeyer, H.-U., ed.). Verlag Chemie/Academic Press, pp. 1205-1227 (1974).
- Bernt, E. and Bergmeyer, H.-U., *Methods of Enzymatic Analysis*, 2nd ed., pp. 2246-2248 (1974).
- Andrews, J. et al., Plant Cell Environment, 23(4), 415-422 (2000).
- Nakajima, T. et al., Cell Transplant., 10(4-5), 423-427 (2001).
- Zhou, J. et al., Anal. Chem., 74(16), 4007-4010 (2002).
- Wagschal, K. et al., Appl. Environ. Microbiol., 71(9), 5318-5323 (2005).
- Zucca, P. et al., Chem. Cent. J., 6(1), 161 (2012).
- 11. Krainer, F.W. et al., Sci. Rep., 3, 3279 (2013).
- 12. Rhee, H.W. *et al.*, *Science*, **339(6125)**, 1328-1331 (2013).

- 13. Krainer, F.W. *et al.*, *J. Biotechnol.*, **233**, 181-189 (2016).
- 14. Pekarsky, A. *et al.*, *Microb. Cell Fact.*, **17(1)**, 183 (2018).
- Regidi, S. *et al.*, *BMC Res. Notes*, **11(1)**, 596 (2018).
- Bredow, M. *et al.*, *J. Vis. Exp.*, (147), doi: 10.3791/59437 (2019).
- 17. Humer, D. *et al.*, *Biomed. Pharmacother.*, **142**, 112037 (2021).
- Ivanov, Y.D. et al., Sci. Rep., **11(1)**, 9907 (2021).
- Connolly, Christopher Neil Peter, "Development of an enzymatic probe for the exocytic pathways of animal cells". University of London, Ph.D. dissertation, pp. 123, 140 (1993).
- Geyer, Cyril, "The coupling of biotechnology with electrochemistry: Electrodenitrification of water using the enzyme nitrate reductase". University of Northern Iowa, Ph.D. dissertation, p. 69 (2001).
- McFetridge, Peter Stuart, "The use of porcine carotid arteries as a matrix material for tissue engineered small diameter vascular grafts". University of Bath, Ph.D. dissertation, p. 54 (2002).
- 22. Sesay, Adama Marie, "Towards a remote portable bio-affinity surface plasmon resonance analyser for environmental steroidal pollutants". Cranfield University, Ph.D. dissertation, p. 138 (2003).
- Powers, Daryl E., "Effects of Oxygen on Embryonic Stem Proliferation, Energetics, and Differentiation into Cardiomyocytes". Massachusetts Institute of Technology, Ph.D. dissertation, p. 22 (2003).
- Németh, Edina Csibrikné, "Investigation of the efficiency and mode of action of lactic acid bacteria suitable for probiotic food production". Central Food Research Institute (Budapest), Ph.D. dissertation, p. 35 (2007).
- Burkhardt, Nicole, "The role of glutamate synthesis and release mechanisms in hormone secretion of pancreatic islet cells". Université de Genève, Ph.D. dissertation, p. 106 (2009).
- Bechoff, Aurélie, "Investigating carotenoid loss after drying and storage of orange-fleshed sweet potato". University of Greenwich, Ph.D. dissertation, p. 154 (April 2010).
- Heijnis, Walter Hendrik, "Peroxidase-mediated cross-linking of bovine a-lactalbumin". Wageningen University, Ph.D. dissertation, pp. 28, 42, 58, 70 (2010).



- Griffin, Carly, "Evaluation of the Selective Anti-Cancer Activity of Natural and Synthetic Alkaloids". University of Windsor, Ph.D. dissertation, p. 179 (2011).
- 29. Paper, Janet Marie, "Secreted Proteins from the Plant Pathogenic Fungus *Fusarium Graminearum*: Identification, Role in Disease, and Biotechnological Applications". Michigan State University, Ph.D. dissertation, p. 126 (2011).
- Williamson, David, "The role of LAT clusters in early immune cell signalling". University of New South Wales, Ph.D. dissertation, p. 75 (2011).
- Ungureanu, Felicia, "Sensing with colors". University of Twente, Ph.D. dissertation, p. 73 (2012).
- Guo, Chuan, "Biosensitive Functionalised Silicon Surfaces: towards Biosensitive Field-Effect-Transistors". University of Sydney, Ph.D. dissertation, p. 136 (February 2013).
- Lochbühler, Bernd Christoph, "The potential of yeast proteins to substitute for traditional fining agents – technological and sensory aspects". Justus-Liebig-Universität Gießen, Dr. agr. dissertation, p. 102 (2014).
- Krainer, Florian Wolfgang, "Production and characterization of recombinant plant peroxidases in *Pichia pastoris*". Graz University of Technology, Dr. rer. nat. dissertation, p. 46 (2015).
- Peters, Bridget A., "Experimental woodsmoke exposure during exercise and an evaluation of blood oxidative stress". Auburn University, Ph.D. dissertation, p. 105 (2015).
- Nanadikar, Maithily, "Application of redox biosensor mouse models to study redox processes in cardiomyocytes". Georg-August-Universität Göttingen, Dr. rer. nat. dissertation, p. 35 (2019).
- Nowak, Claudia, "Characterization of Biomimetic Nicotinamide Cofactors and Optimization of Enzymatic Regeneration Processes". Technischen Universität München, Dr. rer. nat. dissertation, p. 36 (2019).
- Abernathy, Breann Elise, "Prebiotic Dietary Fibers Must Achieve a Threshold of Beneficial Gut Bacteria in Order to Prevent Adiposity and Fatty Liver in Rats Fed High Fat Diets". University of Minnesota, Ph.D. dissertation, p. 234 (2020).

- Chen, Xi, "Carbon Nanotube Porin Based Biosensors". University of California Merced, Ph.D. dissertation, p. 65 (2020).
- Fischer, Hillary Donna, "The Influence of Singlet Oxygen and Loss of Function of Fatty Acid Desaturase 7 in the Chloroplast on Aphid Resistance in *Arabidopsis thaliana*". University of Arkansas, Ph.D. dissertation, p. 84 (2021).
- 41. Welinder, K.G., *Eur. J. Biochem.*, **96(3)**, 483-502 (1979).
- 42. Delincée, H. and Radola, B.J., *Eur. J. Biochem.*, **52(2)**, 321-330 (1975).
- 43. Schomberg, D., Salzmann, M., and Stephan, D., Enzyme Handbook 7, EC 1.11.1.7, 1-6 (1993).
- Zollner, H., Handbook of Enzyme Inhibitors, 2nd ed., Part A. VCH Verlagsgesellschaft, pp. 367-368 (1993).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at <u>SigmaAldrich.com/techservice</u>.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at <u>SigmaAldrich.com/terms</u>.

Contact Information

For the location of the office nearest you, go to <u>SigmaAldrich.com/offices</u>.



Table 1. Peroxidase Substrates

Substrate	Cat. No. or Cat. Nos.	Color Reaction	End Product	Applications
2,2'-Azino-bis(3- Ethylbenzthiazoline-6- Sulfonic Acid; ABTS)	A3219, A9941	Green	Soluble	ELISA
o-Phenylenediamine (OPD)	P9187	Orange	Soluble	ELISA
3,3',5,5'- Tetramethylbenzidine (TMB)	T8665, T3405	Blue	Soluble	ELISA
	T0565	Deep Blue	Insoluble	Blotting
o-Dianisdine	D9154	Yellow-Orange	Soluble	ELISA
5-Aminosalicylic Acid (5AS)	A79809, A3537	Brown	Soluble	ELISA
3,3'-Diaminobenzidine (DAB)	D7304, D5905, D4168, D4293, D4418, D7679	Brown	Insoluble	Blotting, Histochemistry
	D0426	Blue-Black		
4-Chloro-1-Naphthol (4C1N)	C6788	Blue	Insoluble	Blotting
3-Amino-9-Ethylcarbazole (AEC)	AEC101, A6926	Red	Insoluble	Blotting
CPS-1	CPS160, CPS1A120, CPS1A300	Chemiluminescent	Soluble	Blotting
CPS-3	CPS350, CPS3100, CPS3500			
CPS-2	CPS260	Chemiluminescent	Soluble	ELISA

Table 2. Other Grades of HRP Available

Cat. No.	RZ value	Specific Activity (*)	
P8250	≥ 1.8	150 – 250 units/mg solid	
P2088	2.6 - 3.4	200-300 units/mg solid	
P8415	≥ 3.0	≥ 250 units/mg solid	
P8125	≥ 1.0	50-150 units/mg solid	
P8375	2.5 - 4.0	≥ 250 units/mg solid	
P6140	2.5 - 3.5	≥ 225 units/mg protein	

(*) Specific activity is reported in terms of purpurogallin units.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. P6782dat Rev 02/22 RBG,CS,MAM,GCY

