

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

# ProductInformation

#### Potassium phosphate monobasic Plant Cell Culture Tested

Product Number **P 8416** Store at Room Temperature

## **Product Description**

Molecular Formula: KH<sub>2</sub>PO<sub>4</sub> Molecular Weight: 136.1 CAS Number: 7778-77-0 Synonyms: monopotassium phosphate, potassium dihydrogen phosphate

This product is plant cell culture tested and is suitable for use in plant cell culture applications.

Potassium phosphate is a reagent with very high buffering capacity that is widely used in molecular biology, biochemistry, and chromatography. Potassium phosphate occurs in several forms: monobasic (KH<sub>2</sub>PO<sub>4</sub>), dibasic (K<sub>2</sub>HPO<sub>4</sub>), and tribasic (K<sub>3</sub>PO<sub>4</sub>). Most neutral potassium phosphate buffer solutions consist of mixtures of the monobasic and dibasic forms to varying degrees, depending on the desired pH. A table for preparation of 0.1 M potassium phosphate buffer at 25 °C using various proportions of potassium phosphate monobasic and potassium phosphate dibasic has been published.<sup>1,2</sup>

Some limitations of the usefulness of phosphate buffers include their precipitation of Ca<sup>2+</sup> and Mg<sup>2+</sup>, their inhibition of restriction enzyme activity, and their interference in protocols related to DNA ligation and bacterial transformation.<sup>1</sup> A study of the effect of freeze-thaw storage cycles on proteins in potassium phosphate and sodium phosphate buffer solutions has been reported.<sup>3</sup> The use of high concentrations of potassium phosphate in the immobilization of affinity ligands onto epoxide-activated stationary phases has been reviewed.<sup>4</sup> A two-phase system of aqueous potassium phosphate and poly(ethylene glycol) for the isolation of *E. coli*  $\beta$ -galactosidase and  $\beta$ -galactosidase fusion proteins has been published.<sup>5</sup> The quantitation of nonionic surfactants in buffered solutions using strong cation and anion exchange HPLC guard columns and potassium phosphate solution has been investigated.<sup>6</sup>

## **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

# **Preparation Instructions**

This product is soluble in water (100 mg/ml), yielding a clear, colorless solution.

#### References

- Molecular Cloning: A Laboratory Manual, 3rd ed., Sambrook, J. F., et al., Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY: 2001), p. A1.5.
- Green, A. A., and Hughes, W. L., Protein Fractionation on the Basis of Solubility in Aqueous Solutions of Salts and Organic Solvents. Meth. Enzymol., 1, 67-90 (1955).
- Pikal-Cleland, K. A., et al., Protein denaturation during freezing and thawing in phosphate buffer systems: monomeric and tetrameric betagalactosidase. Arch. Biochem. Biophys., 384(2), 398-406 (2000).

- Wheatley, J. B., and Schmidt, D. E. Jr., Saltinduced immobilization of affinity ligands onto epoxide-activated supports. J. Chromatogr. A., 849(1), 1-12 (1999).
- Enfors, S. O., et al., Combined use of extraction and genetic engineering for protein purification: recovery of beta-galactosidase fused proteins. Bioseparation, 1(3-4), 305-310 (1990).
- Pardue, K., and Williams, D., Quantitative determination of non-ionic surfactants in protein samples, using ion-exchange guard columns. Biotechniques, **14(4)**, 580-583 (1993).

GCY/RXR 2/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.