

High Sensitivity Binary Gas Integrity Test for Viresolve[®] Pro Devices The Innovative Parvovirus Clearance Solution

ABSTRACT

A novel non-destructive integrity test for microporous and ultraporous membranes has been developed that is far more sensitive to detecting defects than a conventional gas-liquid diffusion test. The test uses a binary gas mixture and is based on the principle of differing gas permeabilities through the liquid layer that results in a concentration enhancement of the faster permeating gas. In an integral membrane, the permeate composition can be predicted based on the transport properties of the gases permeating through the liquid layer and the known operating conditions. A deviation from the expected concentration is an indication of the presence of a defect or open pores. Unlike the gas-liquid diffusion test, the binary gas test has low sensitivity to membrane porosity, liquid layer thickness, and membrane area. Consequently, integral devices will exhibit a relatively narrow range of integral values resulting in a superior defect signal-to-noise ratio. The binary gas integrity test method was applied to the Viresolve Pro device and was demonstrated to provide a significantly higher level of virus retention assurance compared to the air-water diffusion test. Binary gas testing has been implemented as a 100% quality assurance test for Viresolve[®] Pro devices.

BACKGROUND

Unique Technology



Novel Dev

Formats

Scalable

Asymmetric PES		
Membrane		
0	Patented	
	asymmetric	
	membrane	
0	Robust viral	
	clearance	
0	High capacity	
	and flux	
0	Caustic stable	

1	ALL REAL PROPERTY AND INC.
lovel Device	Integrity and
ormats	Performance
Scalable devices	• Testing
design	 LRV/capacity assays
Disposable flow	using model viruses
path	and proteins
Designed to enable	 Liquid-liquid
integrity testing	porometry tests
Manufacturing	 Patented Binary
process design	Gas Test

Building Assurance Test by Test

No one does more than Millipore to help assure virus filter performance and compliance

• End-User Validation and Testing • Device 100% Tests • Device Release Tests • Membrane Release Tests • Viresolve Pro Process & Product Validation

Integrity Test Objective

Purpose of Integrity test is to detect the presence of oversized pores or defects that can compromise the retention capability of a filter

- Must correlate to bacteria/virus retention or retention shortfall
- Must be non-destructive (except post-use tests)
- Ideally provides Log Reduction Value (LRV) assurance consistent with product claim

THEORY

Gas

EXPERIMENTAL VERIFICATION

Conventional Gas-Liquid Diffusion Test



• Air/water test used commonly to test virus filters

- Sufficient and appropriate as an end-user test
- Test sensitivity limitation: convective flux vs. background diffusive flux noise

Experimental Verification of Defect Sensitivity

• 142-mm virus membrane discs (127 cm^2) - Laser drill holes of 2,5,10 micron sizes • 0.7 m2 microporous membrane device - Insert orifices of 0, 2, 5, 10 micron sizes • Perform integrity tests - Air/water diffusion test $-90\% CO_2/C_2F_6$ binary gas test • For virus filter, retention test to 75% fouling with lgG and Θ X174



Binary Gas Test Offers Higher LRV Assurance



[•] Data obtained from double layer membrane discs • Discs drilled with laser holes • Binary Gas Test far more sensitive to defects than air/water diffusion test

Viresolve Pro Solution Includes a Novel Binary Gas Test



Binary Gas Test Defect Sensitivity -Double Layer Membrane Discs



Double layer 142-mm discs of virus filter membrane

Binary Gas Test Offers Higher LRV Assurance



• Data obtained from devices manufactured with various size orifices • Binary Gas Test detects defects as small as 2-3 microns • These defects cannot be detected by current technology (air/water diffusion test)

Binary Gas Test Offers Higher LRV Assurance



Binary Gas Pair Selection



• Selected gas pair and liquid: 90/10 CO_2/C_2F_4 in water • High sensitivity to defects • Easy and accurate composition measurement • High diffusive flow rate of CO₂ through water enables integrity

 $\Pr = \frac{P_p}{D}$

Air Diffusion

1.0

Air Diffusion Permeability (sccm/m²-psi)

¥-1

Integrity Testing of Magnus Devices

0







Effect of Pressure, Temperature, and Purge Ratio on Permeate Composition



🗖 Data

— Model

Effect of Staging on Membrane Gas Separation

50 PSIG

< 100 ppm C₂F₆

0 PSIG

Feed composition = $90/10 CO_2/C_2F_6$

• Concentration depends only on pressure ratio (not additive)

• Flow rate depends only on pressure difference (additive)

Feed conc. = 10.1% C₂F₆/CO₂ Sweep/Feed ratio = 4.7

Feed pressure = 50 PSIG

Temperature (C)

180 ppm C₂F₆ 0 PSIG

• Flat sheet virus filtration membrane • Perfect mixing gas separation model (Weller-Steiner)

Feed Conc = 9.93% C₂F₆₊/CO₂

Temperature = 22 C Permeate Pressure = 0 PSIG Purge/Perm Flow Rate = 4.6 CO_2/C_2F_6 Selectivity = 1075

🗖 Data

- Model

60

Feed Pressure (psig)

 $<500 \text{ ppm } C_2 F_6$ ~25 PSIG

• Retention data on 24 Magnus 2.1 and 10 Magnus 2.2 devices At least 3 lots of each size • Defect detection sensitivity of air diffusion and binary gas tests verified • Binary gas test provides added level of retention assurance

*This point project from 5 min data

TESTING BY MILLIPORE DELIVERS HIGHER LRV ASSURANCE

- o Binary Gas Test is performed by Millipore as part of QC testing
- Viresolve Pro devices have been designed specifically to accommodate the Binary Gas Test



Binary Gas Test system dedicated for Viresolve Pro Micro devices

Binary Gas

Binary Gas Value (ppmv C₂F₆)

• Multi-function integrity test system dedicated for Viresolve Pro Modus/ Magnus devices

– Air diffusion

• Binary Gas Test provides superior sensitivity in detecting small defects that can affect LRV - Supported by theoretical and experimental data

SUMMARY

- Based on the principle differing gas permeabilities through liquid layer of a wetted membrane
- Key measurement is gas composition, not gas flow
- Demonstrated superior defect detection sensitivity compared to conventional air diffusion test
- Theoretical model closely predicts measured performance
- Enables the detection of defects in a wide range of device sizes from $<3 \text{ cm}^2$ to $> 1.5 \text{ m}^2$
- Used as 100% integrity test for Viresolve Pro devices

Acknowledgements

Nick Harrington, Steve DiTullio, Mike Doty, Mike Mills, Mani Krishnan, Greg Straeffer



- Pressure hold - Water flux
- Binary Gas

