3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Immunohistochemistry Procedure

Preparation Instructions

Citrate Buffer - 10 mM sodium citrate buffer, pH 6.0

 $1\times$ Tris/EDTA Solution - 10 mM Tris base with 1 mM EDTA solution and 0.05% TWEEN $^{\!@}$ 20, pH 9.0

2 M HCl Solution

Procedure

This immunohistochemical procedure is for deparaffinized sections, rehydrated with PBS.

Pretreat the sample using **one** of the following options:

- No treatment at all
- Place sample in Citrate Buffer, pH 6.0, microwave at 750 W for 20 minutes, and then cool the sample.
- Place sample in 1x Tris/EDTA Solution, microwave at 750 W for 20 minutes, and then cool the sample.
- Place sample in 2 M HCl Solution at room temperature for 10–20 minutes.
- Place sample in 0.1% trypsin and shake for 25 minutes at 37 °C.

Immunohistochemical Procedure

- 1. Incubate deparaffinized, pretreated sections in 3% (v/v) H_2O_2 in $1\times$ PBS at room temperature for 10 minutes and then wash the sections again.
- Incubate sections in blocking solution for 10 minutes.
- 3. Add primary antibody diluted in blocking solution and incubate the sections overnight at 2–8 °C, then wash sample with 1× PBS.
- 4. Incubate sections with peroxidase labeled polymer conjugated to a secondary antibody for 30 minutes followed by washing the sections with PBS.
- 5. Application of substrate solution (DAB or other suitable peroxidase substrate). Wash sample thoroughly under running tap water.
- 6. Counterstain the samples in Mayer's hematoxylin.
- 7. Dehydrate and mount samples.

TWEEN is a registered trademark of Croda International PLC.

JX.LPG.MAM 12/08-1