



Product Information

Monoclonal Anti-Human CD29

Clone P4C10

Purified Mouse Immunoglobulin

Product Number **C 9062**

Product Description

Monoclonal Anti-Human CD29, clone P4C10 (mouse isotype IgG1) is derived from the hybridoma produced by the fusion of mouse myeloma cell line NS-1/FOX-NY and splenocytes from RBF/Dn mice immunized with non-trypsinized HT1080 cells.

Monoclonal Anti-Human CD29 antibody may be used for the inhibition of cell-cell adhesion and for flow cytometry/immunocytochemistry identification and localization of CD29 expression.

CD29 is expressed on resting and activated leukocytes and on non-lymphoid cells. CD29 is the β_1 subunit of integrin. Integrins share the common structural feature of two non-covalently associated subunits, α and β . There are at least eight β -subunits which show strong homology to each other and 12 α -subunits which exhibit less sequence conservation. The amino acid sequence homology of the β -subunits is between 40 and 48%. All β -subunits contain 56 conserved cysteines arranged in four repeating units. The N-terminus is in the extracellular domain. The cytoplasmic domain is generally short, 40-50 amino acids long, except for β_4 . Phosphorylated serine and tyrosine have been reported in the cytoplasmic domains of the integrin β subunits. CD29 polypeptide β_1 has potential tyrosine phosphorylation sites and an apparent MW of 110 kDa in non-reduced SDS-PAGE. In contrast to earlier classification schemes, individual α -subunits have been shown to associate with more than one type of β -subunit. CD29 is known to associate with α_{1-6} CD49a-f. The non-covalent complexes of CD49a-f/CD29 are also known as very late antigens (VLA). Interaction of the cytoplasmic domain of CD29 and other β subunits with the cytoskeleton has been proposed as a mechanism for linking the cytoskeleton with the extracellular matrix. Integrin molecules bind a wide variety of extracellular matrix proteins such as collagen, laminin, fibronectin, thrombospondin, vitronectin and von Willebrand factor. Integrins also function in cell-cell adhesion. The cell surface ligands

for integrins are members of the immunoglobulin superfamily, Cell Adhesion Molecules (CAM).

Reagents

The product is provided as purified immunoglobulin in 0.01 M phosphate buffered saline, pH 7.4, containing 1% BSA and 0.1% sodium azide as a preservative.

Precautions and Disclaimer

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

Storage/Stability

Store at 2-8 °C. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Procedure

Indirect Immunofluorescent Staining

Reagents and Materials Needed but Not Supplied

1. a. Whole human blood collected by standard clinical blood evacuation tubes with EDTA, ACD-A, or heparin anticoagulant **OR**
b. Human cell suspension (e.g. peripheral blood mononuclear cells isolated on HISTOPAQUE[®] (Product Code 1077-1)).
2. Diluent: 0.01 M phosphate buffered saline (PBS), pH 7.4, containing 1% BSA, and 0.1% NaN₃.
3. Fluorochrome (FITC, PE, or Quantum Red[™]) conjugated anti-mouse secondary antibody diluted to recommended working dilution in diluent (FITC-Sheep Anti-Mouse IgG (whole molecule), F(ab')₂ fragment of Affinity Isolated Antibody, Product No. F 2883). Aggregates in conjugates should be removed by centrifugation immediately prior to use.
4. Isotype-matched, non-specific mouse immunoglobulin (negative control, Product No. M 5284).
5. 12 x 75 mm test tubes.
6. Adjustable micropipette.
7. Centrifuge.

8. Counting chamber.
9. 0.2% Trypan blue (Product No. T 0776) in 0.01 M phosphate buffered saline, pH 7.4.
10. 2% paraformaldehyde in PBS.
11. Whole blood lysing solution.
12. Flow cytometer.

Procedure

1. a. Use 100 μ l of whole blood **OR**
 b. Adjust cell suspension to 1×10^7 cells/ml in diluent. Cells should be >90% viable as determined by dye exclusion (trypan blue). For each sample, add 100 μ l or 1×10^6 cells per tube.
2. Add 5 μ l of monoclonal antibody to tube(s) containing cells to be stained. Vortex tube gently. Incubate the cells at room temperature (18 to 22 °C) for 30 minutes.
 Proper controls to be included for each sample are:
 - a. Autofluorescence control: 5 μ l diluent in place of monoclonal antibody.
 - b. Negative staining control: 5 μ l isotype-matched non-specific mouse immunoglobulin (Product No. M 5284) at the same concentration as test antibody.
3. After 30 minutes, add 2 ml of diluent to all tubes.
4. Pellet cells by centrifugation at 500 x g, for 10 minutes.
5. Remove supernatant by careful aspiration.
6. Resuspend the cells in 100 μ l of the fluorochrome conjugated secondary antibody at the recommended concentration. For the autofluorescence control, add 100 μ l of diluent. Incubate at room temperature (18 to 22 °C) for 30 minutes. Protect from light at this and all subsequent steps.
7. a. If whole blood is used, use lysing solution after incubation according to manufacturer's instructions.
 b. If a mononuclear cell suspension is used, proceed to Step. 8.
8. Add 2 ml of diluent to all tubes.

9. Wash as in steps 4-5.
10. After last wash, resuspend cells in 0.5 ml of 2% paraformaldehyde and analyze in a flow cytometer according to manufacturer's instructions.

It is advisable to run the appropriate negative controls. Negative controls establish background fluorescence and non-specific staining of the primary and secondary antibodies. The ideal negative control reagent is a mouse monoclonal or myeloma protein. It should be isotype-matched to the primary antibody, not specific for human cells and of the same concentration as the primary antibody. The degree of autofluorescence or negative control reagent fluorescence will vary with the type of cells under study and the sensitivity of the instrument used.

For fluorescence analysis of cells with Fc receptors, the use of isotype matched negative control is mandatory. In some systems it may be necessary to incubate the cells in 10-20% normal serum from the second antibody host species (at step 2 before adding monoclonal antibody) in order to decrease non-specific staining with the conjugated second antibody.

Product Profile

When assayed by flow cytometric analysis, using 5 μ l of the antibody to stain 1×10^6 cells, a fluorescence intensity is observed similar to that obtained with saturating monoclonal antibody levels. The percent population positive is also at the maximum percentage positive using saturating monoclonal antibody levels.

References

- Hemler, M. E., *Annu. Rev. Immunol.*, **8**, 365 (1990).
 Carter, W. G., et al., *J. Cell Biol.*, **110**, 1387 (1990).

JWM/KMR 08/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.